


ig N \_\_\_\_\_

Con'd from 3884 NB

2/18/95 wed

## MINIPREP DNA

- cfr 500  $\mu$ l of cells for 1 minute in an eppendorf cfr (centrifuge)
- removed supernatant and resuspended pellet in 100  $\mu$ l of 1X PEBI (SI) <sup>(saved)</sup>
- added 200  $\mu$ l of alkaline - SDS mix
- placed the tubes on ice for few minutes. (3-5 min)
- added 150  $\mu$ l of 7.5 M Ammonium Acetate
- Mixed the tubes by inverting
- cfr the tubes for  $\sim$  7-10 min.
- transferred 400  $\mu$ l supernatant to the new eppendorf tube
- added 800  $\mu$ l of ethanol to supernatant. Mixed tubes.
- incubated the tubes for  $\sim$  2 min. Spin.
- dissolved pellet in 50  $\mu$ l of TE + RNase A.
- applied 5  $\mu$ l to a 1% agarose gel.



SI = 0.9% glucose

25 mM Tris HCl (pH 8.00)

10 mM EDTA

alkaline - SDS mix = 1% SDS

0.1 N NaOH

To Page No. \_\_\_\_\_

Read &amp; Understood by me, \_\_\_\_\_

Date \_\_\_\_\_

Invent d by \_\_\_\_\_

Dat \_\_\_\_\_

Recorded by \_\_\_\_\_

4/12/95

Project No. \_\_\_\_\_

Book No. \_\_\_\_\_

TITLE \_\_\_\_\_

2

From Page No. \_\_\_\_\_

2/9/95 The

## Purification of m13 ssDNA

1. cfg 1.0 ml of \*infected cell culture for 2 min. (1 to 5 min)
2. Transferred 800.0  $\mu$ l to the new tubes  
(Pellet was saved for isolation of RF DNA)
3. cfg supernatant again to remove any residual cells
4. added 200.0  $\mu$ l of 20% PEG + 1.5 M NaCl. vortexed
5. Incubated tubes at room temperature for 5 min.
6. cfg tubes for 5 min. & discarded supernatant (sup.)
7. added 200  $\mu$ l of \*TE & vortexed really good.
8. cfg for ~ 1-2 min. (to remove any residual cell debris)
9. transferred sup. to the new tubes. (RNaseA can be added here)
10. added equal vol. of phenol / chloroform / isoamyl alcohol  
(25:24:1) Mixed well.
11. cfg 5 min.
12. removed the aq. (upper) layer to a new tube (be very careful)
13. added  $\frac{1}{10}$  vol. of 3M NaAc +  $2\frac{1}{2}$ -3 vol. of 95%  $\downarrow$
14. Incubated @  $-70^{\circ}\text{C}$  till 2/14/95.

$\left\{ \begin{array}{l} 20.0 \mu\text{l NaCl} \\ 600.0 \mu\text{l Et} \end{array} \right.$

TE ( $T_{10}E_1$ ) = 10 mM Tris-HCl pH 8.0 + 1 mM EDTA pH 8.0

infected cell culture = ① grew an *E. coli*  $F'$  strain to an OD of 0.4 in 2xYT   
 $\downarrow$  next pag  
 $F'$  = Fertility Factor: codes for *tra* genes & pilis to allow infection of the   
 m13 Phage. (transfer of DNA)

cont'd To Page No

Witnessed &amp; Understood by me,

Date

Invented by

Date

R c r d b y

4/12/95

ig N \_\_\_\_\_

2/9/95.

Cont'd \_\_\_\_\_

- ② Inoculated 1.0 ml of the cells ~~not~~ with the phage
- ③ Incubated the phage infected cells at 37°C for 5 hours. Now supernatant <sup>was</sup> ~~can be~~ processed for isolation of ssDNA & cells (pellet) were ready for isolation of RF (Replicating Form) dsDNA.

2/14/95 Tues.

1. Poured 0.8% Agarose gel (250 ml Volumes) in 1X TAE Buffer

2g Agarose

250.0 ml 1X TAE Buffer

- weight the flask
- boiled for 4.20 min. (brought up the)
- weight the flask & adjusted Volume to before boiling with • distilled water
- poured it on the plate.

2. 2X YT

added: 19.3g Tryptone } brought total Volume (TV)  
 12g Yeast Extract } to 1200 ml with water  
 12g NaCl }

- made 5 aliquots (1) 500.0 ml (2) 250 ml (3) 100.0 ml  
 (4) 100.0 ml  
 (5) 100.0 ml
- autoclaved at low pressure for 20 min.

3. 2X YT Top (soft) Agar

added: 0.35 g Agar } made 3 different  
 50.00 ml 2X YT } aliquots.

- autoclaved at low pressure for 20 min. (same as 2X YT)

To Page No. \_\_\_\_\_

Read &amp; Understood by me,

Dolan

Date

4/12/95

Invented by

Recorded by

Date

4/12/95

From Pag No. \_\_\_\_\_

4. started 2, 10 ml cultures

#1 CJ236 in 2XYT + Cm 5  $\mu$ g/ml +/- phage from T. nea / m13  
cloning at 37°C air-shaker

10 ml CJ236 in 2XYT

50  $\mu$ l Cm 5  $\mu$ g/ml  $1000 \mu\text{g} = \text{mg/ml} = 1000 \mu\text{l}$

#2 10 ml culture of T. nea / pTTG and T. nea / pTTG 118 in LB  
+ 100  $\mu$ g/ml Ampicillin at 30°C air-shaker

added: 10 mL

5  $\mu$ l.

5. centrifuged (cf) ssDNA from DH5 $\alpha$  F'IQ (tubes left at -70°C, 2/9/95) for 10 min. at room temperature.

• discarded supernate

(a) Rinsed the precipitate (ppt.) with 70% EtOH. (removed any residual EtOH with another quick spin.

(b) Dried the DNA pellet at 55°C heat block

(c) Dissolved the DNA in 50.0  $\mu$ l of TE

(d) applied 5.0  $\mu$ l on a 0.8% Agarose gel (made today).  
m13 ssDNA was used as a control

T Pag N

Witnessed &amp; Understood by me,

Date

Inv nted by

Dat




4/12/95

R c rd d by

4/12/95

ig N . \_\_\_\_\_



Run 140V ~ 2 hrs.

amp 4/12/95

To Pag N . \_\_\_\_\_

sed & Understood by me,

*J. Polansky*

Date

4/12/95

Invented by

Record d by

*[Signature]*

Date

4/12/95

6

Project No. \_\_\_\_\_

Book No. \_\_\_\_\_

TITLE \_\_\_\_\_

From Page No. \_\_\_\_\_

2/15/95 Wed.

(+) strand (ssDNA) lot # ED5702 260  $\mu$ g/ml  
RF strand (dsDNA) lot # CC3111 5  $\mu$ g/18.4  $\mu$ l

calculation: ssDNA = 260  $\mu$ g/ml = ng/ $\mu$ l

$$\frac{260 \cancel{\mu\text{g/ml}}}{1000 \cancel{\text{ng}/\mu\text{g}}} \cdot 1000 \text{ ng}/\mu\text{g} \cdot \text{ml} \cdot 1000 \cancel{\mu\text{l}} = 0.260 \mu\text{g}/\mu\text{l}$$

$$\frac{1000 \text{ ng}/\mu\text{g} \cdot (0.260 \mu\text{g})}{260 \text{ ng}/\cancel{\mu\text{g}}} = 260 \text{ ng}/\mu\text{l}$$

$$\frac{260 \text{ ng}/\cancel{\mu\text{g}}}{x \mu\text{l}} = 100 \text{ ng}$$

$$\left\{ \begin{array}{l} 260 \left( \frac{1}{2.6} \right) = 100 \text{ ng} \\ \text{or} \\ \frac{260}{2.6} = 100 \text{ ng} \end{array} \right.$$

for 2.6 total or final volume  
you need 1.0  $\mu$ l DNA

$$\therefore \frac{1 \mu\text{l DNA} (260 \text{ ng}/\mu\text{l})}{1.6 \mu\text{l TE}} = 2.6 \mu\text{l}$$

for 100 ng/ $\mu$ l, } 2.0  $\mu$ l DNA (260 ng/ $\mu$ l)  
multiply by 2 } 3.2  $\mu$ l TE

dsDNA = 5  $\mu$ g/18.4  $\mu$ l.

$$1000 \text{ ng}/\mu\text{g} \times 5 \mu\text{g} = \frac{1000 \text{ ng} (5 \cancel{\mu\text{g}})}{\cancel{\mu\text{g}}} = 5000 \text{ ng}/18.4 \mu\text{l}$$

$$\frac{5000 \text{ ng}}{18.4 \mu\text{l}} = \frac{272 \text{ ng}/\mu\text{l}}{x \mu\text{l}} = 2.72 \text{ ng}$$

for 2.7 total volume you need  
1.0  $\mu$ l DNA

Total Volume (TV)

T Pag N

Witnessed & Understood by me,



Date

4/12/95

Invent d by

R c rded by

Dat

4/12/95

Pag No. \_\_\_\_\_

	RF (ds)	Tube # 1	Tube # 2
DNA	1.0 $\mu$ L $\times 3 = 3.0 \mu$ L		(+) ssDNA DNA = 1.0 $\mu$ L $\times 3 = 3.0 \mu$ L
TE	1.7 $\mu$ L $\times 3 = 5.1 \mu$ L		TE = 1.6 $\mu$ L $\times 3 = 4.8 \mu$ L
TV	2.7 $\mu$ L $\times 3 = 8.1 \mu$ L		TV = 2.6 $\mu$ L $\times 3 = 7.8 \mu$ L

Tube # 1, 2, 3, 4 of RF (dsDNA)

	(1) Alu I	(2) Hind III	(3) Sau 3A I	(4) Bam HI	
H <sub>2</sub> O	16.0 $\mu$ L				(all 4 tubes w/ 16.0 $\mu$ L)
10x Buffer	2.0 $\mu$ L				
DNA	1.0 $\mu$ L				
	(React 1; React 2; React 4; React 3)				
Alu I	+	-	-	-	
Hind III	-	+	-	-	
Sau 3A I	-	-	+	-	
Bam HI	-	-	-	+	

Tube # 1, 2, 3, 4 of + (ssDNA) same order as RF

2 tubes were set-up for uncut, 1 with RF &amp; 2nd with (+)

- each tube added 16.0  $\mu$ L H<sub>2</sub>O  
2.0  $\mu$ L React 2 10x buffer  
1.0  $\mu$ L DNA

Put all 10 tubes in

- ran the sample (all 10) on a gel next morning.  
(0.8% agarose gel, 147 volts)
- picture of the gel is on the next pg (pg # 8)

To Pag No. \_\_\_\_\_

Used &amp; Understood by me,

J. Olamp

Date

4/12/95

Inv. nted by

Record d by

J. Olamp

Date

4/12/95

Project No. \_\_\_\_\_

Book No. \_\_\_\_\_

TITLE \_\_\_\_\_

B

From Page No. \_\_\_\_\_


tube# 1 T. nea / pTTC

1.0 ml

tube# 2 T. nea / pTTC

1.0 ml

- Cfg. for 1 min. at room temperature
- discarded supernate and added: 100  $\mu$ l S1 to the pellet. mixed  
200  $\mu$ l S2 lysis put both tube  
ice.  
150  $\mu$ l S2 with RNASE A
- Cfg. for 5 min. at 4°C
- transferred 400  $\mu$ l of supernatant to the new tubes.
- added 800.0  $\mu$ l E+OH to the supernatant
- put both tubes in the fridge till tomorrow (2/16/95)



arp 4/12/95

(+) Sal3A1 - gel shift (did not cut  
but binded)

To Page No

Witnessed &amp; Understood by m ,

Date

Investigated by

Date

4/12/95

Recorded by

4/12/95



ag N .

2/21/95 TUE

DIGEST T.nea/pSPORT with SstI & SphI

DIGEST M13 mp18 & M13 mp19 w/ SstI & SphI

M13 mp 18 RF (0.44 ug/ul) } cut 500.0 ng  
 M13 mp 19 RF (350.0 ug/ml)

$$\rightarrow 1000 \text{ ng/ug} \times 0.44 \text{ ug} = 440 \text{ ng}$$

$$\frac{500 \text{ ng}}{440 \text{ ng}} = 1.1 \text{ ul}$$

$$\rightarrow 1000 \times 0.350 \text{ ug/ul} = 350 \text{ ng}$$

$$\frac{350}{350} = 1.4 \text{ ul}$$

mp 18

mp 19

H <sub>2</sub> O - 35.0 ul	H <sub>2</sub> O - 35.0 ul
10x buffer - 2.0 ul ← REact 2 →	10x buffer 2.0 ul
500 ng DNA - 1.1 ul	DNA 1.4 ul
1 ul SstI - 2.0 ul	SstI 2.0 ul
40.0 ul	40.0 ul

T.nea / pSPORT

H<sub>2</sub>O - 81.0 ul  
 10x buffer - 10.0 ul  
 ng/ul DNA - 4.0 ul  
 SstI - 5.0 ul  
 100.0 ul

- Incubated all 3 tubes @ 37°C for 1/2 hour
- Made 0.8% agarose gel  
 250.0 ml TE buffer  
 2.0 g Agarose
- boiled for 4.0 min.
- added 12.0 ul E. Bio mide
- Poured the gel.

To Page No. \_\_\_\_\_

s d & Understood by me,

Date

Invented by

Dat

50 okup

4/12/95

Recorded by

4/12/95

From Page N \_\_\_\_\_

added: 2.0  $\mu$ l of 1M KCl  $2 \mu$ g = 2000 ng  
to 40.0  $\mu$ g (mp 18 & 19)

5.0  $\mu$ l of 1M KCl  
to 100.0  $\mu$ g (pSPORT)

added ~~Sph~~ Sph I - 2.0  $\mu$ l mp 18  
2.0  $\mu$ l mp 19  
5.0  $\mu$ l pSPORT

- Incubated @ 37°C for 1/2 hour
- put the tubes in the fridge till
- ran samples on the gel ~~to~~ on 2/22/95



arp 2/22/95 ②

M13mp18 & mp19 RF D  
are ds, supercoiled forms  
the DNAs of phages M13  
& 19. Using this vector  
foreign DNA can be  
inserted into the mul  
cloning site in an  
oriented fashion.

To Page N \_\_\_\_\_

Witnessed &amp; Understood by me,

Date

4/12/95

Inv nt d by

R c rd d by

Date

4/12/95

19 N \_\_\_\_\_

2/22/95

1. grow cells overnight (O/N) 10.0 mL

= 9.0 mL (1.0 mL in ea. nine tubes). Each tubes labelled DH10B

• Quick freeze all nine tubes in a powdered  
dry ice.

Plz C/T. nea

2/22/95 BJS

LB + AP100

### GENE CLEAN

2) Did electrophoresis of yesterday's DNA (2/21/95)

M13 mp 18 and M13 mp 19 and pSPORT

b) Took the picture of the gel

c) cut off mp 18 fragment, mp 19 fragment & pSPORT fragment from the gel & transformed the gel w/ <sup>pa</sup> DNA into the separate eppendorf tubes.

d) added 700.0  $\mu$ L NaI to each <sup>2</sup> tubes. Vortexed mp 18 & mp 19 tubes.

e) Incubated both tubes @ 55°C to melt agarose. mixed ~~after~~ after incubation.

f) added 5.0  $\mu$ L glass milk to both tubes.

g) Incubated both tubes on ice for 5 min.

h) Cfg. both tubes (quick spin)

i) discarded supernate

j) added 500.0  $\mu$ L New Wash buffer

k) discarded supernate & again added 500.0  $\mu$ L New Wash buffer.  
washed both tubes 3 times.

l) added 10.0  $\mu$ L dH<sub>2</sub>O to the tubes. mixed well by vortexing. 55°C for 2-5 min

m) set up Ligation

ligation

H<sub>2</sub>O = 12.0  $\mu$ L

ligase) 5x Buffer = 4.0  $\mu$ L

mp 18 DNA = 2.0  $\mu$ L

(1  $\mu$ L/ $\mu$ L) Ligase = 2.0  $\mu$ L

TV = 20.0  $\mu$ L

H<sub>2</sub>O = 12.0  $\mu$ L

5x buffer = 4.0  $\mu$ L

mp 19 DNA = 2.0  $\mu$ L

Ligase = 2.0  $\mu$ L

TV = 20.0  $\mu$ L

n) Incubated both tubes overnight @ room temperature (cond)

T Pag No.

ed & Understood by m ,

Dat

Inv nted by

Dat

2/22/95

4/12/95

R cord d by

4/12/95

From Page No. \_\_\_\_\_

(cont'd.)

T.nea/Ptarc E. pttc

- 1.0 mL of ea.
- Cfg.
- discarded supernate
- added 100.0  $\mu$ l SI mixed well
- Incubated on ice for few min.
- added 200.0  $\mu$ l S2 14SIS
- Incubated on ice for few min.
- added 150.0  $\mu$ l S3 w/ RNAase A
- Cfg. for 7.0 min. @ 4°C


T Pag

Witnessed &amp; Understood by me,

Dat

Invented by

Date



4/12/95

R c rd d by

4/12/95

Project No. \_\_\_\_\_

Book No. \_\_\_\_\_

TITLE \_\_\_\_\_

20

From Page N \_\_\_\_\_

2/28/95 TUE

I. set up digest DNA ppt.

① M13mp 18, ② mp 19 and ③ T. nea/pSPORT

1. - To ea. 3 added 100.0  $\mu$ L TE } to ppt.  
- " " 10.00  $\mu$ L NaAc } DNA  
- " " 300.00  $\mu$ L EtOH }

2. Incubated on dry ice for ~5 min.

3. Cfg. for 10 min. @ room temp. (no ppt.)

4. no ppt., added 2.0  $\mu$ L (carrier molecule) Yeast tRNA. Vortexed

5. incubated on dry ice for ~5 min.

6. Cfg for 10 min. @ room temp. (Supernate saved) Pellet was saved on mp 18

7. added 200.0  $\mu$ L 70% EtOH to the pellet

8. Cfg. discarded supernate, air dried by putting tubes on the heat block.

II. DIGEST set-up. (to map Bam HI site)- cut T. nea/pSPORT with Hind III, Bam HI, Xba, NOT I, Sst, Eco R  
SeparateH<sub>2</sub>O - 13.0  $\mu$ Lbuffer - 2.0  $\mu$ LT. nea/pSPORT DNA - 3.0  $\mu$ Lenzyme - 2.0  $\mu$ LTV = 20.0  $\mu$ LEnzymes - Hind III, Xba, Sst had REact 2 1  
buffer.- Bam HI, NOT I, Eco RI had REact:  
- buffer.Control: H<sub>2</sub>O - 13.0  $\mu$ L(REact 2) buffer - 2.0  $\mu$ LDNA - 3.0  $\mu$ Lfor  
separate  
enzymes.

- Incubated @ 37°C

- ran on the gel on 3/1/95 (Wed)

Picture shown pg 21

T Pag N

Witnessed &amp; Understood by me,

Date

Invent d by

Dat

4/12/95

R cord d by

4/12/95

ag N . \_\_\_\_\_

# I DIGEST Set-up


- Cut ptra / T. nea w/ Sst I enzyme

H<sub>2</sub>O - 11.0  $\mu$ l  
 (2) buffer - 2.0  $\mu$ l.  
 T. nea DNA - 2.0  $\mu$ l.  
 Sst I - 5.0  $\mu$ l.  
 TV = 20.0  $\mu$ l.

Control: H<sub>2</sub>O - 11.0  $\mu$ l  
 buffer - 2.0  $\mu$ l.  
 DNA - 2.0  $\mu$ l


- Incubate @ 37°C

- Lam DNA on a gel on 3/1/95 (wed) picture shown below



from pg. 20

GRAPH 1:



GRAPH 2: pg 23 of this book.

To Page No. \_\_\_\_\_

s d & Understood by m ,

*Polansky*

Date

4/12/95

Invented by

Recorded by

*Swan*

Date

4/12/95

Project No. \_\_\_\_\_

Book No. \_\_\_\_\_

TITLE \_\_\_\_\_

22

From Page No. \_\_\_\_\_

3/12/95

① 1 kb ladder ② T-neal/pSPORT uncut, ③ Sst, ④ Sst/sph, ⑤ sph, ⑥ 1 kb ladder (from 263 added loading dye, electrophoresis @ 190 v

- digested double digested BamHI/sphI (to map the Bam site T-neal

H<sub>2</sub>O - 14.0  $\mu$ l

(REACT) buffer - 2.0  $\mu$ l

(T-neal/pSPORT) DNA - 2.0  $\mu$ l

(Bam/sph) enzyme - 1.0  $\mu$ l ea.


TV = 20.0  $\mu$ l.

control: - H<sub>2</sub>O - 14.0  $\mu$ l

(uncut) buffer - 2.0  $\mu$ l

DNA - 2.0  $\mu$ l

Incubated @ 37°C for 30 min. (15 min.)



T Page 1

Witnessed & Understood by m ,

Date

Inv nt d by

Dat

*[Signature]*

4/12/95


R c rd d by

4/12/95

*[Signature]*

Tag N. (from pg. 21)

GRAPH 2:




Bam I/Sph I

- (A)
- 4.0 kb
- 700 bp
- 1 kb
- 1.5 kb

Bam I/Sph I

- (B)
- 4.0 kb
- 2.2 kb
- 300 bp
- 700 bp




Cloning mp 18 w/ T. nea / pSPORT & mp 19 w/ T. nea / pSPORT


3/2/95

~~3/2/95~~

Thurs.



Picture before cutting the DNA fragment.



DID GENE CLEAN

To Page No. \_\_\_\_\_

Used & Understood by me,

*[Signature]*

Date

4/12/95

Invented by

Record d by

*[Signature]*


Date

4/12/95



Proj ct No. \_\_\_\_\_  
 B k N . \_\_\_\_\_

ag N (from pg 21)  
 RAPH 2:



*Bam*/sph:  
 (A)  
 4.0 kb  
 700 bp  
 1 kb  
 1.5 kb

*Bam*/sph:  
 (B)  
 4.0 kb  
 2.2 kb  
 300 bp  
 700 bp


3.2  
 1.7  
 1.5

Cloning mp 18 w/ T.nea/psport & mp 19 w/  
 T.nea/psport


3/2/95

~~3/2/95~~

Thurs.



Picture before  
 cutting the  
 DNA  
 fragment.



Picture after cutting  
 the DNA fragment

DID GENE CLEAN

To Page No. \_\_\_\_\_

sed & Understood by me,

*[Signature]*

Date

4/12/95

Invented by

*[Signature]*

Date

4/12/95

Recorded by

*[Signature]*

From Page No. \_\_\_\_\_

GENE CLEAN

- (out) mixed mp 18 with T.nea pSPORT cut w/ Sst I/sph } 1 tube
- " mp 19 with " " " " } 2 tube

- added 700.0  $\mu$ l NaI to each 2 tubes. Vortexed
- put the tubes in  $55^{\circ}\text{C}$  heat block to melt agarose
- after agarose melted, added 5.0  $\mu$ l glass milk to both tubes
- incubated both tubes on ice for 5.0 min.
- c.f.g. both tubes (quick spin)
- discarded supernate & washed pellet 3 x with New Wash b
- added 14.0  $\mu$ l d  $\text{H}_2\text{O}$  to each tube
- quick spinned, discarded pellet & saved supernate

Set-up Ligation

- (mp 18) (mp 19) DNA - 14.0  $\mu$ l
- (ligase) 5x buffer - 4.0  $\mu$ l
- ligation - 2.0  $\mu$ l.
- TV - 20.0  $\mu$ l.

- incubated both

Transformation Cells

- (1) 100.0  $\mu$ l Competent
- 3.0  $\mu$ l DNA (from ligation)

- (2) incubated on ice for 30 min.

- (3) heat shocked @  $42^{\circ}\text{C}$   $\text{H}_2\text{O}$  bath for 35 sec.

- melted 0.7% 2x YT top agar, added 4.0  $\mu$ l to 6 different glass tubes & put the tubes @  $55^{\circ}\text{C}$  heat block

Didn't work

T Page No. \_\_\_\_\_

Witnessed &amp; Understood by m ,

Dat

Invented by


Dat

Rec rded by

4/12/95


ag N \_\_\_\_\_  
 Repeat

3/7/95 TUE



3/7/95

After taking picture on looking @ the gel, ~~M±~~ \* M±3 mp18 and M±3 mp19 is @ the 7.2 Kb ~~cut~~ which was cut with Sph I. ~~we decided to cut~~ We planned on cutting mp18 and mp19 with Sst I. The gel <sup>picture</sup> below shows mp18 & mp19 before & after cutting the DNA fragments. After cutting the fragments performed Gene CLEAN



3/7/95



To Page No. \_\_\_\_\_

Used & Understood by m ,

*Bokey*

Date

4/12/95

Invented by

Recorded by *Duan*

Date


4/12/95

From Page N \_\_\_\_\_

labelled 2 tubes, 1 w/ mp 18, & 2<sup>nd</sup> w/ mp 19

1. to the DNA w/ agarose gel, added 100.0  $\mu$ l NaI
  2. put the tubes @ 52°C heat block to melt agarose. vortexed constant
  3. added 5.0  $\mu$ l glass milk to both tubes - mixed
  4. incubated on ice for 5 min.
  5. (fg. (quick spin) @ room temp.
  6. discarded supernate, added 500.0  $\mu$ l New wash buffer
  7. discarded supernate, washed pellet 3x with New wash buffer
  8. after washing 3x, added 14.0  $\mu$ l dH<sub>2</sub>O to the pellet (discarded supernate) (mixed)
  9. incubated @ 52°C for 5 min.
  10. discarded pellet & saved supernate for ligation.
- (could this on 3/8/95 wed.)

Purification of m13 ssDNA (T-neo 2kb [Sph1]/mp19) from pg. 1



To Page N \_\_\_\_\_

Witnessed &amp; Understood by me,

Date

4/12/95

Invited by

Recorded by

Date

4/12/95

ag N .

cell growth & infection

- Grew an E. coli F' strain to an OD of 0.2-0.4 in 2x YT
- Inoculated 1-2 ml of the cells w/ the phage. (added 10.0  $\mu$ l from a liquid phage stock & added to cells)
- Incubated the phage infected cells @ 37°C for 5-7 hours.
- The supernate can now be processed for isolation of ssDNA & the cells can be processed for the isolation of Replication Form (RF) dsDNA.

Purification of ml3 ssDNA

- transferred 1.0 ml culture of infected cell to 4 different eppendorf tubes
- cfg 4 tubes for 2 min.
- transferred supernate to the new tubes & saved pellet from 1 tube (out of 4 tubes) for isolation of RF DNA
- Spinned the supernate again & transferred the supernate to the new tubes (done to remove any residual cells remained behind)
- passed the supernate through a 0.45  $\mu$  filter as to remaining cells (done when performing site-directed mutagenesis)
- added 200.0  $\mu$ l of 20% PEG + 1.5 M NaCl. Vortexed
- Incubated tubes for 15 min @ room temperature (or overnight @ 4°C)
- cfg for 10 min in a  $\mu$ cfg. @ room temp.
- discarded supernate & briefly spinned the tubes to remove the residual soln from the side of the tube (removed as much <sup>supernate</sup> as possible)
- added 200.0  $\mu$ l TE. Vortexed
- cfg for 2 min. to remove any residual cell debris.
- Transferred supernate to the new tube. (added 5.0  $\mu$ l RNase I to remove any residual nucleic acid from the prep. Beryonase will remove both RNA & DNA very efficiently.)
- added equal volume of phenol / chloroform / isoamyl alcohol mixed well
- cfg for 5.0 min.
- transferred the upper layer to a new tube (BE CAREFUL NOT TO DISTURB WHITE INTERFACE OR REMOVE ANY PHENOL)
- added 20.0  $\mu$ l NaAC & 600.0  $\mu$ l EtOH
- Incubated @ -70°C for 5-15 min. (we left @ -70°C overnight)

To Page No. \_\_\_\_\_

ssed &amp; Und rst od by me,

Date

Invent d by

Date

Recorded by

J. J. J. J.

4/12/95

4/12/95

Project No. \_\_\_\_\_

28

Book No. \_\_\_\_\_ TITLE \_\_\_\_\_

From Page No. \_\_\_\_\_

3/8/95 wed

- cftg the samples for 10-15 min.
- discarded the supernate & rinsed the pellet w/ 70% EtOH
- dried the pellet @ 55°C heat block or @ room temperature
- dissolved the DNA in 50.0  $\mu$ l TE.

Shm Rxn

3/8/95 wed

Annealing Rxn.

+ Primer (2899)

- Primer (2899)

H<sub>2</sub>O - 3.0  $\mu$ l4.0  $\mu$ l5x Buffer 2.0  $\mu$ l2.0  $\mu$ l19.5% SS DNA 4.0  $\mu$ l4.0  $\mu$ l(200mg/4l Kinasol) Oligo 1.0  $\mu$ lTV 10.0  $\mu$ l10.0  $\mu$ l

Incubated @ 70°C - 75°C for 2 min. (to eliminate non-spf. bin)  
" @ 37°C - 40°C for 2 min.

Synthesis RxnAnnealing Rxn - 10.0  $\mu$ l.5mL 10x buffer - 2.0  $\mu$ lH<sub>2</sub>O - 6.0  $\mu$ lT<sub>4</sub>/T<sub>7</sub> DNA pol - 1.0  $\mu$ lT<sub>4</sub> DNA ligation - 1.0  $\mu$ lTV - 20.0  $\mu$ l.

Incubated @ 37°C for 10 min.

Synthesis rxn - 2.5  $\mu$ lTE - 8.0  $\mu$ lloading dye - 1.0  $\mu$ l

- ran the sample on the gel
- the picture on the next page, # 29.

T Pag No

Witnessed &amp; Underst d by m ,

Date

Invent d by


Dat

4/12/95

Rec rd d by


4/12/95

Tag N \_\_\_\_\_



(con'd from pg. 28)  
 +2899 (w/ primer) oligo forms a ds DNA.  
 +2899 fragment looks brighter because  
 Et. Bromide binds to it better. -2899  
 primer binds but <sup>(GWA)</sup> it does not hold strongly  
 ∴ the DNA fragment looks fainter or  
 light, less Et. Bromide is able to bind.

(con'd on pg. 41)



14, 19, 20, 21, 22, 23, 24, 25, 26, 27, 28, 29, 30, 31, 32, 33, 34, 35, 36, 37, 38, 39, 40, 41, 42, 43, 44, 45, 46, 47, 48, 49, 50, 51, 52, 53, 54, 55, 56, 57, 58, 59, 60, 61, 62, 63, 64, 65, 66, 67, 68, 69, 70, 71, 72, 73, 74, 75, 76, 77, 78, 79, 80, 81, 82, 83, 84, 85, 86, 87, 88, 89, 90, 91, 92, 93, 94, 95, 96, 97, 98, 99, 100

Ligation from 3/1/95 (pg. 26)

H <sub>2</sub> O -	8.0 uL	H <sub>2</sub> O -	8.0 uL
5X buffer -	4.0 uL	5X buffer -	4.0 uL
mp 18 -	2.0 uL	(vector) mp 19 -	2.0 uL
insert -	4.0 uL	insert -	4.0 uL
ligation -	2.0 uL	ligation -	2.0 uL
TV -	20.0 uL	TV -	20.0 uL

- Incubated both samples for 1 hour @ room temp.

100.0 uL Competent cells }  
 3.0 uL DNA } xfection cells.

ran mp 18 on 3/10  
 (used DNA from  
 3/10/95 again on  
 3/15/95 pg. 32)

xfection

10%	mp 18	mp 19
90%	mp 18	mp 19
Control		

To Page No. \_\_\_\_\_

Read & Understood by me,

Date

Inv nted by

Date

*[Signature]*

4/12/95

Recorded by

*[Signature]*

4/12/95

From Page No. \_\_\_\_\_

— 10% mp 18 / mp 19

added → 4.0 mL 2x YT TOP Agar

100.0  $\mu$ L X-Gal 4%

5.0  $\mu$ L IPTG 200 mM (inducer = repressor gives tighter affinity)

60.0  $\mu$ L lawn cells

10.0  $\mu$ L x fraction cells (after heat shock for 35 sec.)

— 90% mp 18 & mp 19.

Same way as 10%

— Control

100.0  $\mu$ L X-Gal

5.0  $\mu$ L IPTG

60.0  $\mu$ L lawn cells.

no 3:30 10%

T Pag N

With ss d &amp; Und rsto d by me,



Dat

4/12/95

Inv nt d by

R c rd d by



Dat

4/12/95



ge No. \_\_\_\_\_

3/14/95 TUE


~~Repro~~: miniprep

1.0 ml culture of T-nealmp 19 grown for 5 hours @ 37°C in  
 6 different glass tubes  
 transferred 1 ml cell to the 6 different labelled eppendorf tubes.  
 cfg all 6 tubes for 2 min. @ room temp.  
 removed supernate & saved in different tubes  
 added 100  $\mu$ l S1 mixed well  
 added 200  $\mu$ l S2. put the tubes on ice (mixed by inverting)  
 added 150  $\mu$ l 7.5 M  $\text{NH}_4\text{OAc}$   
 incubated on ice for 5 min.  
 cfg for 7 min. @ room temp (4°C) NOTE: cfg in 4°C room was taken away for repair: used @ RT  
 transferred supernate (400.0  $\mu$ l) to the new 6 labelled tubes  
 added 800  $\mu$ l of EtOH to the 400  $\mu$ l of supernate (mixed well)  
 incubated @ -70°C for 30 min.  
 cfg for 2 min. @ room temp (discarded supernate)  
 rinsed w/ 70% EtOH (removed supernate)  
 added 50.0  $\mu$ l TE to the pellet.

$\text{H}_2\text{O}$ - 7.0 $\mu$ l	x 6	= 42.0 $\mu$ l
buffer - 2.0 $\mu$ l	x 6	= 12.0 $\mu$ l
Eco47III - 6.0 $\mu$ l	x 6	= 6.0 $\mu$ l
TV		60.0 $\mu$ l

added 10.0  $\mu$ l DNA<sup>+</sup> to each 6 tubes.

the map is on next page # 32. Fragments  
 on all 6 tubes are in still parent,  
 y haven't gone into the mutant.  
 \* I tried miniprep again next day.  
 (started)



ed &amp; Understood by m ,

Date

Invent d by

Date

Polans

4/12/95

R corded by

4/12/95

32


Project No. \_\_\_\_\_

Book No. \_\_\_\_\_

TITLE \_\_\_\_\_

From Page No. \_\_\_\_\_

3/15/95 Wed




parent  
Eco 47 III  
~ 8-9 kb  
0.3 kb

→ (most probably won't see fragment because too small & too light)

mutant  
4 kb  
4.7 kb  
0.3 kb

DNA from date 3/10/95



H<sub>2</sub>O = 6.0 mL.  
React6 buffer = 2.0 mL.  
mpi8 DNA = 10.0 mL  
Sst/Sph 1.0 mL ea.  
TV 20.0 mL

H<sub>2</sub>O = 6.0 mL.  
buffer = 2.0 mL  
mpi8 DNA = 10.0 mL  
Sst/Sph = 1.0 mL ea.  
TV 20.0 mL

from pg 29  
& ran again  
on 3/15

arp  
4/12/95

T Pag

Witnessed & Understood by me,

Date

Inventor by

Date

4/12/95

Recorded by


4/12/95

*[Signature]*

*[Signature]*

ag No. \_\_\_\_\_

- Incubated both tubes @  $37^{\circ}\text{C}$  for 30 min.
- added 2.0  $\mu\text{L}$  loading dye to each tube
- ran both on a gel
- took picture



3/15/95 T. nea/mp

1.0 ml T. nea (sph I) / mp 19 + 2899 + Sau 3AI grown for 5 hours @  $37^{\circ}\text{C}$  in 10 different glass tubes  
 after 5 hours transferred 1.0 ml culture to the 10 labelled eppendorf tubes  
 cfd all 10 eppendorf tubes @ room temperature for 2 min.  
 removed supernate & saved  
 put all 10 tubes w/ pellet & all 10 tubes w/ supernate @  $-70^{\circ}\text{C}$  overnight or until 3/16/95 Thursday.

TE: Brian had to leave @ 4:30 pm & this was a point to stop @.

To Page No. \_\_\_\_\_

Read &amp; Understood by me,

Date

4/12/95

Inv. nted by

Recorded by

Date

4/12/95

Project No. \_\_\_\_\_

Book No. \_\_\_\_\_ TITLE \_\_\_\_\_

From Page No. \_\_\_\_\_

3/16/95 Thurs.

con'd from page 33 3/15/95 wed. MINIPREP

- took the pellet out from  $-70^{\circ}\text{C}$  (10 eppendorf tubes)
- added 100  $\mu\text{L}$  S1 mixed well
- added 200  $\mu\text{L}$  S2 put all 10 tubes on ice mixed
- added 150  $\mu\text{L}$  7.5 M  $\text{NH}_4\text{OAc}$
- incubated on ice for 5 min.
- cfg all 10 tubes for 5 min. @ room temp. ( $4^{\circ}\text{C}$ )
- transferred 400  $\mu\text{L}$  of supernate to the new 10 labelled tubes
- added 800  $\mu\text{L}$  EtOH Mixed well
- incubated all 10 tubes for 30 min. @  $-70^{\circ}\text{C}$ .
- cfg & discard for 2 min. @ room temp.
- discarded supernate & washed pellet with 70% EtOH.
- added 50  $\mu\text{L}$  TE to all 10 tubes w/ pellet

	tubes	
$\text{H}_2\text{O}$	1.0 $\mu\text{L}$ x 10	= 10.0 $\mu\text{L}$
buffer	2.0 $\mu\text{L}$ x 10	= 20.0 $\mu\text{L}$
ECO 47 III	1.0 $\mu\text{L}$ x 10	= 10.0 $\mu\text{L}$
TV		= 100.0 $\mu\text{L}$

- added 10.0  $\mu\text{L}$  from TV to all other 9 tubes
- added 10.0  $\mu\text{L}$  DNA to each 10 tubes
- incubated @  $37^{\circ}\text{C}$  for 30 min.
- added 2  $\mu\text{L}$  loading dye
- ran all 10 samples on a gel for 1 hour @ 190 V
- took a picture

picture on pg 35

To Pag No

Witnessed &amp; Understood by me,



Date

4/2/95


Invented by

Recorded by

Dat

4/12/95

ag N \_\_\_\_\_

parent & mutant should look like  
ECO 47 III

Parent

mutant

8.9 kb

4.5 kb

0.3 kb

4.4 kb

0.3 kb

3/16/95

1/12/95

may probably be  
too light to seeNOTE: In this ~~we~~ we could see parent & some mutant. mutant  
is seen on # 5, 6, 7, 8

Po

To Page No. \_\_\_\_\_

sed &amp; Und rstood by me,

Date

Invented by

Date

Dolans

4/12/95

Recorded by

Duan

4/12/95