

Project N _____
 Book No. _____

11/21/94 101

PMC9 / Tag

Tag No. _____

reaction: PMC9 amplification, using dNTPs & primers # 2722 & 2729.

reaction: buffers < Klenow
 Deepvent

cycling:

use dNTP

use primer each

use template PMC9 / Act II

use Mg

80 (94° 30", 56° 30", 72° 30")

prepared 15 x w/o enzyme → added separately in 1x buffer

x buffer (D.V)	75	(K7) 75	5	} 50/1	1 ml
dNTP	15	15	2.5		.5
primer 1	7.5	7.5	2	} 10/1	2
primer 2	7.5	7.5	1.5		1.5
Mg	—	15.0	1	} 10/1	1
template	3.0	3.0	.5		.5
H ₂ O	642.0	627.0	0		0

distributed 50 µl / tube added enzyme.

0	1	2	(tube #)	16	17
1	3	4		18	19
2	5	6		20	21
3	7	8		22	23
4	9	10		24	25
5	11	12		26	27
6	13	14		28	29

0.01 15 30
 Deepvent mix

1x

Deepvent buffer

Klenow buffer

To Page No. _____

Used & Understood by me,

Date

Invent d by


Date

Recorded by

11/22/94

K. Schramm

From Page No. _____

Tag titration

← D.V. buffer

← K.T. buffer

Result: more product with increasing amount of Tag expected.

- K.T. B / 1U better than D.V. / 1U

- 1:0.01 better than 1U Tag alone.

- K.T. B more product than D.V. buffer

- But bit of misreading - adjust the cycling condition

T Page 1

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

K. Subraman

11/22/94