

The over Heparin 40 mL column.

Proj ct No. _____

B ok No. _____

Page No. _____

Spin $(\text{NH}_4)_2\text{SO}_4$ sol'n - @ 18,000 x g 40 minutes -
in GSA rotor -

Save supernatant -

Save pellets -

Store one pellet in -20°C - process the other
pellet 2

pellet 1 slightly greater than half - $\sim 3/5$

pellet 2 slightly less than half - $\sim 2/5$

Resuspend pellet in 20 mL of Buffer 1 -

Buffer 1

5 mM Tris pH 7.5

3.1 glycerol

20 mM KCl

5 mM Bme

.1 mM PMSF

dialyze - against Buffer 1 for ~ 8 hrs -

Exchange buffer 4 times -

Heparin column - use prepacked Heparin from A.G. -
 ~ 40 mL column - bump w/ Buffer + KCl -
wash w/ H_2O -

Previously A.G. used 3 mL Heparin a
5 gram crack

Direct scale up = $\frac{3}{5} = \frac{4}{5} = 30$ mL Heparin
 ~ 50 g

Equilibrate w/ Buffer 1 \rightarrow (Note: made 20 mM KCl -)

To Page No. _____

Used & Understood by me,

Date

Invented by

Date

May Forgo

4/5/95

Recorded by

08/30/95

From Page No. _____

Conductivity of Load - 2.8 mS - after ~ 8 hrs of dialysis

Notice a small precipitate matter in dialysis tube -
Spin down in SS-34 - 18,000 x g - 10 minutes -
same pellet - small + white -

① Load - 21 mL of sample - 75 ~~mm~~ mL/min - collect FT -

② Wash - 2 V_t of Buffer 1 - collect ^{7.5} 8 mL fractions
1 mL/min

③ Gradient - Buffer 1 to Buffer 2 - 25 mM Tris pH 7.5
8% glycerol
5 mM BME
1 mM PMSF
2 M KCl

10 V_t - 400 mL gradient - linear - 1 mL/min -
collect 7.5 mL fractions -

④ Wash w/ 2 V_t Buffer 2 1 mL/min - 7.5 mL fractions -

Let column run O/N -

Note Next time gradient should be much shallower - 1 M KCl -

To Page N

Witnessed & Understood by me,

Date

Invented by

Date

- Mary Long

4/5/95

Recorded by

8/31/95