

## Annealing:

	Kinased Oligo 2680	Kinased oligo 2723	Control
Ann. Buffer	1 $\mu$ l	1 $\mu$ l	1 $\mu$ l
U-SS DNA (gene cleaned)	1.5 $\mu$ l	1.5 $\mu$ l	1.5 $\mu$ l
H <sub>2</sub> O	6.5 $\mu$ l	6.5 $\mu$ l	7.5 $\mu$ l

Anneal at 70°C for 2 min  $\rightarrow$  let cool to  $\sim$ 35°C over 30-40 min.

## Synthesis:

Add to tubes.

1  $\mu$ l 10x Syn. Buffer  
0.3  $\mu$ l T7 DNA Pol (USB)  
1.0  $\mu$ l T4 ligase (LT1)

S' on ice, S' at room temp — 45 min at 37°C.

Add 90  $\mu$ l TE (10mM Tris pH8.0, 10mM EDTA).

Transfect DH5 $\alpha$  F' IR.

Transform DH5 $\alpha$  F' IR 100  $\mu$ l with 1.5  $\mu$ l ligation/synthesis mix.

30 min on ice

42°C for 40 Sec.

$\rightarrow$  Split 10% and 90% for 2680 and 2723

Keep 100% for the control.

Add 6  $\mu$ l lawn cells + 3 ml Soft agar

- let the agar solidify

$\rightarrow$  Incubate plates at 37°C.

To Page No.

Issued &amp; Understood by me,

Date.

Invented by

Date

Recorded by