

Kinasing

oligo	2679	2680	2723
-------	------	------	------

✓ μ l	2.5 μ l	2.5 μ l	2.5 μ l
-----------	-------------	-------------	-------------

✓ 1 M Tris pH 8.0	3 μ l	→	
✓ 0.15 M MgCl ₂	2 μ l	→	
✓ 0.5 M DTT	0.3 μ l	→	
✓ 10 mM ATP	1.3 μ l	→	
✓ H ₂ O	20 μ l	→	
✓ T4 Kinase (100/ μ l)	0.5 μ l	→	(T4 Kinase from J. Wang).

30 min at 37°C. Heat killed at 65°C for 10 min.

Add 30 μ l TE - Freeze in dry ice → Save at -70°C until further use.

∴ Final conc. of oligo \approx 250 pmol/60 μ l \approx 4 pmol/ μ l.

To Page No. _____

essed & Understood by me,

Date ,

Invented by

Date

Recorded by

Deb H. Bhatia