

Project No. _____

Book No. _____

TITLE

SOS gel on 80ml Heparin

From Page No. _____

I' DHIOBRK2
P190mv
157 200mv
152

I' TFI P171

200

Am 504 resur 362 μ A

1.5

5200.50 μ A

105

Blue

62

Heparin frns

105

150

106

150

107

150

108

150

109

150

110

150

111

150

112

150

113

150


114

150

115

150

TFI 10 + #TF31010A-502

H₂O

150

100 100 - - 238 150 →
K_f = 300

15% TCA

300 → - - 300 →

P 72, 9 and P 50, 7)

30' ice

1 x sample buffer

35 50 - - 400 → 35

Witnessed & Understood by me,

Deanna Polys

Date

5/1/95

Invented by

Record d by

Dat

4-28-95

To Page 1

PAGE 5 OF NOTEBOOK WAS BLANK

storage buffer

From Page No. _____

buffer F, G as per rTag 91342. PRP

make F, 4L first as follows.

make vol up to $\frac{3200}{2700}$ ml (i.e. 80% of Vf of

remove	160 ml	and add
	20 ml	Tween 20 (Pierce)
	20	NP40 (Pierce)
ml	$\frac{20}{200}$	

= buffer G

Take the remaining buffer up to Vf = 3800
for 1X buffer Fpool frn 105-114 of Heperan (P. 1-4)
= 90 ml vol (actually measured 87 ml)Dialyze against 2 buffer F, 5 hr
2 ml C O/N

Recovered 33 ml after Dialysis -

4-30

combine with 33 ml buffer G = Vf 66.

labeled:

TFI DNA pol in
storage buffer 4-30-95

stock at -20

4.33u/l
all PB

To Page N

Witnessed & Understood by me,

Deena R. Rupp

Date

5/1/95

Initiated by

Recorded by

Date

4-29-95

4-20-95

Project No. _____

Book No. _____

TITLE

Storage buffer

6

From Page No. _____

buffer F, G as per rTag 91342. PRP

make F, 4L first as follows.

make Vol up to $\frac{3200}{2700}$ ml (i.e. $80\frac{1}{2}\%$ of 4L of
remove 160 ml and add
20 ml Tween 20 (Pierce)
20 ml NP40 (Pierce)
ml $\frac{200}{200}$
= buffer G

Take the remaining buffer up to Vf = 3800
for 1X buffer F

pool frn 105-114 of Heparin (P. 1-4)
= 90 ml vol (actually measured 87 ml)

Dialyze against 2 buffer F, 5 hr
2 ml L O/N

Recovered 33 ml after Dialysis - 4-30
combine with 33 ml buffer G = Vf 66.

labeled: TFI DNA pol in
storage buffer 4-30-95 stock out - 20

4.33u/
ACP 8

Witnessed & Understood by me,

Deena R. Rump

Date

5/1/95

Invented by

Recorded by

Date

4-29-95
4-30-95

T Pag N

Unit assays for TH1

ge N	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27
storage buffer																											
1-30.95																											
ten*																											
.50	2				2				2				2														
.00		2				2				2				2													
.000			2				2				2				2												
.0000				2				2				2				2											
30209B																											
5 (ACP181)																											
50																2											
.00																	1										
.00																		1									
.000																			2								
31010A-502																											
4-25.51 PFD, 9																											
75 gel																											
70																					2						
.00																						2					
.00																							2				
400																								2			
4-25.95																											
SPS gel																											
70																									2		
.00																										2	
.00																											2
.000																											2
Exu mix 48 μ l																	48										
0, 9)																											
unit assay																	48										
mix																											

74°C, 10', → 10 μ l 0.5M EDTA → spot 40 μ l on GFCtubes # 13-16 are same serial dilution as tubes # 1-4
all 3 independent serial dilutions in Tag dip buffer

To Page No. _____

ed & Understood by me,

Date

Invented by

Date

Serecia Polak

5/1/95

Recorded by

5-2-95

Endonuclease Qc

Pr j ct N _____

Book N _____

Exhibit 84

Appl. No. 09/558,421

9

eN follow + Tag QC: 10342 • QCP

	22 Rxns	SS DNA Rxn mix	ds DNA Rxn mix	
CP buffer		110 μ l	110 μ l	✓
M Tris HCl pH 8.4				
10 mM KCl				
1 mM MgCl ₂		110	110	✓
74 (+) ssDNA 0.2 μ g/ μ l		110 μ l		✓
174 RF 0.33 μ g/ μ l			66.7	✓
Qc H ₂ O		666 990	703.3 990	✓

7.30.55
3 μ /l

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17
.48	.92	1.39	1.8	2.3												

micr
0A-502
1/1 (P8)

.57	1.14	1.7	2.3	2.9												

rudd xpt
I 4-13-95
K

10,000
1,000
100
10
0 dil

4.5 4 3.6 3.15 2.7 4.4 3.9 3.3 2.7 2.1 4 4 4 4 5 - ✓
5

17 is ⁴⁵ μ l SS DNA Rxn mix and #18-34 is ds DNA mix ⁴⁵ μ l Rxn mix

To Page No. _____

ed & Understood by me,

blanco

Date

5/15/95

Inv nted by

R c r d b y

Date

5-3-95

Project No. _____


Book No. _____

TITLE _____

From Page No. _____

units:
 LTI TFI 2 4 6 7 8 9 10
 Epicenter 1 2 3 4 5 6 7 8 9 10
 crude ext 1 2 3 4 5 6 7 8 9 10
 Blank 1 2 3 4 5 6 7 8 9 10
 buffer blank

(relaxed) RFII
 ssDNA →



} ssDNA
 end

} ssDNA
 end

note for RF II substrate with positive control (FrT)
 there is a slight increase in relaxed circles (ie nick
 and immediate conversion to small fragments.
 no nicking seen for either ssDNA or dsDNA
 substrates by any TFI or Epicenters.

To Page N

Witnessed & Understood by me,

[Signature]

Date

5/15/95

Invented by

[Signature]

Date

5-3-95

Recorded by

Solutions of Tag for PCR functional assay


Project N _____
Book N _____

Exhibit 85
Appl. No. 09/558,421

11

Tube No.	#	1	2	3	4	5	6	7	8	9	10
CF u/ml		.0625	.125	.188	.25	.5					
Lot #											
EM7414		5	5	5	5	5	3	3	3	3	4
Solution buffer		395	195	128.3	95	45	237	117	77	57	27
Vf =		400	200	133.3	100	50	240	120	80	60	36

R401	#	6	7	8	9	10
------	---	---	---	---	---	----



To Page No. _____

ed & Understood by me,

Polamp

Date

5/15/95

Invented by

Record d by

Date

32P "33 correct" (as done P146, 9 and 136, 9) and at 5' end primer of so

(P-2131) 3' end reaction can be Cf = 200 nM

33 mer correct (P13P, 9)	20 pmol/μl (20 μM)	7.5	15 μl	✓	✓	~ 300 pmol 33 mer
32P & ATP 10 mCi/μl		5	30 μl	✓		1.1 x 10 ⁶ cpm/μl
5-12-95		6	12 μl	✓		~ 90
5 x Kinase buffer		1	3 μl	✓		
PNK 1/2						
H ₂ O		30	60 μl			
					37°C, 30'	
					55°C, 5'	

Plan for fidelity assay for pol ± 3' end

(33 correct P13P, 9)

(-) dCTP (+) dATG-TP

CCAGTG A A T T G A G T G G T A
 C T T A A G E T C G A G C C A T G G G C C C C

↑
 same 5' end as
 23 mer on mp19+

↑
 only 3' into
 for rescue will
 have to do
 n-1 m
 3' into downstream
 from primer
 3' end


↑
 most run through
 stopped here for
 quantitation
 at 49-51

Unit array on 1he
(in Tag unit array map)

Project N. _____

Block No. _____

13



Understood by m ,

Signature

Date

5/15/95

Invnted by

Signature

Date

5-5-95

To Page No. _____

5' exp in pools: TFI, Tne

From Page No. _____

see P 136, 9 137, 9 for procedure: 200 nm primer (so its like a real PCR)

	(1)	(2)	(3)	(4)	(5)	(6)	(7)	(8)	(9)	(10)	(11)	cocktail
10x PCR buffer	5	5	5	5	5	5	5	5	5	5	5	60 μ l ✓
50 mM MgCl ₂	1.5 μ l	1.5 μ l	1.5 μ l	1.5 μ l	1.5 μ l	1.5 μ l	1.5 μ l	1.5 μ l	1.5 μ l	1.5 μ l	1.5 μ l	18 μ l ✓ (1.5 mM)
3'P 33mer correct"	2	2	2	2	2	2	2	2	2	2	2	24 μ l ✓
P12 (5 μ M)												102 μ l - use 8.5 μ l/
H ₂ O												
TFI FrI 4-1355	10											for # (11)
TFI FrI'/PEI 4-1355	10											✓ 10 μ l of 5x buffer
(2.45 μ /l)												Nin (2-14.55 on
TFI in storage buffer	10											has dNTPs and 1
4.33 μ /l (P8)												
TFI epimenter 6ATF31010A-502	12.4											20 μ l TFI/vant
(3.5 μ /l P.7)												0.745 μ l from 1
rTag 5 μ /l EKBT1	10											✓ 2 μ l 3'P 33 a
The 3-245 1.72 μ /l												✓ 18 μ l H ₂ O
(P.13) dilute to 0.5 μ /l	2/0											50 μ l ✓
The 5-7-95 36 μ /l												2/10 -
(according to Lig F and see P13												
where I got 34 μ /l)												
dilute to 0.5 μ /l												
1% Tween 20/PP4												
H ₂ O												
Tag dil buffer												
74°C												
remove 10 μ l to 5 μ l vial seq stop												
at 0.33, 1, 2 hr												
run on 8% PAGE												

see analysis of TFI/vant exp rate and TFI loss of full length 33mer on P 46.

To Page N

Witnessed & Understood by me,

Date

5/15/95

Invented by

Date

5-10-95

Recorded by

PAGE 15 OF NOTEBOOK WAS BLANK

From Page No. _____

10 PCR buffer

Mix A
 110 μ l

✓

CP = 1X

acc
 to
 QC1

50 mM MgCl₂

110 μ l

CP = 5 mM

103'

3 (JS) substrate

44 μ l

(1 pmol / reaction)

0.5 pmol / λ
 lot EFE73

autoclaved, filtered H₂O

726 μ l

✓

VF = 990 μ l

winds 0 2 4 6 8 10

tube # 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15

4-30-95

TF1 (P. 8)
 4.33 μ l / λ

0.46 0.92 1.39 1.85 2.3

TF1 epicenter

lot TF31N10A502

3.5 μ l (see P. 8)

0.57 1.14 1.71 2.28 2.86

TF1 Fr I

4-13-95

dil none 1/100

1/10 1/10

1/100 none

2
 2
 2

H₂O 3 4.54 3.63 2.75 4.43 93.3 2.7 2.1

✓

Tox storage buffer 2

Mix A 45 μ l

Tox dil buffer

✓

✓

✓

VF = 50 μ l

5 3 →

74°C, 60'

37°C, 60'

put tubes on ice

To Page N

Witnessed & Understood by me,

J. Polansky

Date

5/15/95

Invented by

Record d by

Date

5-12-95

PET

Project N _____


Block No. _____

17

1 N — cut plate in half → get 20 x 10 cm plates

use 0.9 cm wide lanes
origin at 2 cm

score lanes
every 0.9 cm



spot 5 μ l, dry, spot another 5 μ l, dry
resolve in 2 N HCl

bring solvent front to top of plate

Dry by heat lamp, not more than 7 min

count bottom — 3.2 cm (= ori)

3.2 cm — 7.6

7.6 — top

for each

3 x 12 = 36 tubes

add 3.5 ml fluor

count 3H

To Page No. _____

Read & Understood by me,

Date

3/15/95

Invented by

Recorded by

Date

5-12-95

201amp

Project No. _____

Book No. _____

TITLE no released

From Page No. _____

		CPM	
TFI 4-30-85	0	4802.00	
	1	673.00	
	2	92.00	
	3	4782.00	BKGD
	4	709.00	
	5	123.00	31
	6	4935.00	0
	7	702.00	
	8	114.00	22
	9	4670.00	0
TFI 4-30-85	10	743.00	
	11	129.00	37
	12	4732.00	0
	13	669.00	
	14	123.00	31
	15	4788.00	0
	16	767.00	
	17	123.00	31
	18	4661.00	0
	19	613.00	
TFI 4-30-85	20	84.00	
	21	4677.00	
	22	636.00	
	23	145.00	53
	24	4185.00	0
	25	767.00	
	26	140.00	48
	27	4586.00	0
	28	774.00	
	29	112.00	20
TFI 4-30-85	30	4136.00	0
	31	887.00	
	32	111.00	19
	33	4202.00	0
	34	793.00	
	35	97.00	1
	36	1404.00	0
	37	396.00	
	38	1801.00	1709
	39	70.00	50%
TFI 4-30-85	40	133.00	
	41	2536.00	2461
	42	27.00	95%
	43	208.00	
	44	2466.00	2374
	45		91%

Conclude no ds 3' end activity is detected
in either LTI or Epicenter primed TFI pol

To Pag Nc

Witnessed & Understood by me,

Date

Invented by

Dat

5/15/95

Record d by

5-12-95

ge No. _____

Epicenter TFI
storage buffer

LTI rtag
storage buffer

is HCl

50 mM pH 7.6
 (pH at Room Temp
 based on Tale call to
 Epicenter Tech line)

20 mM pH 8

need to increase
 molarity from
 20 mM to 50 mM
 and pH from
 8 down to 7

all

100 mM

0

need to add
 100 mM NaCl to
 LTI SB

glycerol

50%

"

LTI Buffer

DTT

1 mM

"

as

EDTA

0.1 mM

"

epicenter

rem 20/MP40

0.5% each

"

Experiment

1 M Tris HCl

pH (room Temp)

10 ml buffer G

0

8.03

(expect 8.0)

5 λ

7.99

20 λ

7.92

20 λ

7.81

20 λ

7.77

20 λ

7.72

20 λ

7.67

add 125 μ l
 of 1 M Tris HCl
 to 10 ml TFI
 4-30-95

$C_f = 20 \text{ mM Tris (in SB)}$
 $+ 12.5 \text{ mM Tris HCl added}$
 32.5 mM

will follow this exact procedure
 for 10 ml of TFI 4-30-95
 (P.8 4.33 μ l \Rightarrow new $C_f = \frac{10.0}{4.33} = 2.31$

add 182 μ l 1 M Tris pH 7.5 $\Rightarrow C_f = 50 \text{ mM}$ pH = 7.60 pH

$V_f = 182 \mu\text{l} + 125 \mu\text{l} = 307 \mu\text{l}$

mix end over end

plus

307 μ l glycerol

30 to 60 min

V_{total}

= 10.614 ml

plus

62 mg

$C_f = 100 \text{ mM}$ To Page No. _____

Val
 (new)

d & Understood by me,

Date

Inv nt d by

Date

Solans

5/15/95

R cord d by

5-12-95

Stability study of Δ in (m)
TFI / vent mix (5-16-95)

From Page No. _____

Follow p. 84, 9

[A]

5 x Cheng (no ATP)
H₂O
activation DNA 3.7 mg/ml
dATG-C-TP 10 mM ea
32P dATP 10 mCi/ml
5-19-95 w/ date 3000 Ci/mmol

200 μ l
637.9
135.1

5 μ l
2 μ l

✓

✓

✓

✓

Cheng at 1X

glycerol

Tris pH 9.0

K₂CO₃

DMSO

mgOAc

pH

2.0

85

2.0

1.05

0.5 mg/l

CF 50 μ mV_f = 570 μ l

for 10 Rxns

① ② ③ ④ ⑤ ⑥ ⑦ ⑧ ⑨

[A]

981

TFI LTISB

1.891

+ Vent (5-16-95)

TFI Epicenter SB

+ Vent (5-16-95)

2

TFI Epicenter Enzyme + Vent

(5-16-95)

1.891

V_f = 100 μ l

68 °C

remove 10 μ l to 5 μ l 0.2M EDT

(spot 10 μ l on 6 FC) and remove 5 μ l to 5 μ l Kill solution
with cold dAMP (PP4, 9) and spot 2 μ l on PEI (PP4,
at 0, 10, 20, 40 min

Repeated on P. 40

Conclusion: There is no real good way to do this upper
because Vent is present at very low level compared to
so incorporation is saturated at levels where turnover
is barely detected by PEI method. Will have to settle
dDNA 3' end QC assay for units see p 26

To Page 1

Witnessed & Understood by m,

Date

Inv nt d by

Date

D. S. S. S. S.

6/9/95

Rec rd d by

5-11-95

Pr j ct N _____

Bo k N _____

21

ge No.

247.00

5736.00

4880.00

5046.00

2170.00

3243.00

3897.00

5694.00

2774.00

3405.00

4090.00

5613.00

5451.00

5498.00

5034.00

5820.00

5169.00

5604.00

5612.00

580.00

5193.00

5755.00

480.00

542.00

053.00

113.00

817.00

182.00

1925.00

3149.00

2983.00

3230.00

1091.00

1772.00

3000.00

3722.00

1 68601.00

2 118176.00

3 155582.00

4 95344.00

5 125044.00

6 138326.00

7 93269.00

8 119376.00

9 167655.00

10 93777.00

11 116666.00

12 131003.00

13 109619.00

14 126936.00

15 143456.00

16 90599.00

17 103792.00

18 162204.00

19 96493.00

20 124924.00

21 162506.00

22 91166.00

23 125191.00

24 167630.00

25 84292.00

26 127063.00

27 153977.00

28 2135.00

29 975.00

30 110991.00

31 110539.00

pmol

1487

2562

5372

BKGD

BKGD

ave 110765 69.2 cm/yr

To Page No. _____

ed & Understo d by me,

Date

Inv nted by

Date

R corded by

B. A. Polay

5/1/85

5-16-85

Fr m Page N _____

* buffer "S"

1M Tris HCl 7.5 pH
 0.5M EDTA
 glycerol
 3 me
 3M KCl

20 ml ✓
 0.2 ml ✓
 80 ml ✓
 0.357 ✓
16.7 ml ✓
 1 L

cf
 20 mM
 0.1 mM
 8%
 5 mM
50 mM

provided 6 ml S200

wash equilibrated at 1 col vol/hr
 (= 0.1 ml/min) for 2 hr

load 120 ml of Tne 36 u/l 5-7-95
 (= 2% vol/vol) (4320 units total)
 (load before S200)

by gravity.

elute at 1 col vol/hr
~~into 95 fr~~

collect 50 µl frns (30 sec/fr)
 span ~0.5(??)
 A

* note no detergent in buffers
 maybe this is why activity died after
 few weeks at 4°C (see P53)

To Pag 1

With ssed & Understood by me,

Date

Inv nted by

Dat

Polans

6

9/9/95

Rec rd d by

5-17-95

Project No. _____

Book No. _____

TITLE Assay 5200 fractions

24

From Page No. _____

fr #
Silicate 100

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16
10 11 16 19 22 25 28 30 35 40 45 50 55 60 65 70

2 μ l

Top Rxn m
pH 2.0, 5

48

10, 74

Results from
Top of P. 25 μ /ml

fr		1, mol	μ /ml
10	1	9978.00	1.99
13	2	9430.00	
16	3	16574.00	
19	4	12045.00	
22	5	10155.00	
25	6	7256.00	2.2
28	7	6430.00	1.9
30	8	6091.00	1.8
35	9	3680.00	1.1
40	10	2746.00	0.72
45	11	2002.00	0.65
50	12	1476.00	0.44
55	13	1118.00	0.35
60	14	694.00	0.2
65	15	970.00	
70	16	603.00	
172	17	120535.00	75 cpm/pmol

fr			μ /ml
3	1	3216.00	0.97
4	2	3990.00	1.2
5	3	5279.00	1.6
6	4	5148.00	1.55
7	5	6128.00	1.8
8	6	8950.00	2.7
9	7	11386.00	3.4
10	8	16350.00	4.9
11	9	14464.00	4.5
12	10	19127.00	5.8
13	11	23242.00	7.0
14	12	24609.00	7.4
15	13	25276.00	7.6
16	14	24319.00	7.5
17	15	23374.00	7.4
18	16	16929.00	5.1
19	17	14764.00	4.4
20	18	15075.00	4.5
21	19	13028.00	3.9
22	20	10781.00	3.2

Witnessed & Understood by me,

Polamp

Date

6/9/95

Inv nted by

[Signature]
R cord d by

Dat

5-12-15

To Pag 1

N	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	
1/100	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	
Results on P 24																										
unit assay on fr 13-17 : pool 20 µl of each																										
3-17	1	2	3	4	5	6	7	8	9	10																
100	2																									
200	2																									
400	2																									
800	2																									
1600	2																									
3200	2																									
7-95 36 u/l											134625	11	13	1	183.00											
1000											22 Remains	12	14	1	216.00											
2000											22 Remains	13	15	1	115777.77											
4000											22 Remains															
8000											22 Remains															
<div style="display: flex; justify-content: space-between;"> <div> <p>pool u/ml</p> <p>338 5.1</p> <p>7.8</p> <p>6.9</p> <p>2.5</p> <p>2.7</p> </div> <div> <p>ave 7.7 u/l</p> </div> <div> <p>well normalize 5200</p> <p>pool 13-17 against 36 u/l</p> <p>value for Tite 7-95</p> <p>(see P 13) 00</p> <p>pool 13-17 (7.7 u/l / 36 u/l)</p> <p>170.7</p> <p>= 3.9 u/l</p> <p>(3.9 u/l) (7 frms) (15 µl/frm)</p> <p>13.65 units recovered.</p> <p>70.7 ave</p> <p>adjusted only 36 (see P 13)</p> </div> </div>																										
3	1	16254.83																								
4	1	12439.50																								
5	1	5518.00																								
6	1	2909.00																								
7	1	1741.00																								
8	1	866.00																								
9	1	17122.03	355	53.2																						
10	1	10433.3	45	65																						
11	1	5683.00	70	70																						
12	1	3137.00	77	77																						

To Page No. _____

d & Understood by m ,

Date

Invented by

Date

To Pag No. _____

R corded by

5-18-95

5-18-95

ig N	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30
1/100	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30		

Results on P 24

9-95 unit assay on fr 13-17 : pool 20 µl of each

13-17	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30
1/100	2																													
1/200		2																												
1/400			2																											
1/800				2																										
1/1600					2																									
1/3200						2																								
5-7-95 3/4 µl																														
1/1000																														
1/2000																														
1/4000																														
1/8000																														

136025 11 13 1 183.00

12 14 1 216.00

22 Remains 13 15 1 115777.77 72.9

10, 74°C

will normalize 5200
 pool 13-17 against 36%
 value for Tue 5-7-95
 (see P 13) 00

$$\text{pool 13-17} \left(\frac{7.7}{170.7} \right) \left(\frac{1}{36\%} \right) = 3.9 \mu/l$$

(3.9 µ/l) (7 frms) (15 µl/frm)
 1365 units recovered.

70.7 ave
 expected only 36 (see P 13)

	3	1	16254.83	pool	u/ml
	3	1	16254.83	338	5.1
	4	1	12439.50		7.8
	5	1	5518.00		6.9
	6	1	2909.00		2.5
	7	1	1741.00		8.7
	8	1	866.00		
	9	1	17122.00	355	53.2
	10	1	10433.3	45	6.5
	11	1	5683.00	70	70
	12	1	3137.00	77	77

sed & Underst od by me,

Polamp

Date

6/9/95

Invent d by

Record d by

Date

5-18-95
 5-19-95

To Pag No.

Project No.

Book No.

TITLE

3' exp (QC) assay
for TFI / vent

26

From Page No.

(see P20 where turnover on PET didn't work)
follow assay on P 11-17actual units of Vent
added based on 0.09 u vent / 0.945
TA / u

Mix A P. 16

0 1 2 3 4 5 6
4545TFI / vent LIESB
(5-16-95)

1/10

1/100

1/1000

1/10000

1/100000

2

2

2

2

~ 0.2

.02

.002

.0002

.00002

.000002

~~TFI / vent Epinephrine SB
(5-16-95)~~~~1/10~~~~1/100~~~~1/1000~~~~1/10000~~~~1/100000~~~~TFI / vent with
epinephrine TFI (5-16-95)~~~~1/10~~~~1/100~~~~1/1000~~~~1/10000~~~~1/100000~~H₂O

3

BK 6-D

68 °C, 60 min

88.00 21.31

BK 6-D start 2:35 pm

no dil

1

1684.00

4.87

159.6

1/10

2

1323.00

5.49

123.5

1/100

3

448.00

9.44

36.0

1/1000

4

143.00

16.71

5.1

1/10000

5

96.00

20.39

8

1/100000

6

101.00

19.90

13

Result: assay is not very h
in assay range.
with assay at 1/10 and 1/100
in triplicate dilution

Witnessed & Understood by me,

D. Olap

Date

6/9/95

Invented by

Recorded by

Date

5-22-95

T Page 1

Result
5 λ enzyme dil

Project N

Exhibit 95

Appl. No. 09/558,421

Book N

29

TIME			Total Sum		%	
					lost from	
1 front	324.00					
2 middle	726.00	1.00				
3 ori	1780.00	1.00	2830			
4	348.00	1.00				
5	960.00	1.00	3117			
6	1809.00	1.00				
7	458.00	1.00	2991		(43%)	
8	830.00	1.00				
9	1764 1703.00	1.00	2817			
10	483.00	1.00				
11	608.00	1.00	2690			
12	1726.00	1.00				
13	477.00	1.00				
14	515.00	1.00				
15	1698.00	1.00				
16	288.00	1.00				
17	856.00	1.00				
18	168 1670.00	1.00				
19	729.00	1.00				
20	627.00	1.00				
21	1452.00	1.00				
22	632.00	1.00				
23	511.00	1.00				
24	1199.00	1.00				
25	471.00	1.00				
26	686.00	1.00				
27	135 1410.00	1.00				
28	374.00	1.00				
29	664.00	1.00				
30	1398.00	1.00				
31	200.00	1.00				
32	786.00	1.00				
33	1641.00	1.00				
34	128.00	1.00				
35	947.00	1.00				
36	1573 1682.00	1.00				
37	97.00	1.00				
38	477.00	1.00				
39	3007.00	1.00	3581			
40	66.00	1.00				
41	591.00	1.00				
42	3146.00	1.00	3803			
43	60.00	1.00				
44	480.00	1.00				
45	3122 3214.00	1.00	3754			

Result. The sum of all counts when enzyme is present (ave 2830) is less than for the no enzyme blank (ave 3713)

so 22% of counts unaccounted for!

in each case ~ 1500 CPM is lost from origin (ie ^3H removed by expo) but ~ 1/2 of that appears in front and middle. Also,

why does only appear in middle it looks like ^3H gap partly missing in middle.

since ^3H runs in front, it may be quenched by growth contaminants that smear PET at the front.

% label removed from dsDNA substrate

	2 λ enzyme mix	5 λ enzyme mix
LTISB	19%	43%
Epicate SB	15	46
Epicate TFI	26	57
Void pol	21	50

To Page No.

d & Understood by me, Bokamp		Date 6/9/95	Invented by 	Date 5-24-95
			Recorded by	

g N

	SAM	CPM1	pmol	u/ul	
			(-D/GD)		
			127	4.77	(was 4.33 on P8)
Tf1	1	4165.00		5.14	
	2	2395.00		6.22	
	3	1575.00		6.9	
	4	1018.00		9.70	
	5	4108.00		6.7	
	6	3019.00		7.5	
	7	2232.00		7.7	
	8	1205.00		4.3	
LTI SB	9	3788.00		5.0	
	10	2354.00		5.47	
	11	1425.00		7.62	
	12	1090.00		4.7	
	13	4135.00		5.0	
	14	2353.00		6.78	
	15	1688.00		7.07	
	16	1095.00		4.00	
Epimant SB	17	3543.00		5.06	
	18	2336.00		5.51	
	19	1431.00		7.13	
	20	1040.00		3.81	
	21	3388.00		4.64	
	22	2191.00		5.12	
	23	1354.00		5.31	
	24	857.00			
Epimant Tf1	25	3746.00			
	26	2188.00			
	27	1548.00			
	28	866.00			
	29	4053.00			
	30	2474.00	71	5.3	
	31	1456.00			
	32	905.00			
B/KGD 2x Rx mix	33	322.00			
	34	72445.00	45.3	cpm/pmol	

2
2
2 -

will repeat this with 5 duplicates
of the 1/350 ul dil for each
on P. 47

To Page No. _____

d & Understood by me, <i>Polamp</i>	Date 6/9/95	Invented by <i>[Signature]</i>	Date 5-24-95
	R corded by		

Project

Book No.

TITLE

units: 1.1X field test ("old")
new, B and M 2X mix

34

From Page No. 10

Tag 1-31-95

54/pt

1/125

1/250

1/500

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24 25

2 2 2 2 2

1.1X new (may 8, 1995)

no dil

1/2

1/4

2 2 2 2 2

1.1X old

no dil

1/2

1/4

2 2 2 2 2

BM 2X

1/2

1/4

1/8

2 2 2 2 -

Tag 1-31-95

P120, 9)

48/pt

1-500

74°C, 10'

see circled #0 P35

for 25 Tag 1/25 at 2 pt undisturbed of sample
new 6506 + 7080
Tag 10692 + 12403 $(.04 \mu\text{A kg by def}) = .025$ 6729 + 6660
10692 + 12403 025

Witnessed & Understood by me,

Polansky

Date

6/9/95

Invented by

Rec'd by

Date

5-25-95

To Page 1

ig No.

		pmol	u/μl			
1/25	1	10692.00	412	7.7		
1/25	2	5333.00		7.7		
1/25	3	3112.00		8.0		
	4	12403.00		8.9		
	5	6387.00		9.2		
	6	3853.00		10.0		
1/25	7	6505.00	250	0.38	1.78	1.71
1/25	8	4314.00		0.50	2.50	2.27
1/25	9	2381.00		0.47	2.39	2.27
	10	7000.00		0.40	2.0	2.27
	11	4401.00		0.51	2.5	2.27
	12	2364.00		0.47	2.4	2.1
	13	6729.00		0.39	1.9	1.77
	14	3962.00		0.46	2.3	1.73
	15	2247.00		0.41	2.2	
	16	6660.00		0.38	1.9	
	17	3659.00		0.38	1.9	
	18	1940.00		0.37	1.9	
1/2	19	3456.00		0.36	0.77	
1/4	20	1705.00		0.32	0.77	
1/2	21	1368.00		0.47	0.94	
1/2	22	3028.00	117	0.55	0.84	
1/4	23	2005.00		0.39	0.94	
1/2	24	900.00		0.26		
	25	333.00	BKGD			
	26	62215.00	38.9 cpm/pmol			

7.7 is ~ right since this "5 u/l" is normalized to amplitude that was ~ 8 u/l

Total units in 50 μl at 1.1X

* in 50 μl at 1X

average

* 1.76

2.22

* 1.75

1.88

0.86

red is average
total units in 50 μl
at 1X for the array of
undiluted μl of 1.1X mix

* mix is 1.1X
so at 1X mix is
(u/μl in 1.1X)

or use 1.76 units for "new"
or 1.75 units for old "field test"

as my first time points for these samples. (note: expect ~2
nd dilutions above indicate 2.22 u and 1.89 u for old and new).
see next assay on P52 which is 1 month time point
using 2 μl undiluted.

To Page No.

ed & Understood by me,

Polans

Date

6/9/95

Invented by

Record d by

Dat

5-25-95

Project No. _____

Book No. _____

TITLE

unit assay for stability
array of 1.1 x 10⁶ (p121, 9) (at
ave of 1.1 x 10⁶ of p. 122) zero time pre

36

From Page No. _____

1	8878.00			
2	9039.00			
3	9623.00			
4	8220.00			
5	8228.00			
6	8109.00			
7	8855.00			
8	8307.00			
9	8584.00			
10	6857.00			
11	7096.00			
12	6660.00			
13	9295.00			
14	8535.00			
15	8519.00			
16	6780.00			
17	5930.00			
18	5879.00			
19	8250.00			
20	8545.00			
21	9288.00			
22	8605.00			
23	7590.00			
24	7975.00			
25	7909.00			
26	7993.00			
27	7301.00			
28	2863.00			
29	3151.00			
30	3188.00			
31	7926.00			
32	7626.00			
33	7355.00			
34	8180.00			
35	8930.00			
36	9000.00			
37	2662.00			
38	2885.00			
39	2632.00			
40	8091.00			
41	7872.00			
42	7722.00			
43	7664.00			
44	7828.00			
45	8063.00			
46	10091.00			
47	9701.00			
48	10062.00			
49	10476.00			
50	10230.00			
51	101.00			
52	56411.00			
53	57488.00			

all can
result
on p. 1

← (repeat in
del. 10/15/5)

← (repeat)

data from
p. 37, 38

10112 ave ⇒ 42.6 pmol ⇒ 7.99 units / 1 in
(expect up since 5th 11 EKBTI
as unitized & amplifying of w8
(35.6 pmol / pmol)

To Page No. _____

Witness d & Underst od by m ,

Date

6/9/95

Inv nt d by

Record d by

Dat

5-26-95
5-26-95

Stability of 1:1X at room temp

Project N

BKN

Exhibit 98

Appl. No. 09/558,421

37

g No.

0 time on P. 154, 9
1 month P. 174, 9

3-13-95
4-11-95

assay same as P. 121, 9 for 4°C stability study
used same assay mix as P. 5-25-92

amp 12#	Reaction tube #	ul enzyme 12x100	1% Tween 20 NPH40	Top mit assay mix (P. 121, 9)
	1-3	2		48 µl
	4-6	1		
	7-9			
	10-12			
	13-15			
	16-18			
	19-21			
	22-24			
	25-27			
	28-30	↓	0.5λ	
	31-34	3.64		
→ dil 1/2.5	34-36	2		
	37-39	1		
	40-42			
	43-45			
Trilily	46-48 47		0.5λ	
	48-50 49			
5 µl 1-31-95	52-56	↓		
- dil xx	50-52			
	53 2λ of "old mix"			
	54 2λ of second mix			

74°C 10
kill with 10 µl
0.5M EDTA
spit 40 µl
on 6 FIC

10% TCA, 1% NaPP
5' ↓
3x 5' in 5% TC
1x 5' in 95% ETO

5 µl Top dil buffer + 8 µl sample #12
5λ Top + 620 λ dil buffer

To Page No.

sed & Understood by me,

Polamp

Date

6/9/95

Invent d by

Record d by

Date

5-30-95

Project No. _____

Book No. _____

TITLE _____

38

From						
1	1	3864.00				
	2	4350.00	.024	.037	2 months	1 month
	3	3113.00			65%	89
	4	155.00				
2	5	130.00	0.00090	0.033	2.7%	106
	6	133.00				
	7	4121.00				
3	8	4864.00	0.029	0.032	91%	106
	9	4225.00				
	10	3184.00				
4	11	3267.00	0.021	0.029	72%	83
	12	3240.00				
	13	3952.00				
5	14	3728.00	0.025	0.033	76%	91
	15	3817.00				
6	16	2585.00				
	17	2492.00	0.017	0.027	63%	78
	18	2795.00				
7	19	4009.00				
	20	4530.00	0.028	0.033	85%	84
	21	4334.00				
8	22	4199.00				
	23	4281.00	0.028	0.031	90%	89
	24	4326.00				
9	25	4248.00				
	26	3739.00	0.026	0.031	84%	90
	27	3922.00				
10	28	65.00				
	29	86.00	0.00053	0.022	2.4%	0
	30	93.00				
	31	4340.00				
11	32	3827.00	0.026	0.032	81%	88
	33	4038.00				
	34	5155.00				
12	35	5666.00	0.035	0.034	103%	—
	36	5301.00				
13	37	3927.00				
	38	3981.00	0.026	0.031	84%	97
	39	4352.00				
	40	4214.00				
14	41	3962.00	0.027	0.035	77%	83
	42	4287.00				
15	43	3992.00				
	44	3958.00	0.026	0.032	81%	106
	45	4278.00				
40c	46	2014.00				
	47	2016.00	0.013			
40c	48	1504.00				
	49	1509.00	0.0098			
70c	50	5941.00	(0.04 n/s) by definition			
	51	6122.00	617.7 same as P36 when corrected			
70c	52	6468.00	for 30ul on GFC instead of 40ul			
70c	53	49708.00	and 4 day decay of 72p			
new	54	40962.00	46819			

T Pag N

Witnessed & Understood by me,

Date

Invented by

Dat



6/9/85

Record d by

5-30-85

g N _____

sp act

$$\frac{46819 \text{ CPM}}{100 \text{ pmol dCTP} (4)} = 29.3 \text{ CPM/pmol(nt) DN}$$

$$\left(\frac{6177 \text{ avg Tag CPM}}{29.3 \text{ CPM/pmol}} \right) \cdot \left(\frac{60 \lambda}{30 \lambda} \right) = 421 \text{ pmol DNA synthesis}$$

1 unit 10 nmol / 30' at 74°C

$$\left(\frac{.421 \text{ nmol}}{10 \text{ mol}} \right) \left(\frac{30'}{10'} \right) = 7.89 \mu/\mu$$

$$\frac{1}{2} \text{ of sample} = \left(\frac{.04 \mu/\mu \text{ in Rm}}{7.89 \mu/\mu \text{ for Tag}} \right) \frac{\text{samples CPM}}{\text{Tag CPM}}$$

100 μM each dNTP

To Page No. _____

ed & Understood by m ,

Dolap

Date

6/9/95

Invented by

R corded by

Date

5-30-95

Project No. _____
Book No. _____

TITLE

Repeat of P. 20: stability of
TFI / Vent mixes

40

From Page No. _____

H₂O
5X Cheung buffer (no dNTPs)
activated DNA 3.7 mg/ml
dATG-C-TP 10 mM each
2 ³²P dATP 10 mCi/ml (ref 6-2-5T)
3000 Ci/mmol

[A]
446.4 ✓
140 µl ✓
94.6 µl ✓
3.5 µl ✓
1.5 µl

Cf=1

tube # 1-4 5-8 9-12 13-16 17-20 21-24 25-28 29-32
① ② ③ ④ ⑤ ⑥ ⑦ ⑧

[A] 98 µl →
2 µl TFI in epimix units
(its ~ 3.5 units/µl)
TFI LTISB 2
+ Vent (5-16-95) 10 µl →
TFI Epimix SB 2 →
+ Vent 5-16-95

2X 1X Cheung
- (no enzyme)

0.18 unit
Vent in

TFI Epimix enzyme
+ Vent (5-16-95)

2
15 →

mistake: this is
0.03 units / 100 µl
be 0.18 - ie Vin use
per 50 µl PCR + Vent

2 µl
0.45 u/l
Vent diluted
in 1X Cheung
buffer

2 µl of Vin
2 u/l of 2-
in 86.9 µl
1X Cheung

Vf ~ 100 µl

68°C
start with addition of enzyme to preheated mix.
remove 10 µl to 5 µl 0.2 M EDTA (spot 10 µl on GFI
and 5 µl to 5 µl killing solution with cold dAM.
(spot 2) on PET at 0, 5, 10, 15, 20 min.

resolve in 1M LiCl

T Pag N

Witnessed & Understood by me,

Polay

Date

6/9/95

Invent d by

[Signature]

Record d by

Dat

5-31-95

N		$\frac{\text{cpm} - \text{background cpm}}{\text{specific activity}}$		$\frac{\text{ave pmole turned over}}{\text{ave pmole incorporated}}$ 2 replicates of each time pt	
		pmoles	% turnover	should be constant	
AM	CPM1 turnover				
1 5	885.00	$\frac{885-543}{103.31} = 3.31$	0.042	result: 1) turnover began to reach a plateau by 10 min. we expected turnover to continue increasing over time after DNA synthesis stopped. TFI completes at Vent at nick	
2 10	1241.00	4.76			
3 15	1074.00	5.14	0.069		
4 20	1269.00	7.03	0.065		
5	984.00	4.27	0.062		
6	1332.00	7.64			
7	1678.00	11.0			
8	1590.00	10.1			
9	830.00	2.78		2) turnover by Vent alone in Cheng buffer is lower than expected. 25% turnover was observed in another experiment. 453 pmole in sam. 25% 1/25 sm-by vent. The high signal to noise level for the Vent samples makes it difficult to say what the turnover is. Turnover is higher when Vent is mixed w/ TFI. TFI creates mismatches that are targets for vent-exo. do: repeat w/ 1, 2, 3, 4, 5 min time points and more Vent enzyme in the Vent alone samples, use new PEI plates	
0	1213.00	4.49	0.021		
1	1195.00	6.31	0.066		
2	1460.00	8.88	0.055		
3	555.00	0.116	0.062		
4	1228.00	6.63			
5	1225.00	6.60			
6	1425.00	8.54			
7	764.00	2.14			
8	977.00	4.20			
	CPM3		0.039	1 replicate / time pt	
19	1212.00	6.48	0.028		
20	1453.00	8.81	0.043		
21	895.00	3.41	0.061		
22	772.00	2.22			
23	1009.00	4.51			
24	1365.00	7.96			
25	746.00	1.96	0.70		
26	438.00	1.02	0.47	do: repeat w/ 1, 2, 3, 4, 5 min time points and more Vent enzyme in the Vent alone samples, use new PEI plates	
27	757.00	2.07	0.14		
28	609.00	0.64			
29	412.00				
30	578.00	$\bar{x} 543$			
31	488.00	background			
32	693.00	no enzyme			

specific activity: cpm of 2ul spot of mix A 36438 $\bar{x} = 41,324$ cpm
2 replicates $\angle + 40210$

$$\frac{\left(\frac{100 \text{ul rxn}}{2 \text{ul spot}} \right) (41324 \text{ cpm})}{(5000 \text{ pmole}) (4 \text{ bases})} = 103.31 \frac{\text{cpm}}{\text{pmole (nt) DNA}}$$

50 μ M each dNTP in 100ul rxn
50 μ mole/L $\times 100 \times 10^{-6}$ L = 0.005 μ mole = 5 nmole = 5000 pmole

To Page No. _____

I & Understood by me,

Date

Inv nted by

Date

Polansky

4/9/95

Recorded by

Carolyn Combs

5-16-95

Project No. _____

42

Book No. _____ TITLE _____

From Page No. _____

$$\frac{\text{sample cpm}}{\text{specific activity}} \times \frac{100 \mu\text{L rxn}}{10 \mu\text{L spot}} \times \frac{15}{10} \text{ dilution}$$

Incorporation

pmoles

1	49596.00	103.3	$103.3 \times 15 = 7,201$
2	74066.00	10,754	
3	88521.00	12,853	
4	95661.00	13,889	
5	50395.00	7,317	
6	69543.00	10,097	
7	82738.00	12,013	
8	93515.00	13,578	
9	45114.00	4,550	
10	64768.00	9,404	
11	81250.00	11,797	
12	96711.00	14,042	
13	49095.00	7,128	
14	71796.00	10,424	
15	81335.00	11,809	
16	95798.00	13,909	
17	50290.00	7,302	
18	70938.00	10,230	
19	88754.00	12,887	
20	98147.00	14,250	
21	48881.00	7,097	
22	85245.00	12,377	
23	85694.00	12,442	
24	91420.00	13,274	
25	1932.00	281	
26	2581.00	375	
27	3000.00	436	
28	3120.00	453	
29	854.00	123	
30	777.00	113	
31	3183.00	26.6	no enzyme
32	32487.00	70.7	enzyme
33	9.00		
34	6.00		
2ul mix A	36438.00	-	\bar{x} 41,324 cpm
	46210.00	-	for calculation of specific activity

105 cpm/pmol

Synthesis ~~test~~ was almost complete by 10 min
By 20 min ~14 nmoles of the 20 nmoles
had been incorporated - hi in the 1st

To Page N

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

new dilution of rlog wt EKBT1
to 5 u/ml

Proj ct N .

Exhibit 100

Appl. No. 09/558,421

Bo k N .

43

ge N .

EKBT1
323 u/ml (see p91, 9)

157.2 μ l

Tox storage buffer
(Pierce detergent)

10 ml

$V_f = 10^{.157}$ ml

(cf = 5 u/l)

mix end over end 1 hour

storage buffer is from 12-7-94
with Pierce Detergents

To Pag No. .

ed & Understood by me,

Date

6/9/95

Invented by

Recorded by

Dat

5-31-95

Polansky

Project N

Book No. _____

TITLE

primer degradation using TFI//
(Can see P 17, the reaction # 11)

44

From Page No. _____

① ② ③ ④ ⑤ ⑥

5X Cheng
(no dNTPs)20 μ l \longrightarrow 32P "33 mer correct"
same as (P12 and 14)4 \longrightarrow 5 μ MTFI/Vent 5-16-95
LTI SB

10

TFI/Vent 5-16-95
Epicenter SB

10

(0.9 units va
in 10 μ l rxn)TFI/Vent 5-16-95
Epicenter TFI

10

TFI LTI 4.33 μ l (P.8)

10

Vent .09 μ l
lot #17 (opened 2-24-95)

10

LTI SB P.6 (same stock as in TFI P.8)

10

H₂O

66

V_f = 100 μ l* 2 μ l Vent.
diluted into
LTI Tag SB
2 μ l Vent
42.4 μ l Tag SB
44.4 μ l

68°C.

remove 10 μ l to 5 μ l aple rxn
stop rxn at 0, 5, 10, 20, 40, 60, 80, 100 madd enzyme on ice \rightarrow take 0 time point,
start timing when thin walled tube
put in prewarmed 9600 at 68°C

1.6 % PAGE

~ 44 watts (Volts range from 1600
got ~ 12 cm/hr for BPS

Witnessed & Undersigned by me,

Date

Initialed by

Date

T Page N

Polansky

6/9/95

Recorded by

6-2-95

check 5433mer (P44) on P5


Project No. _____
Book No. _____

Fig N _____

3433mer is $\sim 10 \times 10^6$ cpm/ μ l
 \downarrow dilute $1/1000$

1.5×10^{-6} M
55 μ mol of 23mer ~ 20 μ mol nt
5 μ l H₂O
Full

10 mM dATP stock



100 mM stock
1:5 dilution
10 μ l stock, 40 μ l H₂O
TE

10 mM ATP
10 mM dATP

10 mM 23mer
red

5mer
1/1000

1	2	3
2		2
	2	2
		2

RT 12.5 min
dATP 5.5 min

Read & Understood by me,

BooLap

Date

6/9/95

Initiated by

Recorded by

Date

6-2-95

To Page No. _____

From Page No. _____

Result:

1. ice does not shut down 3' exo (start with Mig. rept (?))
2. there is very rapid loss of first 1-10 uts, then very slow degradation. ie next time, need less enzyme
(and/or shorter time points)
3. No apparent ends or 5' exo activity for TFI alone as seen on P14. Here its 5' TFI (4.33 u/l) per 50 μ l compared to 2 μ l on P14 but still should see plenty of loss of full length here based on rate seen on P14. also Wt. 1.5 but on P14 try \pm JNTPs
also P14 has JNTPs present but not here. Maybe some kind of primer extension involved in loss of full length primer - ie extension to many long contributes to apparent loss of primer or to production of a 5' exo target - could be primer primer extension or hairpin within 33 mer.

To Page No. _____

Witnessed & Understood by me,

D. Olamp

Date

6/9/95

Invented by

Recorded by

Date

6-2-95

Primer degradation by Tfi/Vent

Project No. _____

Exhibit 102

Book No. _____

Appl. No. 09/558,421

1

Fig N. _____

95-6/6/95

Purpose: To measure 3'→5' exonuclease activity of Tfi/Vent using the primer degradation assay.

Background: An earlier trial of this experiment (NB 10 page) was done by removing aliquots of the rxn at time points 0, 5, 10, 20, 40, 60, 80, 100 min. The primer was degraded almost to the maximum amount by 10 min. Since we want to determine the 3'→5' exo activity rate, we need to find the linear range of the assay. This can be done by taking shorter time points or by taking a single time point on a series of enzyme dilutions (doubling [enzyme] should double extent of degradation in the linear range of the assay). We'll do this trial expt w/ just 1 enzyme sample - Tfi/Vent in LTI SB - and 13 different dilutions. Once the linear range is found, we can repeat the exp. just on that range.

Materials: ^{32}P dATP for end labeling primer
primer = 33mer correct
Taq dilution buffer - cc aliquot
LTI storage buffer - RL aliquot
Tfi/Vent enzyme mix - from stability study
9600 PCR machine & tubes
5X Cheng buffer - cc aliquot
PNK = T4 Kinase & buffer - ~~set~~ ^{LTI} - new
8% sequencing gel & buffer - LTI premade
① Stop buffer
sterile H₂O

To Page No. _____

Read & Understood by me,

Polarep

Date

6/9/95

Invented by

Recorded by
C. M. C. C.

Date

6/9/95

Project No. _____

Book No. _____

TITLE _____

2

From Page No. _____

Procedure

4/5 ① end-label the 33-mer primer w/ ^{32}P γ ATP = Kinase reaction

reference = notebook 10 p. 12

mix:	H ₂ O sterile	25ul	+ 25ul	added after 1st 5' co. were added
	5x Kinase buffer	12ul	+ 12ul	
20uM	33mer constant	15ul	+ 15ul	
AC9521	^{32}P γ ATP 10mCi/ml ref	5ul	+ 5ul	
	PNK 1u/ul new F5419	3ul	+ 3ul	
		60ul	120ul	

incubate 37°C 30min ✓

~55°C, 5min ✓

store labeled DNA + unincorporated label at -20°C - run some on

4/5 ② make a ~~8%~~ denaturing-sequencing gel + buffer → both pre-made by

- 1, 75ml bottle of 8% mix (cold room) + 450ul ^{32}P AP (made fresh)

0.0868g AP

$$\frac{0.0868g}{0.868g/mL} \Rightarrow 0.100mL$$

- after pouring gel, shake up remaining gel mix so it can be used to fill leaks, ect.

- store gel upright, ON, at RT w/ H₂O-soaked towels and saran ✓

T Pag Nt

With ss d & Underst d by m ,

Dat

6/9/95

Invent d by

R c rded by

Raymond P. Smith

Date

6/9/95

D. Polansky

Fig N .

↑ up enzyme reactions, 14 runs

For 10 runs:

89.64 H₂O ✓

laster mix A = ~~87~~ 88 ul H_2O

6. gal	³² P	3	3	✓
--------	-----------------	---	---	---

rx rxn \Rightarrow 5 μ l H_2O 3 μ l ^{32}P labeled 33-mer from 6/5
 0.4 μ l ^{32}P 35-mer 30 μ l 5X Cheng buffer

32nd, 51 Cheng

0.4 ul 10³ 35 ul 30 ul 5x Chery buffer
2 ul 5x Chery 20 ul, enough for 15 rns
8 ul store on ice

1256


enzyme dilutions: TFI/vent in SB LTI diluted in Tqg dilution buffer

2. add Tag dilution buffer in tubes 1-13, big tubes

- add 20 μ l TFI/Vent to tube #1 = 1:2

- ~~vortex~~^{flick} to mix and store on ice

- remove soul from tube #1 into tube #2 = 1:4 ext



- store dilutions on ice

dd sul Mix A to 9600 tubes # 1-14 * radioactive and prewarm to 68°C in 9600 machine

Tube #	1	2	3	4	5	6	7	8	9	10	11	12	13	14
--------	---	---	---	---	---	---	---	---	---	----	----	----	----	----

Fix A 8ul — equilibrated to 68°C

[illegible]

add enz. to tube, flick to mix, stagger rxns by 1 min

for 15min, remove tube and add 5ul stop-pipet up + down 3x
and keep on ice until gel is run




gel, run
~~15, 1-15~~
~~14, 1-13, 14~~ (3ul)
~~old labeled 32 P labeled 33-met from 6/1/95~~

at before loading, 95°C 5 min in 9600
in down 11:20 am - 1700 V \Rightarrow 65 W

- load
1. 6
9.5
2. 6
3. 6
4. 6
5. 6
6. 6
7. 6
8. 6
9. 6
10. 6
11. 6
12. 6
13. 6
14. 6
15. 6
16. 6
17. 6
18. 6
19. 6
20. 6
21. 6
22. 6
23. 6
24. 6
25. 6
26. 6
27. 6
28. 6
29. 6
30. 6
31. 6
32. 6
33. 6
34. 6
35. 6
36. 6
37. 6
38. 6
39. 6
40. 6
41. 6
42. 6
43. 6
44. 6
45. 6
46. 6
47. 6
48. 6
49. 6
50. 6
51. 6
52. 6
53. 6
54. 6
55. 6
56. 6
57. 6
58. 6
59. 6
60. 6
61. 6
62. 6
63. 6
64. 6
65. 6
66. 6
67. 6
68. 6
69. 6
70. 6
71. 6
72. 6
73. 6
74. 6
75. 6
76. 6
77. 6
78. 6
79. 6
80. 6
81. 6
82. 6
83. 6
84. 6
85. 6
86. 6
87. 6
88. 6
89. 6
90. 6
91. 6
92. 6
93. 6
94. 6
95. 6
96. 6
97. 6
98. 6
99. 6
100. 6

300 d.l. load g.c.
2nd + 598 ul H₂O
3ul

To Page No. _____

Read and Understood by me, 	Date 6/9/95	Invented by  6/9/95	Date 6/9/95
		Recorded by  6/9/95	

Project No. _____

Book No. _____

TITLE _____


From Page No. _____

- B.G.D. ran down to bottom glass clamp on gel rig, ~ 1 hr, ^{constant} 1700 V, ~ 6.5W
- transfer to ^{drain buffer} Whatman paper
- cover w/ saran & cut to size
- dry - 2 pieces whatman under gel, saran over gel, dry ice in trap
- set vacuum & heat for 1 hr = 12:45 - 1:45 PM
- set in phosphorimager cassette - bottom to bottom w/ saran ~ 2 PM for ON exposure

Result: ON exposure on phosphorimager

csc

C:\DATA\CC.GEL 1995:06:07 07:57:48, Range = 0.11-10000.00 Counts, 0.50x



Conclusion: The 1/4, 1/8, 1/16 dilu
gave span the linear range
the primer degradation
Now, we'll do a cou
time pts of each
dilution to gather
better ^{data} from the lin
range. The data w/
used to show ³⁻⁵ exo ac
for stability study

Witnessed & Understood by m

Polansky

Date

6/9/95

Investigated by

Polansky 6-9-95

Recorded by

Paula Condo

Date

6/9/95

T Page No

Page N. — 6/6/95

purpose: To determine the relative mobility of the ^{32}P -33mer correct primer on a PEI plate, developed in LiCl.

Background: Originally we wanted to determine the specific activity of the ^{32}P -33mer primer that was used in the primer degradation assay w/ Tfl/Vent (NB 10 p). However, we later decided that it is not important to find the specific activity since we can do a no eny. control each time the assay is done. Now we want to determine the mobility because we observed that the ^{32}P oligo did not run as expected on the PEI plate, and we just are curious about how where the oligo ran.

materials: cold 33-mer correct
 ^{32}P 33mer correct - labeled on 6/9/95
 ATP
 PEI plate
 1M LiCl
 scint vials: cocktail

procedure: > spot on PEI plate →


20mM ATP 2ul 2ul
 20mM 33mer 2ul 2ul
 1/1000 ^{32}P -33mer 2ul

- > set plate in developing chamber w/ 1M LiCl - 1hr - ran 1/2 way
 > circle control spots (ATP & cold 33mer) in lane 3
 under UV light → sketch of how plate looked

8 ATP ^{32}P
 ADP
 ^{32}P -33mer*
 cold 33mer
 PPi ^{32}P

cut

- > lane 3 into 8 pieces and
 count in scintillation counter



expected to stay
 at base, but ran
 near the moving
 front

To Page No. _____

Designed & Understood by me,

Polansky

Date

6/9/95

Invented by

R. Corded by
 Paulen & P. M. K.

Date

6/9/95

Project No. _____

Book No. _____

TITLE _____

8

From Page No. _____

CSC

Result:

PAGE:

USER: 1 ID:32P 1.0 CPM PRESET TIME: 1.00 TUE 06 JUN 1995 15:36
 SAMPLE REPEAT: 1 CYCLE REPEAT: 1 SCR:N RS232:N
 H#: 0 ABC:N BCF:N RCM:N
 CHANNEL 1-LL: 0 UL:1000 2SIGMA: 0.05 BKG SUB: 0.00 BKG 2SIG: 0.00 LSR:
 DATA CALC: CPM. UNKNOWN REPLICATES: 1 NORM FACTOR: 0 1.00000
 HALF LIFE(DAYS):N

SAM	CPM1	TIME	EF
1	111.00 - baseline	1.00	
2	307.00 - area around where cold ATP	1.00	
3	136.00 standard ran	1.00	
4	117.00	1.00	
5	60.00	1.00	
6	215.00 - area around where cold 33mer primer ran	1.00	⇒ primer runs like dAMP
7	27.00 - near solvent front	1.00	

The labeled 33-mer primer ran like the dAMP runs on a PEI plate, up near the solvent front, not at the origin.

This result was expected from the information we heard from the chemistry group: DNA stays at the origin because of its large size.

oligonucleotides run like dAMP because they have same charge:mass ratio as dAMP

To Page 1

Witness d & Und rsto d by m ,

Dat

6/9/95

Invent d by

6-12-95

Dat

6/9/95

R c rded by

Paulen Poub

1/Vent primer degradation assay - time course on serial dilutions for a time point of stability study

pose: To measure 3'-5' exonuclease activity of TFI/Vent mixes in 3 buffers at the zero time point of stability study

background: linear range of assay determined in previous expt, NB II page 1 - dilutions 1/4, 1/8, 1/16 looked good

- linear range $\hookrightarrow 0.045$ units

- now do time course of these dilutions

1/4 dilution $\rightarrow 0.09 \frac{\mu\text{g}}{\mu\text{L}} \times 2 \mu\text{L} = 0.18 \mu\text{g}$ in 100L reactions
1/8 dil $\rightarrow 0.0225$ units
1/16 $\rightarrow 0.01125$

materials: ^{32}P 33mer correct - labeled on 6/5/95

TFI/Vent in LTISB, epicenter SB, epicenter TFI

8% gel
mix A

* Vent dilution 2000 u/mL
200 uL vent Lot 17 2/24/95
42.4 uL TAE 513
44.4 uL at 0.09 u/uL
22.2 x dil

procedure:

make mix A, enough for 14 rxns - 90 uL per rxn

per rxn $5.6 \mu\text{L} \times 14 = 78.4 \mu\text{L}$ $20 \mu\text{L} \times 14 = 280 \mu\text{L}$ 5x Cheng

TFI/Vent in LTISB

" in epicenter SB

" in epicenter TFI

4 uL $\times 14 = 56 \mu\text{L}$

66 uL $\times 14 = 924 \mu\text{L}$

90 $\times 14 = 1260 \mu\text{L}$

^{32}P 33mer 500 uL stock

H₂O

mix A

vent

each enzyme/buffer mix

mix A

H₂O

undiluted enzyme (0.09 u/u)

90 uL

90 uL

90 uL

8

6

4

2

4

6

take 3, 6, 9, 12 min time points by removing 10 uL of rxn to 5 uL stop in small tubes, keep on ice

start rxns 1 min apart time on clock - 0 1 2 start

1. no enzyme control

2. 33mer

3. H₂O

4. 33mer + H₂O

5. 33mer + H₂O + 2 uL stop

6. 33mer + H₂O + 2 uL stop

7. 33mer + H₂O + 2 uL stop

8. 33mer + H₂O + 2 uL stop

9. 33mer + H₂O + 2 uL stop

10. 33mer + H₂O + 2 uL stop

10 uL Cheng

2 uL ^{32}P 33mer

33 uL H₂O

30 uL

30 uL

30 uL

30 uL

30 uL

30 uL

30 uL

3

4

5

6

7

8

9

10

11

12

13

14

15

16

17

18

19

20

21

22

23

24

25

26

27

28

29

30

31

32

33

34

35

and Und rstood by me,

Date

Invented by

Date

Polansky

6/9/95

Recorded by

4/19/95

4/19/95

To Page No. _____

Project No. _____

Book No. _____

TITLE _____

From Page No. _____

Note: samples Vent/Tf1 in L.T.I. 3 min 2nd > only 5th of can was stopped → load mass 4.5 instea
3 min 4th

4:50 PM - ~6:15 PM

1700 V constant, ~45 W, gel was dried & put in P.I.

* 39 & 40 may be underloaded due to problem: expelling full vol. from

order: control, 1 → 48 where 1-12 are Tf1/Vent in L.T.I. SB

of samples on gel

13-24 are Tf1/Vent in epicenter SB

25-36 are Tf1/Vent in epicenter Tf1

37-48 are Vent alone

T Pag N.

With ss d & Und rsto d by m ,

D Polans

Dat

6/9/95

Invent d by

R Cord d by

Randyn Comb

Date

6/9/95

Project No. _____

Book No. _____

TITLE _____

48

From Page No. _____

Repeat unit assay for TF1/V.

incorp of p.30 for stability of
Units TF1/Vent mixes of A J-
with 5 duplicate dilutions
for optimized signal/noiseUse 2 μ l of 1/250 dil
and linearity vs (units)

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20

TF1 Chem unit
assaying (same
as P. 30)

47 →

TF1 4.33 μ l (P. 8)
1/250 dil *

2 →

TF1/Vent LTISB

5-13-95 (Ning Guan &
1/250 notebook)

2 →

TF1/Vent Epicent SB

5-16-95 1/250 dil

2 →

TF1/Vent (Epicent TF1)

5-16-95 1/250 dil

2 →

Vf 50 μ l

74°C, 10'

1 μ l 0.5M EDTA → spot 40 μ l on* all dilutions are done as 5 separate dilutions
of 2 μ l Enzyme into 498 μ l Tag dilution buffer

T Pag

Witness d & Understood by me,

Dat

Inv nted by

Dat

R c rded by

O O Olamp

6/9/95

6-9-95

Project N _____

Book N _____

49

Page N _____

1	4163.00
2	4429.00
3	4636.00
4	4646.00
5	4349.00
6	4550.00
7	4529.00
8	4623.00
9	4350.00
10	4315.00
11	3995.00
12	4339.00
13	3732.00
14	4695.00
15	4428.00
16	3975.00
17	4584.00
18	4541.00
19	4297.00
20	4412.00
21	259.00
22	84613.00
23	87557.00

T Page No. _____

Received & Understood by me,

Date

Initiated by

Date

Recorded by

Test run of PEI plates - prior to turnover exp

No. 195

Purpose: To test how well a fresher batch of PEI plates can resolve dAMP from dATP and how tight/clean the spots are. This is being done prior to using this batch of plates for another TFI/Vent turnover experiment

Background: Last time the dAMP spot did not resolve well from the yellow "junk" that runs near the 1M LiCl solvent front, making it difficult to cut out + count just the dAMP (without Pi) for accurate turnover results. We'll try washing a different batch of plates - from Jesse - in dH₂O & drying them 1st, before running samples and compare to unwashed plate. Also we see if running a whole plate gives better resolution than a 1/2 plate.

Materials: PEI plates - from Jesse ← Mackerey Nagel Polygram cell 30
Kill soln = ~20mM dATP
20mM dADP
20mM dAMP
100mM EDTA
1M LiCl - Fresh, see recipe on next page
Aldrich cat # Z12255-2
PEI/UV

To Page No. _____

Read & Understood by me, _____	Date _____	Inv. noted by <u>[Signature]</u> <u>6/19/95</u>	Date <u>6/12/95</u>
		Recorded by <u>Carolyn E. [Signature]</u>	

see P14 apparent exo present
in PCR buffer Project No. _____
50 P44, no exo in Book No. _____
From Page No. _____

Repeat apparent exo result for TFI
on P14 with different primer
~~to home + hand check~~

Exhibit L-107
Appl. No. 09/558,421

① ② ③ ④ ⑤ ⑥ ⑦ ⑧

10x PCR buffer
50 mM MgCl₂

5 →
10.5 μl →

5 ✓
15 ✓

5x Cheung (no dNTPs)
(see P200)

10 μl → 10 ✓

32P 33mer correct (P51) 5 μm
32P 23mer "AC" 5 μm
32P 42mer "fidel" 5 μm
see P12, 14 for method

2 2 2 2 2 ✓
2 2 2 ✓

Tag storage buffer

10 10

TFI 4.33 μl (P.8)

X10 →

H₂O

31.5 → 28 → 31.5 28 ✓

V_f = 50 μl

74°C. remove 10 μl to 5 μl cycle seq stop sol
at 15 min 60 min

run on 8% PAGE

* Zero time point:

1. mix buffer^{mix}, 32P primer (MgCl₂ if needed) and H₂O. Volume =
2. remove 8 μl to 2 μl Tag storage buffer and 5 μl cycle seq for 0 (no enzyme) time point.
3. now have 32 μl of reaction^{mix} left. preheat to 74°C, add 8 μl TFI so V_f = 40 μl again and remove 10 at 15 and 60 min to 5 μl cycle seq stop sol.

To Page 1

Witnessed & Understood by me,

Deanne Boland

Date

6/1/96

Invented by

Recorded by

Dat

6-13-97

32P

oligos (follow P12, 14 for 5' end)

Project No. _____

Book N. _____

51

tag N

5x kinase buffer
"correct" 20 μ M
13P9

1

4 μ l
5 μ l

2

3

✓
✓mer AC(P1365)
100 μ M1 μ l

✓

← 23mer has termin
A instead of G
its called "AC".mer "fidel" 6-13-95 "old Temp" This
is different from
Fidel Temp P541 μ l

✓

ATP 10 μ M (ref 6-16-95)
VR
H₂O2 μ l
1 μ l
8

✓

← (5x less ATP)
(than on P12)V_P = 20 μ l

12

12

✓

37°C, 35', 5', 55°C

2 mer fidel

51351 CAC (012)

56.89 nmoles primer
56.89 μ l H₂OF = 100 μ M

To Page No. _____

Read & Understood by me,

Cecilia Polansky

Date

6

19/95

Investigated by

Recorded by

Date

6-13-95

Incorporation and turnover by Vent alone

Time course on 0.1, 0.2, 0.3 units Vent

ag No. _____

Purpose: To measure both incorporation and turnover by varying amounts of Vent over a 1 hr time course. These values will be compared to those of TFI/Vent mixes. Vent alone is a control for the stability assays of TFI/Vent mixes.

Background: NB 10 p. 41, turnover was not above background for Vent alone, when 0.09 units were used in a 20 min time course. However, ^{this time well by more exp.} Roger had observed turnover over background when 0.15 units were used 100ul rxn. Turnover is expected to increase linearly w/ time even after incorporation has stopped, if Vent can bind at the nick. New PET plates may give better results than last time.

Materials: PET plates from Jesse, tested on NB 11, p. 15, Vent 2u/ul
activated DNA = gapped DNA made w/ DNase I, Kill sol'n = NB 11 p. 14
^{80 mM dAm, DT, 100 mM EDTA}
³²P dATP - ref date 6/16/95
mix A: per 1, 100ul rxn: { 63.786 ul H₂O
20ul 5x Cheng buffer
13.5ul activated DNA, 3.7 mg/ml, Cf = 0.5
0.5ul dATGC-TP, 10mM each, Cf = 50u
0.214 ul ³²P dATP 10mCi/ml

make enough mix A for today + the next expt - TFI/Vent turnover
16x 98ul/rxn = 1568ul

(for 16 reactions)

A { 1020.6ul H₂O go 1ml + 20.6ul ✓
320ul 5x Cheng ✓
216ul DNA ✓
8ul dATGC-TP - BRL lot FBH001 ✓
3,424ul ³²P dATP
VF 1569ul use 98ul/100ul rxn

7.5×10^7
 4.7×10^6 CPM total
294-250/100ul reaction
 7×10^6

Vent dilutions: - 2/14/95 Lot 13 spch + mix of Vent
LTISB = (1005B)

1st dilute 5x Cheng → 10ul 5x Cheng + 40ul H₂O - mix

① dilute Vent stock (2u/ul) to 0.15u/ul → 2ul Vent + 24.66ul 1x Cheng LTISB

② dilute ① to 0.1u/ul → 10ul ① + 5ul 1x Cheng LTISB

③ dilute ② to 0.05u/ul → 5ul ② + 5ul 1x Cheng LTISB

Read & Understood by me,

Ernest Polary

Date

6/19/95

Invented by

Recorded by

Date

6/13/95

Page No. _____

Project No. _____

Book No. _____

TITLE _____

18

From Page No. _____

procedure: - deliver 10ul of kill soln to 1.5mL Cppendorfs x 1-
 - label tops 10' 20' 30' 40' 50' 60' & set on ice until

rxn x

1

2

3

4

stop tube x

1- ~~7~~67-12
~~7-14~~13-18
~~15-21~~19-24
~~20-28~~

mix A

98ul / rxn

prewarm to 68°C

→

0.15 u/l Vent
 Lot # 17
 opened 2-24-95

0.10 u/l Vent

0.05 u/l Vent

2ul

start rxns by adding enz., keep at 68°C in 9600

2ul

thin in
 n the mix level in
 of vent as TPI/Vent
 was not

2ul

thin
 0.2 u/l
 and TPI
 in 0.1
 Vent


no enz

2ul + change
 LTI SB
 = Tag SB


100ul

- remove 10ul rxn to stop tubes at each time pt (10', 20', 40', 50')
- spot 10ul / GFC & 2ul PET
- 10ul of mix A x 3 (for determination of specific activity)

Test of old (Baker) PET plates - 2ul, cold kill soln spotted



- old plate → dAMP runs w/
 solvent front junk



EC

- new plate → dAMP runs ~1/2
 between solvent front & or

T Pag N

Witness d & Und rst od by me:

Deborah Polak

Date

6/19/95

Invent d by

R cord d by Paulm Pomb

Dat

6/13/95

6-19-95

Project No. _____

Book No. _____

TITLE _____

Incorporation (pmoles)

From Page No. _____

$\left(\frac{\text{cpm}}{\text{specific activity}}\right) \left(\frac{100 \mu\text{l rxn}}{10 \mu\text{l spot}}\right) \left(\frac{20 \mu\text{l}}{10 \mu\text{l}}\right)$

0.3 mls	55	21752.00 - 936 pmoles
	56	40494.00 - 1742
	57	50701.00 - 2745 2181
	58	63810.00 - 2723 2745
	59	63423.00 - 2723
	60	61923.00 - 2663
1.3 mls	61	12710.00 - 547
	62	21727.00 - 934
	63	32040.00 - 1378
	64	39939.00 - 1718
	65	43064.00 - 1852
	66	51401.00 - 2211
0.1 mls	67	9060.00 - 390
	68	14810.00 - 637
	69	19948.00 - 858
	70	24421.00 - 1050
	71	31940.00 - 1374
	72	30490.00 - 1311
no engr.	73	420.00
	74	540.00
	75	299.00 $\bar{x} = 348$
	76	310.00 $n = 6$
	77	197.00
	78	323.00
	79	923719.00 10ul spot of mix A
	80	973931.00 $\bar{x} = 929,46.2$ for specific activity
	81	890737.00

The ^{observed} specific activity of mix A is 2x higher than anticipated by the following calculation:

$$\frac{10 \mu\text{L Li}}{\mu\text{L}} \times \frac{2.2 \times 10^6 \text{ cpm}}{\mu\text{L}} = \frac{2.2 \times 10^7 \text{ cpm}}{\mu\text{L}} \times \frac{3.4 \mu\text{L in mix A}}{1568 \mu\text{L}} = \frac{8.7 \times 10^7 \text{ cpm}}{1568 \mu\text{L}} = 5.6 \times 10^4 \text{ cpm}$$

0.925 2 days to ref. date

$$\frac{5.6 \times 10^4 \text{ cpm}}{\mu\text{L A}} \times 10 \mu\text{L A spotted} = 5.6 \times 10^5 \text{ cpm expected}$$

$$\frac{9 \times 10^6 \text{ cpm}}{9.3 \times 10^5 \text{ cpm observed}}$$

I don't know where the error came from, but the results should still be consistent within this experiment

To Page No

With ss d & Understood by m ,

Demetrius Polaris

Dat

6/19/95

Inv nted by

[Signature] 6-19-95
 R c rd d by
Pauline P. [Signature]

Dat

6/14/95

purpose: To ^{more} accurately determine turnover by T₄I/Vent in LTT SIB

- 5 replicates of the 5min and 10min ^{± 5min} time points, which are within linear range of the assay. This data will be used for the stability study.

materials: new PET plates
mix A from 4/13/95
TFI/Vent in LTI SB
kill soln from NIB 11 page 14 = 20mM each dA-MDT-P
100mM EDTA

	1	2	3	4	5	6
2400 tubes						
100 tubes						
50 tubes						
25 tubes						
12.5 tubes						
6.25 tubes						
3.125 tubes						
1.56 tubes						
0.78 tubes						
0.39 tubes						
0.19 tubes						
0.09 tubes						
0.04 tubes						
0.02 tubes						
0.01 tubes						
0.005 tubes						
0.002 tubes						
0.001 tubes						
0.0005 tubes						
0.0002 tubes						
0.0001 tubes						
0.00005 tubes						
0.00002 tubes						
0.00001 tubes						
0.000005 tubes						
0.000002 tubes						
0.000001 tubes						
0.0000005 tubes						
0.0000002 tubes						
0.0000001 tubes						
0.00000005 tubes						
0.00000002 tubes						
0.00000001 tubes						
0.000000005 tubes						
0.000000002 tubes						
0.000000001 tubes						
0.0000000005 tubes						
0.0000000002 tubes						
0.0000000001 tubes						
0.00000000005 tubes						
0.00000000002 tubes						
0.00000000001 tubes						
0.000000000005 tubes						
0.000000000002 tubes						
0.000000000001 tubes						
0.0000000000005 tubes						
0.0000000000002 tubes						
0.0000000000001 tubes						
0.00000000000005 tubes						
0.00000000000002 tubes						
0.00000000000001 tubes						
0.000000000000005 tubes						
0.000000000000002 tubes						
0.000000000000001 tubes						
0.0000000000000005 tubes						
0.0000000000000002 tubes						
0.0000000000000001 tubes						
0.00000000000000005 tubes						
0.00000000000000002 tubes						
0.00000000000000001 tubes						
0.000000000000000005 tubes						
0.000000000000000002 tubes						
0.000000000000000001 tubes						
0.0000000000000000005 tubes						
0.0000000000000000002 tubes						
0.0000000000000000001 tubes						
0.00000000000000000005 tubes						
0.00000000000000000002 tubes						
0.00000000000000000001 tubes						
0.000000000000000000005 tubes						
0.000000000000000000002 tubes						
0.000000000000000000001 tubes						

1/Vent mix	2ul	- add enzyme to start rxns
TSB	2ul	

Zul

Zul

zul

~~2nd no enzyme control~~
~~2nd LTI SB = Tag SB~~

incubate at 68°C

100ul

> At 5, 10, 15 min remove 200 μ l rxn to the 2nd kill soln in stop tubes, mix well & keep on ice.

> spot row/GEC filter $(1-15) \times 2 = 30$ spot rxn G, no enz control 4x per time pt
 > 2w/PEI plate $(1-15) \times 2 \times 2 = 54$ = 12 spots

4 = k

2nd / 1st plate

16-1-2-3, 4
17-1-2-3, 4 18-1-2-3, 4

= 102

To Page No.

ed & Understood by me,

Dat

Invented by

Date _____

ruel a Polak

6/19/95

Recorded by

6-15-95

Date 6/14/95

22 Results:

Book No. _____

TITLE

BOOK NO. _____
Turnover (cont)

$$= \frac{(cpm - \text{background cpm}) / \text{specific activity}^*}{\left(\frac{40}{20}\right) \left(\frac{100}{2}\right)} = 109$$

SAM		CPM1
odd = dAMP		
even = dADP + dATP		
1	> 5	950.00
2		72035.00
3	>	1727.00
4		71656.00
5	> 5	2489.00
6		67554.00

replicate 2	7	1176.00
	8	72213.00
	9	1768.00
	10	69650.00
	11	2189.00
	12	69544.00

replicate 3	13	5	1219.00
	14	7	1519.00
	15	10	1784.00
	16	7	0384.00
	17	15	2314.00
	18		69324.00

19	1189.00
20	71061.00
21	1914.00
22	72483.00
23	2096.00
24	67757.00

25	1032.00
26	70489.00
27	1714.00
28	64939.00
29	2289.00
30	64174.00

31	612.00
32	75976.00
33	494.00
34	79452.00
35	421.00
36	79068.00

37	587.00
38	78840.00
39	460.00
40	80087.00
41	481.00
42	78697.00

43	419.00
44	78893.00
45	404.00
46	76101.00
47	450.00
48	68672.00

49	549.00
50	67776.00
51	498.00
52	73381.00
53	559.00
54	72680.00

55	620.00
56	1714.00

Witnessed & Understood by me, _____

Date 6/19/95

Inv nt d by *K. M. L.* 6-19-55

Dat
6/14/95

Record by *David R. Smith*

Page No. _____	Incorporation (pmol)
	(cpm/specific activity) $\left(\frac{100\mu\text{L rxn}}{10\mu\text{L spot}}\right) \left(\frac{40\mu\text{L}}{20\mu\text{L}}\right)$
5' 93072.00	93072/419 $\left(\frac{100}{10}\right) \left(\frac{40}{20}\right) = 4443$
10' 107957.00	-5153
15' 140583.00	-6710
5' 107888.00	-5150
10' 116159.00	-5545
15' 157153.00	-7501
5' 89224.00	-4259
10' 129878.00	-6199
15' 158185.00	-7551
5' 86678.00	-4137
10' 129770.00	-6194
15' 146342.00	-6985
5' 71757.00	-3425
10' 127388.00	-6081
15' 158825.00	-7581
no. 285.00	background
5' 355.00	
no. 291.00	$\bar{x} = 301 \pm 31 \text{ cpm}$
10' 300.00	214.4 pmoles
no. 310.00	$n = 6$
15' 262.00	
839570.00	oul mix A
831885.00	spotted 3X
840299.00	$\bar{x} = 837,251$

$\bar{x} \pm 1SD$ Incorporation

$n = 5$

5' $\Rightarrow 4283 \pm 620 \text{ pmol}$

10' $\Rightarrow 5834 \pm 467 \text{ pmol}$

15' $\Rightarrow 7266 \pm 395 \text{ pmol}$

$$\text{specific activity} = \frac{837251 \text{ cpm} \times \frac{100\mu\text{L rxn}}{10\mu\text{L spot}}}{(5000 \text{ pmol}) 4} = 419 \text{ cpm/pmol (nt)}$$

To Page No. _____

Read and Understood by me,

Michael Polansky

Date

6/19/95

Inv. nted by

Recorded by

Carlynn Lamb

Date

6/14/95

Project No. _____


Book No. _____

TITLE _____

24


From Page No. _____

PROOFREADING: TFI/VENT



- 6/15/95
cc
- Background was 118 pmd, so the best signal to noise occurs at 15 min 13.6x background. However, by 15 min, incorporation is slowing down as gaps are filled in. At 15 min some turnover is occurring at nicks - not a good model of a PCR reaction.
 - Trade off between good model and linearity of the time points is good.
 - A 10% drop in activity would be detected by this assay, using 5 replicates.

POLYMERIZATION: TFI/VENT



- Incorporation falls off after 5 min, because gaps are filled by high TFI polymerase activity. After that, turnover occurs at nicks.
- % turnover increases because DNA synthesis is slowing down while turnover keeps going at the same rate.

T Page 1

Witnessed & Understood by me,

Date

Inventor by

Date

Danica Pokany

6/26/95

Recorded by

Paula Pomb

6/25/95

sig N	Turnover (pmoles) dAMP	% Turnover
u115	$\frac{(\text{cpm} - \text{background cpm})}{\text{specific activity}} \left(\frac{30}{10} \right) \left(\frac{100}{2} \right)$	$\left(\frac{\text{pmoles dAMP}}{\text{pmoles incorporation} + d.t.} \right) \times 100$
619 Turnover CPM1 = dAMP, dATP = dAMP	ex. $\frac{(1776 - 807)}{465} \left(\frac{30}{10} \right) \left(\frac{100}{2} \right) = 20.8$	ex. $\frac{208}{208 + 936} \times 100 = 18.2$
710 1776.00 86470.00		
720 2994.00 86141.00	470	21.2
730 4209.00 91512.00	732	25.1
740 4983.00 85588.00	898	24.7
750 6822.00 94359.00	1290	32.1
760 7013.00 85869.00	1330	33.3
770 1216.00 79679.00	88.0	13.9
780 2179.00 83426.00	295	24.0
790 2954.00 81631.00	462	25.0
800 3716.00 83944.00	623	26.6
810 4469.00 85258.00	788	29.8
820 5283.00 87259.00	963	30.3
830 1223.00 85430.00	89.5	18.7
840 1807.00 90067.00	215	25.2
850 2316.00 88894.00	325	27.5
860 2953.00 84914.00	462	30.6
870 3572.00 90268.00	595	30.2
880 3815.00 92711.00	647	33.0
890 792.00 81173.00	\bar{x} dAMP = 807 cpm is background = 174 pmoles $\frac{929,462 \text{ cpm} \times 100 \mu\text{l}}{(5000 \text{ pmoles})(4)} = 464.7 \text{ cpm}$ pmole nt DNA	
900 749.00 87079.00		
910 655.00 87371.00		
920 880.00 86383.00		
930 785.00 86929.00		
940 978.00 90674.00		
950		

sed & Understood by m	Date	Invented by	Date
revised Polamp	6/19/95	6/19/95	6/19/95
Recorded by			
Condon Condon			

To Page No. _____

TITLE

R c r d d by

Primer degradation by TFI/Vent-Cpicenter
time course & 5 replicates on 12% gel

B k No. _____

eN. _____

pose: To measure 3'-5' exonuclease activity of TFI/Vent mix using the primer degradation assay. The data will be part of the stability study on the mix. The following changes to the assay will be made in hopes of obtaining more accurate data on the rate of primer degradation:

- 1) 12% gel will be used instead of an 8% gel - may give better peak resolution. Last time, peak shoulders and double peaks were a problem during quantitation.
- 2) 44mer^{Fidel} will be used instead of 33mer - correct - the 33mer could form primer dimers &/or a hairpin that may have altered the degradation rate.
- 3) By doing 5 replicates, we can assess the accuracy of this assay as compared to the turnover assay.

Primer degradation by Vent alone will also be measured.

Background: Note - although we tested TFI/Vent mixes on p. 9 NB11, only the TFI purchased from Epicenter is TFI LTI's TFI is really Tth. That is why we are not doing any stability tests on LTI's enzyme, until a new TFI alone is obtained.
- the amount of enzyme and time course of this experiment are known to be in linear range of assay from the earlier expt on p. 9 NB11.

Materials: γ -³²P-end labeled 44-mer Fidel
fresh mix A

12% denaturing sequencing gel

stop sol'n

TFI (Vent) mix made w TFI purchased from Epicenter
Vent diluted w/ TAE Epicenter 5B, Vent 10 \times 17

recipe for 12% gel, 100mL: 48g urea

30mL, 40% acrylamide:bis mix
10mL, 10 \times TBE

dissolve by stirring & low heat
+ 600 μ L 10% AP (100mg/L)

qs to 100mL w/ H₂O - squirt bott.

30 μ L TEMED

100mL

To Page No. _____

I & Understood by me,

slamp

Date

6/26/95

Invented by

Recorded by

Carolyn Conk

Date

6/21/95

15

Project No. _____

Book No. _____

TITLE _____

No. _____

A: 90ul per rxn, enough for ¹⁵ rxns

15
60ul x 9 = 540ul H₂O 540ul

20ul x 9 = 180ul 5x Cheng buffer 300ul

4ul x 9 = 36ul 3²P-44mer Fidel, 5uM stock C.F. = 200n

90 x ⁹/₁₅ = 540ul 1350:

procedure:

label the 44-mer Fidel primer - 70.3ul H₂O

24ul 5x Kinase buffer

6ul 44mer Fidel, 100uM stock

13.7ul 3²Pγ ATP, 10uCi/ul

ref. date = 6.

6ul PNK 10uCi

120ul

37°C 30min

55°C 5min

store at -20°C overnight

deliver sub stop soln to 9600 tubes & ¹⁵/₁₋₂₈ label rxn tubes 1-
make mix A

rxn	1	2	3	4	5	6	7	8	9	10	11	12
H ₂ O	8	8	8	8	8	8	8	8	8	8	8	8
mix A	90	90	90	90	90	90	90	90	90	90	90	90

preheat to 65°C in 9600, start by adding

50*
Epi / no
enzyme
2ul

2ul Trf / Vent - Epicenter 2ul Vent.

At each time pt remove rxn rxn to the sub stop soln - heat 90°C

no enzyme 0', 10' replicas: 4', 6' at 75°C
time course 2', 4', 6', 8', 10', 20' w/ Trf / Vent

& Understood by me,

Dat

Invented by

Date

6/26/95

Record d by

4/21/95

Darwin Rom

g N _____

center SB for no enz control: final 100mM NaCl, add solid to LTI SB
and to dilute vent with.

final 50mM Tris

now at 20mM

5 mL Tag SB = LTI SB ✓

91 µL 1M Tris 7.5 - premade by LTI ✓

62.5 µL 1M Tris-HCl ✓ → 157.649/mole

153.5 µL glycerol ✓

0.031 g NaCl (58.449/m) ✓

$$1 \text{ mole} \times 10 \times 10^{-3} \text{ L} = 0.01 \text{ mole}$$

$$\frac{g}{157.649/m} = 0.01 \text{ mole}$$

$$= 1.5764 \text{ g Tris-HCl} \\ + 10 \text{ mL dH}_2\text{O}$$

Vent dilution in Epicenter storage buffer = 22.2x dilution to 0.09 u/L
2 µL Vent stock (2 u/L) w/ p2 pipetman → 1st spin down + vortex Vent stock
+ 42.4 µL Epicenter storage w/ p200 pipetman
44.4 µL vortex to mix

mples were heated to 90°C, 5' in 9600 prior to loading.

sample #'s 1 = no enz 0 min

2 = no enz 10 min

3 = TFI/Vent 2 min - not preheated = may be off

4 = " 4' 21 4 > 2

5 = " 6' 22 6 > 2

6 = " 8' 23 4 > 3

7 = " 10' 24 6 > 3

8 = " 20' 25 7 > 4

9 (1) TFI/Vent 4' 26 4 > 1

10 " 6' 27 4 > 5

11 (2) 4' 28 6 > 2

12 6' 29 6 > 3

13 4' 30 6 > 4

14 6' 31 6 > 5

15 4' 32 6 > 1

16 6' 33 6 > 2

17 4' 34 6 > 3

18 6' 35 6 > 4

19 Vent 4' 36 6 > 1

20 6' 37 6 > 2

2nd load 1, 3, 28, 12
2, 3 - 28, 8
1st load 1, 2, 4 - 28, 2
1, 19 - 28, 2 no times course
2nd run longer

To Page No. _____

d & Understood by me,

Date

Invented by

Date

Solamp

6/26/95

Recorded by

6/21/95

Paula Lomb

Project No. _____

Book No. _____

TITLE _____

om Page No. _____

sample #

1st loading: 1, 9-28, 2 (no time course) ~3:15^{pm}, rate $\frac{0.7 \text{ cm}}{8 \text{ min}} \approx 0.1 \text{ cm/min}$

2nd loading 2, 3-28 w/ time course ~5:15^{pm}, run until $\beta\beta$ reaches bottom. At this time ~10 bases will have run. This will be our whole gel loading to see the whole gel most of the products.

Gel was run at 1700V constants for a total of 5.5 hrs.

To Page N. _____

Used & Understood by me,

Date

6/26/95

Invented by

Recorded by

D. M. M. M.

Date

6/21/95

J. Polak

From Page No.

PCR mix:

LTI 10X PCR buffer
 50 mM MgCl₂
 4 dNTPs, 10 mM each
 M13 RF 1 pg/ μ l Δ see page 42 for dilution
 M13-6301 (anchor), 20 μ M
 H₂O

(40 Rxns)

400 μ l
 120 μ l (1.5)
 80 μ l (2.0)
 40 μ l (1.0)
 80
 3160
 3880 μ l (400 primers)

(18 Rxns)
 (in falcon tube)

[1]

1746 μ l

18

1764 μ l

[2]

1746 μ l15.5 μ l

2.5 μ l 18 μ l
 1764

(5 μ M / 100 μ l)

H₂O
 rTag 5 μ l/ μ l
 Tne 3.6 μ l/ μ l 5 μ l/ μ l

[1]

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16
 98 —————>

[2]

98 —————>

20 μ M primer see page 42-43

M1366R1

2

2

7069

2

2

407

2

2

806

2

2

1491

2

2

2506

2

2

3972

2

2

5464

2

2

100 μ l

make 2 sets

1-16 gets elongation

1'-16' gets elongation

RF

1 pg M13 / 100 μ l
 = 3×10^{15} molecules
 = 2×10^{19} molecules
 = $(2 \times 10^{19}) (6.0 \times 10^{23} \text{ molecules})$
 = 125,000 molecules
 anchors

94°C 1 min \rightarrow 30 cycles: 94°C 15 sec, 53°C 30 sec
 elongation is 6 min for PCR 1-16 and 2 min for 1'-16' \rightarrow cont. on,

T Page 1

Witnessed & Understood by m

Dat

Inv nt d by

Date

6/26/95

R c rd d by

6-23-95

0000101010

Incubate Tag with Cheng buffer
in gap DH P

Project No. _____
B kN. _____

41

N	1	2	3	4
5 x Cheng buffer (containing dNTPs) (all P20, 10)	10	10	10	10
10 mM dNTPs	1	1	1	1
Human spleen genomic DNA 80 ng/μl	1.25	—	—	—
2112 10 μM	1	—	1	—
2113 10 μM	1	—	1	—
5 μl Tag	0.5	0.5	—	—
2115 5 μl	—	—	0.5	0.5
H ₂ O	35.25	33.25	35.25	33.25
	50 μl	—	—	—

T Page No. _____

Read & Understood by m ,
B. Camp

Dat
6/26/95

Inv nt d by
R cord d by

Dat
6-23-95

Fig No. _____

SigmaPlot regression lines:

1) for Proofreading: TFI/Vent the slope = $28.676 \frac{\text{pmole}}{\text{min}}$, $r^2 = 0.978$

2) for Polymerization: TFI/Vent the slope = $417 \frac{\text{pmole}}{\text{min}}$, $r^2 = 0.873$

Units of TFI/Vent proofreading: $28.68 \frac{\text{pmole}}{\text{min}} \times \frac{30 \text{ min}}{10,000 \text{ pmole}/\mu\text{l}} = 0.043 \mu\text{l}$

Units of TFI/Vent polymerization: $417 \frac{\text{pmole}}{\text{min}} \times \frac{30 \text{ min}}{10,000} = 0.03 \mu\text{l}$

expect: 0.045 μl

conclusions

1) The turnover assay can detect a 10% loss of 3' exo activity in TFI/Vent mixes. By repeating the assay more frequently and/or with more replicates the error may decrease so that a 5% loss of activity could be detected.

2) Early time points, before 10 min, reflect turnover during DNA synthesis - the best model of PCR. Later time points reflect turnover during DNA synthesis plus turnover at DNA nicks - not such a good model of PCR. However, the later time points give better data because the signal to noise ratio is higher (3.6x versus ~1.3x early on). Both all 3 time points should probably be done during the stability study.

3) Turnover by TFI/Vent mix is about 3x higher than by Vent alone. This result was observed in an earlier experiment too. TFI may create more mismatches for Vent to turnover than when no TFI is present.

↑
not true: exo in TFI/Vent is only ~2x higher than for Vent alone (see mismatch on p 40 10). In this experiment (P17.11) the Vent alone of 2 μl of 0.1 μl is the one we should compare to the TFI/Vent mix.

To Page No. _____

I & Understood by me,

Wade C. Polay

Date

6/26/95

Invented by

Record d by

Cawlyn Lamb

Date

6/24/95

PLATES

Both the dATP+dADP peak and dAMP peak become more spread out on the PEI plates as the solvent front runs further. The distance between the two peaks becomes greater as the solvent front runs further: ~0.5 cm vs 2 cm

Conclusion: For best resolution of dAMP from dATP+dADP on PEI plates, run the Lill solvent front to the top of the plate, ~16 cm from origin.

To Page No. _____

d & Understood by m ,

Carol Polansky

Date

6/26/95

Inv nted by

Recorded by

Carolyn Combs

Date

6/24/95


Project No. _____

Book No. _____


TITLE _____

from Page No. _____

Dilution of M13 RF to 1 pg/ul in TE: stock lot CN4132, 5 μ g/18.4 μ l
do serial dilutions to 1:2.72 x 10⁵



6/20/95



(M13mp19 prime)

→ 20mers

→ 9G+C

→ 11A-T

→ Tm = ~58°C

primer name

PCR product length (bases)

sequences

M13-6301 anchor (= +strand)

5' GTTTTACAAC

M13-6681 (= -strand)

380

5' TTCC TGTAGCCAGCTTTC

M13-7069 "

768

5' ATGCC TGA GTAATGTGTAGG

M13-407 "

1356

5' GAAGCAAAG CGG ATTGCA

M13-806 "

1755

5' TTA TACCA GTCAGGACGT

M13-1491 "

2440

5' AGCTTGATA CCG ATAGTTGC

M13-2566 "

3455

5' CGACAG AATCAAG TTTGCC

M13-3972 "

4921

5' AATCG CAAGACAAAGAACG

M13-5464 "

6413

5' GTATAACG TGC TTT CCTCC

Witnessed & Understood by me,

Dat

Invented by

Date

T Pag N

6/26/95

6/26/95

Record d by

6/24/95

Pamb Pamb

e No. _____

variation of primer stocks and dilutions:

Primers were made by Gibco BRL Custom Primers order# 510790A

Each primer was ~~first~~ resuspended in sterile dH₂O at a CF=100µM
- spun down, H₂O added, 2 min RT, vortex, invert

primer	nmols/tube	volume of dH ₂ O added to make CF=100µM
3-6301 anchor	54.7	547 µL
3-6681	44.09	440.9
3-7069	39.14 57.29	391.4 572.9
3-407	39.14	391.4
3-806	51.92	519.2
3-1491	69.49	694.9
3-2506	66.34	663.4
3-3972	34.68	346.8
3-5464	42.45	424.5

Each 20 µM aliquots of each primer were made from the 100 µM stocks:

⇒ 1:5 dilution, 40 µL of 100 µM primer stock

+ 160 µL dH₂O-sterile

200 µL for each primer except m13-6301 anchor

1:5 dilution of m13-6301 anchor, 200 µL of 100 µM stock

+ 800 µL dH₂O-sterile

1 mL

variation of PCR expt from p 402 µL of 1 Kb ladder and 10-18 µL of PCR products were run on a 1% TAE
agarose gel at 190 V (~180 mA, 40 W)

Gel recipe: 220 mL 1x TAE

buffer: 2 L 1x TAE

2.2g agarose

w/ 70 µL CTBr at each end

wt = 4.35.6g, boil, reweigh & add H₂O back to orig. wt
stir & cool

add 15 µL, 10 mg/mL CTBr & pour into rig. w/ quarters

To Page No. _____

d & Und rst od by m ,

Date

Inv nted by

Dat

6/26/95

R corded by

4/26/95


Ensign Lamb

Project No. _____

Book No. _____


TITLE _____

From Page N. _____

taped into
book 6/26/95 cc

> Since 10ul of the rTa
PCR rxns contained suc-
low level of product,
tried loading 15ul in or
to be able to visualize
products better

> No specific products
seen in the Tne. lat

taped into book
6/26/95 cc

> note false product
1491 primer

- only 1 band of 2.
was expected

- instead there are 2
~ 0.5-1Kb and 2-3

> Note false product
m13 2506 primer -
should be ~3.5Kb

> also a short, false
product w/ m13-5464

To Page

Witnessed & Understood by me,

A handwritten signature in cursive script, appearing to read "David P. Lamb".

Date

6/30/95

Inventor by


Recorded by

David P. Lamb

Date

6/26/95

g N _____



> 6' extension time did not
give more product. - worse
yield than with 2' extension

cc taped into
book 6/26/95

To Page No. _____

d & Understood by me,

SB Kemp

Dat

6/30/95

Inv nt d by

Recorded by

Cawley Comb

Date

6/26/95

113 PCR system: optimizing annealing temperature

N		A		(2 runs)	
LTI	10x PCR PCR buffer	120	✓		
50 mM	MgCl ₂	36	✓		
4 dNTPs	10 mM each	24	✓		
M13 RF	1 pg/μl (CN4132)	12	✓		
M13 6301 "Anchor"	20 μM	24	✓		
rTag	5 μg/μl	12	✓		
H ₂ O		948	✓		
		1.176	ml		

15	15	51, 49, 47°C	annealing temp w/ 1st 3 oligos	program 19, 94°C
				18, 94, 15"
				51, 3"
				70, 2"

A	97 μ l	→			
6671	2	2	2	2	
7065	2	2	2	2	
407	2	2	2	2	

$V_f = 100 \mu l$

14	9600	51°C	annealing temp	1-3
15		49°C		4-6
16		47°C		7-9

% gel: $\frac{2.64g \text{ agarose} + 220mL}{220mL}$ 1 - 9 100
100B 100B
ladder ladder

3.36g agarose + 220mL TAE wt = 480.3g

To Page No. _____

Read & Understood by me,

Bo/ang

Date

6/30/95

Invented by

Recorded by

Emily Pomeroy

Date

6/24/95


Project No. _____

Book No. _____

TITLE _____

om Page No. _____

Result:



6/20/95 cc

primer combo	expected product length (bp)	observed product length (bp)
anchor + 4681	380	✓ looks like 380 relative
anchor + 7069	768	✓ 768 "
anchor + 8407	1356	✓ 1356

> specific products of the right length were made, but the yield was still v. low. Cole thought lowering the annealing temp might result in a higher yield. It did not. The lanes on the right side of the gel just look a bit darker because the light box is brighter on that side (note how the 10kb ladder looks more intense on right side even though 10ul was loaded on left & right sides. Also note

> Next we'll try to increase the yield by using denaturation time (from 15" to 30" - Veri's suggestion), [Tag], 1 cycle & [primer] - in c the anchor primer has a hard time annealing due to 2° struc. If it does, then lowering annealing temp would exacerbate t problem.

T Pag N.

Witness d & Underst d by m ,

Dat

Inv nt d by

Dat

6/30/95

R c rded by

6/24/95

Parker P. Smith

M. J. Olson

Δ denaturation temp, Δ [Tag]
 Δ cycle number, Δ [primer]

Project No. _____

Block No. _____

49

e N _____

10x PCR buffer

10 mM MgCl₂

4 dNTPs 10 mM

413 RF (lot FAST01) 1 pg/μl
 opened 6-26-95

H₂O

100 ✓

30 ✓

20 ✓

10 ✓

740 ✓

900 μl

different lot than on p. 40
 dilution of stock

$\sqrt[3]{\text{stock}} \xrightarrow{2\mu\text{l}} \sqrt[3]{198\mu\text{l H}_2\text{O}} \xrightarrow{2\mu\text{l}} \sqrt[3]{198\mu\text{l H}_2\text{O}} \xrightarrow{10\mu\text{l}} \sqrt[3]{360\mu\text{l H}_2\text{O}}$
 $\frac{1}{100} \quad \frac{1}{100} \quad \frac{1}{37}$
 = 3.7×10^5 fold dilution

1 2 3 4

90 μl

5 1 4

4 8 4 8

5 6 7 8

5 1 4

4 8 4 8

1:1 mix of anchor + 407

30 μl 20 μM stock anchor

30 μl 20 μM 407

60 μl, 10 μM each

(400 or 800 nM oligo each)

2
 301 Anchor
 413-407
 (each)

5 μl 1 1 2 2

1 1 2 2

74°C, 1' initial denaturation

15" denaturation
 94°C

53°C, 30 sec

70°C 2 min

Lab 15, 9600

30 sec denaturation
 94°C

11

11

Lab 16, 9600

30 and 40 cycles 3:15 PM

520 μl rxn 1-8
 + 2.3 μl of LI

left in PCR machine, at 4°C ON

freeze ON

To Page No. _____

& Understood by me,

lump

Date

6/30/95

Invented by

Recorded by

Carolyn Cont

Date

6-26-95

Project No. _____

Results P. 49

Book No. _____

TITLE _____

50

From Page No. _____

15 sec

30 sec

Denaturation time

units ~~1000~~

5

10

5

10

primer (nM)

400

700

400

700

400

700

400

700

cycles

30

38

30

38


30

38

30

38

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16



E 1356

← primer dimer

0. tube 1 compared to tube 3 on 48 h low new M13RF only gives little or no improvement in yield at 30 cycles.
1. lower [primer] (400 nM) is best (tubes 1, 2 vs 3, 4) (where denaturation is for 30 sec, 400 vs 700 nM primer about equal)
2. Lower T_{aq} is best: only primer dimer made for 10 (# 5-8 and 13-16)
3. 38 cycles made more product than 35 cycles (tube 2 vs 1)
4. 30 sec denaturation gave less product than 15 sec.

To Page No. _____

Witnessed & Understood by m,

Date

6/20/95

Invented by

Recorded by

Paulson Pankas

Date

6-26-95

R. Pankas

ig No. _____

Therefore:

keep 15 sec denaturation

Try even less primer (eg 100 200 300 400 nm)

Try even less time (2, 3, 4, 5 min)

Try different cycle number (30 - 40)

Try M13 6301 (anchor) alone and with other primers
with no target to look for primer dimerTry more M13 RF target
eg. 0.1 pmol — 10 nmolTry R1 (or Bam, Hind III etc) "check buffer salt
cut↓
PCR with no purification

Test all primers with test conditions

To Page No. _____

ed & Understood by me,

5/30/95

Date

6/30/95

Invented by

Recorded by

C. M. J. (only)

Date

6-26-95

Project No. _____

Book No. _____

52

TITLE

1.1X PCT stability assay
 4°C storage (see P121, 9) -20°C and -70°C.
 freeze freeze Thaw

From Page No. _____

	Rxn #	array	ul	
# 10 (P121, 9): no det	1-3	2	✓	This is 5 month point for 4°C study
# 11 " 1.1X	4-6	3.64	✓	(same as P121, 9; 154, 9; 174, 9; 37, 10 = 0, 1, 2, 4 m)
untamp # 11: 1.1X	7-9	3.64	✓	(154, 9 is 0 time point, P 38, 10 is 1 month)
Tag # 125 did (same as P121, 9)	10-14	2	✓	
1.1X May 8, 1995	15-17	2	✓	called "new" on P 34, 10
1.1X field test	18-20	2	✓	called "old" P 34
1-27-95 -20°C	21-23	2	✓	Joel's took aliquot from samples on 1-27-95 gave 2 freeze Thaw stored at -20°C, 5 months at with unknown effect from freeze/thaw
-20°C 5/24/95	24-26	2	✓	from 1.1X May 8, 1995 (compare to Rxn # 15-17 above) ~1 month at -20°C with no extra freeze th
-70°C 5/24/95	27-29	2	✓	from 1.1X May 8, 1995 - its ~1 month at -
* 10 freeze Thaw	30 32	2	✓	"new" * used 1.1X May 8, 95 (above) in
20 freeze Thaw	33 35	2	✓	dry ice EtOH → 30°C bath re
30 freeze Thaw	36 38	2	✓	start with 60 µl and take out at 10, 20, 30 freeze Thaw
Fr 13-17 S200 Tnl				
P 25, 10				
1/700	39	2	✓	
1/700	40	2	✓	
1/700	41	2	✓	
The 5-7-95 (~70.7 µl on P 25)				
1/8000	42	2	✓	
	43	2	✓	
	44	2	✓	
48 µl Tox unit array mix (P120, 9) in each				
10, 74°C				
Kill' with 10 µl 0.5M EDTA, spot 20 µl on 6-FC				


Witnessed & Understood by me,



Dat

6/30/95

Inv nt d by



R c rd d by

Dat

6-28-95

To Page 1

Corrells R.J.L

Appl. No. 09/558,421

Project N _____

Book N _____

Percent 5:
of zero time for
in P122/9

e No.		u		u/A relative to log (#10-14)		Percent	
1	1597.00						
2	1752.00						
3	1760.00						
4	7709.00						
5	6150.00	} 6760 u/r		.030		94%	
6	6422.00						
7	4510.00						
8	5085.00			.025		72%	
9	5662.00						
10	8620.00						
11	9351.00	} - 8931 CPM ave					
12	8531.00			.04 (by definition)			
13	8321.00						
14	9832.00						
15	5618.00						
16	5895.00	} 5632 ave	5632	.025	.025	10.9%	
17	5384.00		5631				
18	5128.00						
19	5036.00	} 5215 ave	5215	.023 u/A	.023	10.5%	
20	5481.00		5231				
21	3989.00						
22	4058.00	} 3673 ave	its decay at -20°C for 1 month + 6 freeze thaws			54% recover.	
23	2971.00		above it came from tube #11				
24	5931.00		above (P122/9)				
25	5591.00	} 5921	conclude -20°C for activity			from .55 u/A measure	
26	6242.00		1 month -20°C from 5-8-97 (1.1X #15-17 above)			above (#4-6)	
27	5891.00						
28	5381.00	} 5712	Sample from 1.1X of 5-8-97			no activity but	
29	5865.00		(corrected to 1.1X above)			for 1 month at -20°C	
30	5644.00						
31	5407.00						
32	5271.00	} 5440	started with 1.1X			97%	
33	5362.00		(stored out at 40°C)				
34	5494.00	} 5405	(stored out at 40°C)			96%	
35	5361.00		(1.1X - 17 above)				
36	5556.00					100%	
37	6159.00	} 5638					
38	5200.00						
39	138.00						
40	287.00						
41	137.00						
42	1014.00						
43	960.00						
44	1092.00						
45	395.00						
46	110131.00						
47	110429.00						

The died at 40°C off 5200 col of 5-1PST
got 7.7 u/A on PST

aux-RKGD
= 62.7 CPM = 27.3 pmol (32.8 u/A) agrees with 36 u/A
- eg Flynn got out 8 gpt P13. May 70.7 u/A
BKGD of 25 may be wrong. will conclude LHS
2.1 mix The is stable until more data is available
2.1 mix (68.9 CPM/pmol)

To Page No. _____

I & Understood by me,

stamp

Date

6/30/95

Invented by

Recorded by

Date-30-95

6-28-95

①
17 sample
orgnt bff. formulation -

created
 ↙ ↓ ↓ Date 27/1
 40°C RT -20°C

~~RTG~~
~~NaN₃~~

1.1x
+ ~~long~~ ~~del~~ ~~del~~

The unit area

② Date _____
Field Fest _____

needed

RT (40C)
*

{ then
 make aliquots
 20 µl each
 at -20°C
 (only 1 or 2 freeze + thaw)

notes from
Joe Jones

6-27-95

May 8, 95 5:

(= "new" on P34, 10)

500ml batch

107

centered

40c/

24/5/95 ✓
ad 19 and 1-2

-20°C -70°C π 4°C
* * *

2000

3

70

• • •

25

Turnover by TFI/Vent - 7 replicates of _____ Bo k No. _____
 epicenters TFI/Vent and LTI's TFI/Vent in Epicenters SB - zero time of for stability

19 N —

pose: To establish the 3' exo activity of TFI/Vent by Turnover on gapped DNA at time zero of stability study. 7 replicates to reduce error

Background: The zero time point turnover assay has already been for Vent alone and LTI's TFI/Vent in LTI SB with 5 replicates NB11 page 17 and 21

pro to deliver 100ul stop & 100ul rxn-wipe tips

ials: mixA, enough for 40 rxns, each using 98ul of mixA ✓
 for 23 = 1467ul H₂O ✓
 $40 \times 63.786 = 2551.44 \text{ ul H}_2\text{O}$
 $40 \times 20 = 800 \text{ ul 5x Cheng}$ → 4005x Cheng ✓
 $40 \times 13.5 = 540 \text{ ul activated DNA}$ 310ul act DNA
 $40 \times 0.5 = 20 \text{ ul dATGC-TP, 10mM each}$ 11.5¹ dATGC-TP
 $40 \times 0.214 = 8.56 \text{ ul } \alpha^{32}\text{P dATP}$ 4.92ul $\alpha^{32}\text{P dATP}$
 2253.422 ^{6/30} Amc
 3920ul

5 1-7 8-14 15-21 22-28 29-35
 tube 1-21 22-42 43-63 64-84 85-91
 A 98ul — prewarm to 68°C 15-21 43-49

Vent (LTI SB) 2ul

Vent (epi SB) 2ul

1/Vent (epicenter TFI) 2ul

nt diluted (5-14.55)

enz

2ul LTI SB - Tag ³²P
 To Page No. _____

ed & Understood by m ,

Date

Invented by

Dat

Shulans

6/30/95

R corded by

Carolyn Combs

6-29-95

From Page No. _____

At 5, 10, 15 min remove ¹⁰20ul rxn to ¹⁰20ul stop soln (P.14)
 Spot 2ul on PET plates
 Spot 10ul on GFC filters (2 per rxn) + 2ul mix A 3x

tube #s

1	2	3
4	5	6
7	8	9
10	11	12
13	14	15
16	17	18
19	20	21

epicenter's Tfl enzyme + Vent in ^{epi}LTISB
 (=Tth?)

22	23	24
25	26	27
28	29	30
31	32	33
34	35	36
37	38	39
40	41	42

LTIS's Tfl eng in ~~LTIS~~ epicenters SB + Vent
 =Tth

5	10	15
43	44	

43-49 = no enzyme

GFC'S 1-6 = 2 replicates of epicenter eng + Vent
 7-12 = 2 replicates of LTIS's eng in epicenter SB
 13 = no eng
 16-18 = 2ul mix A

T Page 1

Witnessed & Understood by me,

Date

6/30/95

Invented by

Recorded by

Date

6-29-95

From Pag No. — see. p. 51

purpose: The PCR worked well (gave a large amount of product when 38 cycles were done using 400nM primer 5u/100ul Taq, 15 sec denaturat.

Now we'll try to get the same good plateau yield by optimizing [primer], [enzyme], checking for anchor probe [target] & linearizing target - using just 30 cycles

program 74 mix

expt 1. Δ [primer] from 100, 200, 300, 400 nM w/ 5u Taq, 15" denat 30 cycles

make **A**: 120ul 10x LTI PCR buffer
 36ul 50mM MgCl₂
 24ul 10mM 4dNTP
 12ul 1pg/ul m13RF → dilute w/ TC
 888ul H₂O
 1069ul
 1080ul

run	1	2	3	4
A	90	90	90	90
H ₂ O	9	8	7	6
primer mix				

A: 50ul 10x PCR buffer
 15ul 50mM MgCl₂
 10ul 10mM 4dNTP
 5ul 1pg/ul m13RF
 395ul H₂O
 5ul 5u/ul Taq
 480ul

2ul of 370ug stock + 198ul TE - mix

2ul + 198ul TE - mix

10ul + 360ul TE ⇒ 1pg/ul

make primer mix: 6ul 20uM anchor primer
 6ul 20uM 407 primer
 12ul 10uM each

To Pag 1

Witnessed & Und rsto d by m ,

Date

Invented by

Date

R corded by

g No.	Tube	1	2	3	4
0	0	1	2	3	✓
mer mix	4	3	2	1	✓
□	96	96	96	96	✓

Tag from 1 to 7 units / 100ul rxn → 400nm primers

80ul	10x PCR buffer	✓	dilute rTaq 5x	✓
24ul	50mM MgCl ₂	✓	6ul rTaq	5 ^u /ul ✓
16ul	10mM 4 dNTPs	✓	+ 24ul Tag	dilution buffer ✓
8ul	1 pgul m13RF	✓	30ul	
16ul	20uM anchor primer	✓		
16ul	20uM 407 primer	✓		
584ul	H ₂ O			
744ul				

g No.	5	6	7	8	9	10	11
3	93						✓
20	6	5	4	3	2	1	0 ✓
Tag 1 ^u /ul	1	2	3	4	5	6	7 ✓
	100ul						

target DNA, 800nm primer, 10u Tag:

3	35ul	10x PCR buffer	✓
	10.5ul	50mM MgCl ₂	✓
	7ul	10mM 4 dNTPs	✓
	242.5ul	H ₂ O	✓
	7ul	rTaq (5 ^u /ul)	✓

g No.	12	13	14
er	4ul anchor	4ul 407	4ul anchor
	4ul 407	4ul 407	primers are 20uM
2	94	4	0 ✓
	92	92	92 ✓

To Page No. _____

Read & Understood by m ,

Bolar

Date

6/30/95

Invented by

Recorded by

Cameron P. Miller

Dat

6-21-95

From Page No. _____

Δ [Target], use 400nM primers and 5u Taq 100ul
for 7 rxns

[D]: 70ul 10x PCR buffer ✓
21ul 50mM MgCl₂ ✓
14ul 10mM 4 dNTPS ✓ * dilute m13 RF to 0.5 pg/ul
14ul 20mM anchor ✓
14ul 20mM 407 ✓
462ul H₂O ✓
7ul 5u Taq 5u/ul ✓
20ul 1 pg/ul m13 RF
20ul TC ✓

Tube #	15	16	17	18	19	20
m13 RF (0.5 pg/ul)	1	1.5	2	6	10	14 ✓
H ₂ O	13	12.5	12	8	4	0 ✓
[D]	36					→ ✓

RI digestion of template:

3ul 1 pg/ul m13 RF ✓
5ul 10x PCR buffer ✓
3ul 50mM MgCl₂ ✓
13ul H₂O ✓
1ul CcoRI 10u/ul ✓
25ul

3ul ✓
5ul ✓
3ul ✓
13ul ✓
1ul H₂O ✓

37°C, 30'

37°C 30'

1 pm - 1.50 pm

8.3ul 8.3ul 8.3ul
91.7ul E →

same

21 22 23

24

25

26

To Page 1

Witnessed & Understood by me,

Dat

6/30/95

Invented by


Recorded by

Dat


6-29-95

e No. _____

4.6 ul	10x PCR buffer	x 7	=	32.2	✓
2 ul	10mM 4dNTPs		=	14 ul	✓
1 ul	M13 RF	100 ul	=	7	✓
83.1 ul	H ₂ O		=	581.7	✓
1 ul	Tag	5 ul	=	7	✓
<hr/>				<hr/>	
91.7 ul				641.9	



anchor + 407 are primers



15 16 17 18 19 20 21 22 23 24 25 26

15 16 17 18 19 20 21 22 23 24 25 26
pg target + EcoRI - EcoRI

To Page No. _____

I & Understood by me,

Date

Invented by

Date

olanges

6/30/91

Recorded by

Carolyn P. [Signature]

6-28-91

Proj ct No. _____

Book No. _____

TITLE _____

From Page No. _____

Results -

1. lag is inhibitory for 100 μ l PCR Best results at 1 and 2 units
2. yields improve with increasing target
3. both primers present is required to make primer dimers
3. its not happening substrate or one primer self annealing to another copy of that same primer.

conclude the anchor primer should be OK for most PCR's

Witnessed & Understood by me ,

[Signature]

Date

6/30/95

Investigated by

Recorded by

Date

6-29-95

T Page 1

23mer. mp19 ssDNA

Project N _____

Book No. _____

63

N (see P17, 9)

13 mp19 0.26 μ g/ μ l
79 nmol nt/ μ l
109 pmol circle/ μ l

200 μ l

21.7 pmol circle
0.594 nmol nt/ μ l
43.5 pmol primer
total

3mer 5 ng/ μ l
66 pmol primer/ μ l

66.1 μ l

266 μ l

23mer
mp19
circles
= 2

70°C, 5 min

cool at room temp 40 min

(for 60 μ l 23mer mp19)
conc in 100 μ l rxn

mp19 is 10 mM Tris pH 7.4
5 mM NaCl
0.1 mM EDTA

6 mM Tris pH 7.4
3 mM NaCl
0.06 mM EDTA

note Cheng 1X (P20, 10) =

20 mM Tris pH 9
85 mM KOAc
2% DMSO
1.05 mM MgOAc
8% glycerol

So buffer in 23mer mp19 will alter the reaction conditions a little since 60 μ l 23mer mp19 is needed per 100 μ l rxn

in future need more concentrated DNA

To Page No. _____

I & Understood by me,

Date

Invented by

Date

30 days

6/30/95

Recorded by

6-29-95

Carson Com

Results P.52

Exhibit 120
Appl. No. 09/558,421

Project N _____
Book N _____

Percent
of zero time pt
on P.122/9

g N		u	u/1 relative to 10g (#10-14)	
1	1597.00			
2	1752.00			
3	1760.00			
4	7709.00			
5	6150.00	} 6760 ave	.030	94%
6	6422.00			
7	4510.00		.025	72%
8	5085.00			
9	5662.00			
10	8620.00			
11	9351.00	} 8931 CPM ave	↑ .04 (by definition) ↓	
12	8531.00			
13	8321.00			
14	9832.00			
15	5618.00			
16	5895.00	} 5632 ave	.025	109%
17	5384.00			
18	5128.00			
19	5036.00	} 5215 ave	.023 1/4	54% recover
20	5481.00			
21	3989.00			
22	4058.00	} 3673 ave	it's been at -20°C for 5 months + 6 freeze-thaws it came from tube #11 above (tubes #4-6) conclude -20°C broad activity 1 month -20°C from 5-8-95 1.1X (#15-17 above)	from .054 µl meas above (#4-6)
23	2971.00			
24	5931.00			
25	5591.00			
26	6242.00			
27	5891.00			
28	5381.00	} 5712	Sample from 1.1X of 5-8-95 (corrected 15-17 above) ⇒	no activity lost for 1 month at -70°C
29	5865.00			
30	5644.00			
31	5407.00	} 5440	started with 1.1X (stored out at 40°C) of 5-8-95 (see #15-17 above)	conclude no loss of units for even 30 freeze-thaws
32	5271.00			
33	5362.00			
34	5494.00	} 5405		
35	5361.00			
36	5556.00			
37	6159.00	} 5658		
38	5200.00			
39	138.00			
40	287.00			
41	137.00			
42	1014.00			
43	960.00			
44	1092.00			
45	395.00			
46	110131.00			
47	110429.00			

ave - BKGD

= 6270 cpm

BKGD

2x mix

2x mix

The died at 4°C off 5200 col of 5-1895
got 7.7 u/l on 10/25

32.8 u/l

agrees with 36 u/l

leg Flynn got out & got P.13. May 70. 7 u/l

of 1.25 may be wrong

The is stable

will conclude Lys

until more data is available

68.9 CPM/ppm

To Page No. _____

Read & Understood by me,

10/25/95

Date

6/30/95

Invent d by

R corded by

Date - 30-95

6-28-95

May 8, 95 5:

(3)
("new" on P34, 10)
500ul batch
lot

3 batch:

①
17 sample
original batch formulation -

②
Date
Field test

created
4°C RT -20°C Date 27/1

created

created

* RT & NaN₃ *

1.1x
+ 1.1x
* *

8x
thaw
*

RT (40C) *

4°C
24/5/95
aliquots 1-2ul
-20°C -70°C RT 4°C
* * *

thaw
make aliquots
20ul each
at -20°C
only 1 or 2 freeze thaw

The unit assay

notes from
Joe Jones

Ryker
6-27-95

is
label

ag N — specific activity = $(109,973 \text{ cpm}) \times \frac{100 \mu\text{l rxn}}{2 \mu\text{l spot}} = 275 \frac{\text{cpm}}{\text{pmol (nt)}}$
 background = 188 cpm (5000 pmol) 4

CPM1	Turnover (pmol)	
	$\frac{(\text{cpm} - \text{background cpm})}{\text{specific activity}} \left(\frac{100}{2} \right) \left(\frac{20}{10} \right)$	
5' 774.00	$\frac{774 - 188}{275} \times 100 = 228$	213
10' 1379.00	433	
15' 2170.00	721	
5' 1031.00	307	
10' 1588.00	509	
15' 2241.00	747	
5' 893.00	256	
10' 1731.00	561	
15' 1890.00	619	
5' 788.00	218	
10' 1365.00	428	
15' 1836.00	599	
5' 752.00	205	
10' 1055.00	375	
15' 1732.00	561	
5' 636.00	163	
10' 1140.00	346	
15' 1448.00	458	
5' 854.00	242	
10' 1458.00	462	
15' 2083.00	689	
5' 829.00	233	
10' 1512.00	481	
15' 2124.00	704	
5' 980.00	288	
10' 1612.00	518	
15' 2249.00	749	
5' 1182.00	361	
10' 2028.00	669	
15' 2271.00	757	
5' 1040.00	310	
10' 1816.00	592	
15' 2521.00	848	
5' 944.00	274	
10' 1729.00	560	
15' 2032.00	671	
5' 1087.00	327	
10' 1641.00	528	
15' 2701.00	914	
5' 917.00	265	
10' 2146.00	712	
15' 2530.00	851	
5' 195.00	47	
10' 162.00	47	
15' 169.00	48	
5' 194.00	49	
	186.00	
	189.00	
	223.00	

$\bar{x}_{5'} = 229 \pm 45$ (20%) → large error, next time
 $10' = 436 \pm 86$ (20%) out out dATP + dADP
 $15' = 627 \pm 101$ (16%) spot in order to correct
 for spotting error

$\bar{x}_{5'} = 293 \pm 45$ (15%)
 $10' = 580 \pm 84$ (15%)
 $15' = 785 \pm 88$ (11%)

background for turnover
 $n = 7$ $\bar{x} = 188.3 \pm 19.8$ (~10% error)
 = 6.8 pmoles

sed & Und rstood by m ,
 Polamp

Date 6/30/95
 Invented by
 Recorded by
 Date 6/30/95
 To Pag No.

Projec. _____
 Book No. _____

TITLE Turnover of TF1/Vent on 23. mp

64

From Page No. adjust up
back to 266

23. mp/9 0.46¹⁹⁵ μ g/ μ l
 0.594 nmol at/ μ
 P. 63

Mix A

~~264~~ μ l 264 ✓

(5.4 P.xus)
 23.7
 * 29 nmol at
 * 20.2 μ g/ μ l
 65.6

5 x Chevy (P21, 10)

~~88~~ μ l 108 ✓

* see P. 55
 changed to 6
 of DNA plan

4 JMP₂ 10 mM each

2.2 μ l 2.7 ✓

(50 μ m ea)

α ³²P JADP 10 mCi/ μ l (Amersham)
 H₂O

~~1~~ μ l 1.23
~~78~~ μ l 153.5 ✓

431.2 529.2 (use 98.2/100 μ l)
 see p. 65

Reactions

stop tube *

①

②

③

④

(no enzyme)

1-8

9-16

17-24

25-32

33-40

Mix A

98 μ l

2 μ l Epicentre
 storage by

TF1/Vent = 0.18 units
 (opposite TF1
 5-16-95)

2

Vent lot #17 =
 (opened 2-24-95)

0.09 μ l/ μ l = 0.18 μ l

2

0.5 = 1 μ l

2

2 (no dilution) = 4 μ l

2

VP = 100 μ l

note: come
 to Vent all
 with gap
 on P 20
 } dilute
 opposite
 TF1 into
 buffer

67°C in 9600, remove 8 μ l to 8 μ l killing solution P.
 spot 2 μ l on PET at 10 μ l on GFC.
 * note: 20 μ g/ml 23. mp/9 is ~16% as much total DNA as
 gapped DNA in Tag unit assay (500 μ g/ml) however
 m13 is almost all ssDNA substrate while gapped DNA
 may have ~20 ssDNA gaps.

at 5 10 20 40, 60, 70, 100, 120 min

* spot DuPont on PEI

To Page

Witnessed & Understood by me,

Date

Invented by

Dat

Bolemp

6/30/95

Ryan
 Recorded by
Ryan & others

6-30-95

From Page No.

		Turnover		(-2840 ^{nm}) pmol	percent turnover	percent incorporation (100) (nmol incorporation / nmol input DNA)	
center		min					
TFI	5	1	445.00	83	10.4	24	note: if app there is no T/O at c since T/O stops incorporation stops ~ 20 min
	10	2	1588.00	672	70.8	33	
	20	3	3904.00	1066	13.6	50	
	40	4	6787.00	3352	22.4	49	
	60	5	6818.00	3368	21.3	52	
	80	6	7009.00	3468	21.4	54	
	100	7	7002.00	3462	22.4	50	
	120	8	7164.00	3546	21.4	55	
vent	5	9	333.00	25	—	5	note turnover low and ends in for TFI/vent. conclude most are not full since low per incorp.
	10	10	350.00	34	—	5	
	20	11	381.00	50	—	0.9	
	40	12	521.00	122	37	1.6	
	60	13	832.00	282	43	2.2	
	80	14	1097.00	418	45	2.7	
	100	15	1474.00	613	45	3.2	
	120	16	1928.00	847	53	3.2	
1 unit	5	17	296.00	16	2	1.0	about 25 max incorp 13.8% / 55 (55% seen at 120 min) 3% still miss are not full
	10	18	508.00	115	16	2.6	
	20	1	1322.00	535	25	6.8	
	40	2	3242.00	1524	32	13.8	
	60	3	6116.00	3406	42	17.6	
	80	4	8505.00	4257	45	22.1	
	100	5	11510.00	5787	52	22.4	
	120	6	13872.00	7004	52	27.3	
4 units	5	7	750.00	—	—	—	more conclusions T/O is occurring at an eff high rate (~20-22%) for TFI considering ~30:1 TFI/vent (in units) it appears that substrates after 3' ends accumulate during time. probably at hairpins (on M13) rather than at full length (since T/O stops shortly after incorporation stops even for 1 unit)
	10	8	1728.00	—	—	—	
	20	9	4617.00	—	—	—	
	40	10	9106.00	—	—	—	
	60	11	11531.00	—	—	—	
	80	12	12228.00	—	—	—	
	100	13	12432.00	—	—	—	
	120	14	11890.00	—	—	—	
no enzyme	5	15	269.00	—	—	—	2% are JAMP background
	10	16	278.00	—	—	—	
	20	17	240.00	—	—	—	
	40	18	276.00	—	—	—	
	60	19	293.00	—	—	—	
	80	20	274.00	—	—	—	
	100	21	307.00	—	—	—	
	120	22	331.00	—	—	—	

Incubation plateau for Vent alone (~27% incorp) is lower than for TFI/Vent (~55%)
 at vent alone made at hairpins?

Witnessed & Understood by

Date

Invnted by

Date

T Page N

D. Polansky

7/7/95

Reviewed by

Sawyer Pank

7-1-55

Results of P64

Project No. (23.7 nmol nt input ss M13 DNA)
 Bo k No. Percent substrate
 copied

67

No.	Incorporation	p mol	
23	632.00		
24	37364.00		
25	37879.00		
26	106114.00		
27	104061.00		
28	34.00		
29	57217.00	5298	24
30	76607.00	7098	24
31	114434.00	11797	24
32	112213.00	11587	24
33	120505.00	12425	24
34	123172.00	12698	24
35	115839.00	11942	24
36	125999.00	12890	24
37	1441.00	98	
38	600.00	(98)	
39	1390.00	90	-512
40	2498.00	205	BKGD
41	4176.00	378	
42	5526.00	517	
43	6777.00	646	
44	7820.00	753	
45	2449.00	252	
46	6061.00	624	
47	15599.00	1628	
48	31777.00	3276	
49	40517.00	477	
50	50876.00	5245	
51	52552.00	547	
52	62653.00	649	
53	36563.00		
54#38	260.00		
55	83243.00	8582	
56	58704.00	6052	
57	82715.00	8527	
58	73558.00	7585	
59	69056.00	7119	
60	65008.00	6701	
61	502.00		
62	591.00		
63	450.00		
64	960.00		
65	307.00		
66	563.00		
67	463.00		
68	80862.00		
69	75044.00		
70	81620.00		

ave BKGD
= 512

ave
79175

194 CPM/pmol nt

To Page No.

Read & Understood by me,

Date

Invented by

Date

okays

7/7/95

Recorded by

7-1-95

Carolyn L. Combs

213 PCR system: more optimization of [primer 407], annealing temp. and [target]

g N _____

s A, enough for 28 rxns:			280ul 10x PCR buffer ✓	template primer dilution:
			84ul 50mM MgCl ₂ ✓	to 500pg = 1ul stock m13 RF
15	53°	10:50 am	56ul 10mM dNTPs ✓	370 ug/mL
16	55°	11:02 - 1:15	2094.4 ul H ₂ O ✓	+ 99ul TC
SGI	57'	=	5.6ul rTaq 5uL ✓	mix ✓ ✓
131	21:10:10		2520ul	3.1ul + 7.3ul TC ✓
				→ 7.4ul of 500ug/mL

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18
 3 RF 135 1-3 4-6 7-9 10-12 13-15 16-18
 3 ul of 10 PG/ul 9 ul 10 PG/ul 6 ul of 50 PG/ul 40 ul 10 PG
 mer mix 12 ul 18 24 12.8 ul 24 12.8 24 ul ✓
 um anchor
 um 407
 anchor 407
 0 15 ul 9 ul 3 ul 9, 3 3 ul 0 12 ul 6 ul 0 ul ✓

A 270ul → 270ul → 270ul →
300ul →
 divide into 3 tubes, 1 for each annealing temp.
a, b, c

<p>= 53° 10 pg M13RF 55° w/ 400nM primers 57°</p>	<p>4 a 53° 30 pg M13 b 55° 400nM primers c 57°</p>	<p>7 a 53° 100 pg M13 b 55° 400nM p c 57°</p>
<p>= 53° 10 pg M13RF 55° w/ 400nM primer 57°</p>	<p>5 a 53° 30 pg M13 b 55° 400nM p c 57°</p>	<p>8 a 53° 100 pg M13 b 55° 600nM p c 57°</p>
<p>53° 10 pg M13RF 55° w/ 800nM primers 57°</p>	<p>6 a 53° 30 pg M13 b 55° 800nM p c 57°</p>	<p>9 a 53° 100 pg M13 b 55° 800nM primers c 57°</p>

T Pag No. _____

T Pag No._____

and Understood by me,

Date _____

Inv nted by

Date _____

Polamp

7/7/0

R garded by

7/5/95

Recorded by Evelyn Combs

From Page No. _____

1.7% agarose gels:





6 g agarose
600 mL 1x TAE46 μ L 10 mg/mL EtBr

M13RF: 10 pg 30 pg 100 pg

primer (nm) 400 600 800

cycles

30 35 40 45 50 55 60 65 70 75 80 85 90 95 100



Conclusion

1. no advantage to >400 nm primer
2. more target improves yield and specificity
3. 57°C is most specific


Carlynn Combs
7/4/95

Page N _____


og target

primer (nm)
 cycle #


10			30			100		
400	600	800	400	600	800	400	600	800
130 35	130 35	130 35	130 35	130 35	130 35	130 35	130 35	130 35



anneal
 temp
 57°C



55°C



55°C

To Page No. _____

ed & Understo d by me,

Solamp

Date

7/7/95

Invented by

Recorded by

Date

7-6-95

From Page No. _____

title: m13 PCR system: 4 annealing temperatures and 4 primer sets

53°C, 55°C, 57°C, 59°C

anchor + 6681
7069
407
806

purpose: We have established optimal [template] & [primers] for anchor + 407. Now we'll optimize annealing temp and cycle ~~x~~ for the other primers 6681, 7069, 806. ~~we do it~~ assuming that these other primers will work well with the [template] & [primer] that worked well w/ 407. Later, we'll titrate [Taq] & [Tne] for these primers.

background: • 59°C has not been tried before w/ any primers

- in an earlier expt 57°C worked best for 407
- [template] = 100 pg / 100 μ l rxn - found to be best for 407
- [primer] = 400 nM - found to be best for 407
- [Taq] = 1 μ / 100 μ l rxn, more was inhibitory for 407
2 μ / 100 μ l rxn will be tried for 806 which makes the longest product

• expected product sizes:

anchor + 6681 \rightarrow 380 bp

+ 7069 \rightarrow 768 bp

+ ~~1356~~⁴⁰⁷ \rightarrow 1356 bp

+ 806 \rightarrow 1755 bp

Cowling 7/6/95
Cant

materials: m13RF 50 pg/ μ l, diluted or: ~~etc~~ 7/6/95 NBII p.

Witnessed & Understood by m,

DBLays

Date

7/7/95

Invented by

R c rded by

Date

7/7/95
7/6/95

To Page

procedure:

make a master mix for 21, 100ul rxns - containing everything but the primers:

(A) { 210ul 10x PCR buffer
1654.8ul H₂O
43ul 50mM MgCl₂ Cf = 1.5mM
42ul 10mM dNTP's Cf = 200uM
42ul m13mp19 RF 50pg/ul stock
42ul anchor primer, 20uM Cf = 100pg/100ul rxn
stock Cf = 400nM
4.2ul Taq 5u/ul Cf = 1u/100ul rxn
2058ul

remove 441ul and add 0.9ul Taq (5u/ul) - for 806 primer
w/ 2u Taq/100ul rxn

s (1.5mL)	1	2	3	4	5
ner, 20uM stock	8ul 6681	8ul 7069	8ul 407	8ul 806	8ul 806
	392ul				(B) →
	<u>400ul</u>				

divide into 4, 100ul aliquots in 9600 PCR tubes
and put each tube in different 9600's set to different
annealing temperatures 53° - Lab 15

55°C - Lab 16

57°C - SGI

Tammy 7/6/95
Combs

59°C - Lab 14 * note that 30 cycle

aliquots were taken

during ramp to 94°C

pause the 9600's during later part of the 70°C, 2'
extension to withdraw 25ul samples at 25, 30, 35 cycles.
+ 2.8ul Blue Juice

Read & Understood by m ,

Pokamp

Date

7/7/95

Invented by

Recorded by

Combs

Date


7/6/95

primer
annealing temp


6681	7069	407	1u Taq 806	2u Taq 806
53° 55° 57° 59°	53° 55° 57° 59°			

1/2


25 cycles



30 cycles




35 cycles



Note: 59°
cycles - 1
9600 during
to 94°C. 50
product is,
ssDNA - acc
for the 2.
in 7069 & 407

Dawling
Combs
7/6/95



25 cycles

30 cycles

35 cycles

- 59°C annealing temp for 25 cycles gave a good product yield and the least nonspecific products for primers 6681, 7069, 407
- Try higher annealing temps for primer 806 because product is beginning to come up at 59°C

Witnessed & Und rsto d by me,

Dawling

Date

7/7/95

Invent d by

Dawling

Recorded by

Dawling

Date

7/6/95

m13 PCR system: titration of [Tne] and [Tag] with 6681, 7069, & 407 primers - 3 cycle x 5

g N _____

- purpose:
- 1) To optimize the [Tag] → 100 uM rxn was previously found to be inhibitory, so lower conc. will be tried - 0.25 uM.
 - 2) To see if Tne can synthesize any of the products expected w/ 6681, 7069, & 407 primers + anchor primer - no product was made (p. 44) when 50 Tne/100 uM was tried. The [template] and [primers] that were optimal for Tag will be used in the Tne PCR rxns. 0.1 - 1 uM will be tried.
- 25, 30, 35 cycle samples will be taken
 - 59°C annealing temperature, 400 nM primers, 100 pg/100 uM rxn m13 target

materials:

mix [A] w/ 6681, for 11.5 rxns


885.5 uL H ₂ O	
115 uL 10x PCR buffer	
34.5 uL 50 mM MgCl ₂	Cf = 1.5 mM ✓
23 uL 10 mM dNTPs	Cf = 200 uM ✓
23 uL 20 uM anchor	Cf = 400 nM ✓
23 uL 20 uM 6681 primer	each ✓
23 uL m13RF 50 pg/uL	Cf = 100 pg/rxn ✓
112.7 uL	

see p. 69 for dilution of stock to 50 pg/uL lot FAS 701

mix [B] same as [A] but use 23 uL, 20 uM 7069 primer ✓

mix [C] same as [A] but use 23 uL, 20 uM 407 primer ✓

the dilutions in Tag SB: RT, rinse pipet, vortex 25 sec



Also do same dilutions w/ Tag in Tag SB - only 1 dil will be used

top soln = 100 mM EDTA → 80 uL 0.5 M EDTA - needed to kill exo
 5x Blue Juice + 320 uL 10x Blue Juice ice won't kill it
 400 uL

To Page No. _____

d & Understood by me,

Polamp

Date

7/7/95

Invent by

Recorded by

Date

7/6/95

Project No. _____

Book No. _____

TITLE _____

76

From Page No. _____

run #	1	2	3	4	5	6	7	8 → 14	15-21
	1-98ul [A]			1			1-98ul [B]	1-98ul [C]	1
1.4% Tne	2ul						same as 1-7	same as 1-7	same as 1-7
0.75% Tne	2ul								
0.5% Tne		2ul							
0.375% Tne			2ul						
0.25% Tne				2ul					
0.125% Tne					2ul				
0.05% Tne						2ul			
	100ul								

* remove 25ul rxn after 25, 30, 35 cycles to 2.8ul Blue Juice

note: 59°C annealing temp, 10:40 am - 2
 94°C 15" Lab 15
 59°C 30" program 74
 70°C 2' extension 9600
 4°C final

run #	22	23	24	25	26-29	30-33
	1-98ul [A]			1	1-98ul [B]	1-98ul [C]
0.5% Tag	2ul				2	2
0.375% Tag		2ul			2	2
0.15% Tag			2ul		2	2
0.125% Tag				2ul	2	2
0.05% Tag					2	2

30 cycles only

* * remove 25ul to 2.8ul Blue Juice

* remove 25ul of Tne rxns after 25, 30, 35 cycles during last part of 2min, 70°C elongation into 4ul stop soln in microtiter plate. The final [CDTA] = 10mM
 [Blue Juice] = 1.1X

* * remove 25ul of Tag rxns after 30 cycles + 4ul stop soln

To Page N

Witnessed & Understood by m ,

M Polamp

Date

7/7/95

Invent d by

R. J. ...

R c rd d by


Randy Pombi

Date

7/7/95
7/6/95

Tag N _____ Results

expected product size	380	768bp	1,356bp
primer	0681	7064	407
Tag units	12.5, 1, 75.5, 25.1	some	some



Tag units 1, 75.5, 25.1


1) Tag made the expected specific products.
1 u/100ul rxn Tag was optimal

2) The didn't make any specific products -

← 35 cycles only primer-dimers more The gave more primer dimers
← 30 cycles - try primer extension exp - low pol rate?
← 25 cycles The

← Tag 30 cycles

longer exposure of same gel



To Page N _____

ed & Understood by me,

Dooley

Date

7/13/95

Invented by

Recorded by

Lawson Lamb

Date

7/17/95

Project No. _____

Exhibit 127

Book No. _____

TITLE _____

Appl. No. 09/558,421

78

From Page No. _____

	GAPDH	PCR	
5 x Chemo (w/ dNTPs) (P2, 10)	[A] 100 ✓	[B] 100 ✓	10 R
10 mM dNTPs	10 ✓	10 ✓	(200 μM)
Human spleen genomic DNA (HS #2 19/4) 80 ng/μl	12.5 ✓		(100 ng/5)
GAPDH(+) 2112, 10 μM	20 ✓		Cf (400 nm)
GAPDH(-) 2113 10 μM	20 ✓		Cf mg for
Mg OAC 100 mM	—	2.5 ✓	Cf=1. for [B]
PCR DNA xmn I 25 pg/μl		10 ✓	
2836, 10 μM		20 ✓	
2837, 10 μM		20 ✓	
H ₂ O	317.5 ✓ 480	317.5 ✓ 480	

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18
[A]	48	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—
[B]	—	—	—	—	—	—	—	—	—	48	—	—	—	—	—	—	—	—
The 5-7-95 (Liz)	0.5	2								2								
	1.5		2								2							
	2			2								2						
					2								2					
						2.5								2				
The 5-9-95 Adam 3-2-95							1	2							1	2		
ditto 2 μl																		
Tag 2 μl								1	2									

all dity
made with
Tag
straw
buffer

Witness d & Understood by m ,

[Signature]

Dat

7/14/95

Inv nted by

[Signature]

R cord d by

[Signature]

Dat

7-10-95



T Page N

Sheet3 Chart 1

Sheet3 Chart 1


LT'S TIVent IN EPICENTER'S SB: TURNOVER AT TIME ZERO IN STABILITY STUDY (6/28/95)

LT'S TIVent IN EPICENTER'S SB: TURNOVER AT TIME ZERO IN STABILITY STUDY (6/28/95)



Sheet4 Chart 1

LT'S TIVent: POLYMERIZATION AT ZERO TIME POINT IN STABILITY STUDY (6/28/95)



calculations:

To Pag No. _____

ed & Und rstood by me,

Dat

Invented by

Date

Rec rded by

7/11/95


Carolyn C. Smith

9 N - purpose: Can we make the 380bp m13 PCR product w/ Tne? Earlier p. 77 w Tne failed to make this product in 10x PCR buffer.

4 for rxns, 100ul per rxn: 340 300ul 5x Cheng buffer
 34' 30ul 10mM dNTPs $C_f = 200\mu M$
 34' 30ul 50 99/ul m13mp19 RF in TC $C_f = 100\mu g$
 34' 30ul 20uM anchor primer
 34' 30ul 20uM 6681 primer $C_f = 400nM$
 1190' ul H_2O

1066ul + 470ul

ans of Tag + Tne (5/7/95 Liz) in Tag storage buffer (SB)



1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
93ul A								2ul Tne (5-7) 95 Liz							
2ul Tag								2							
	2								2						
		2								2					
			2								2				
				2								2			
					2								2		
						2								2	
							2								2
100ul															

15', 1 min 94°C
 94°C 30" denat.
 55°C annealing temp 30"
 72°C exten 2"
 remove 25ul aliquots after 25, 30, 35 cycles
 program 76
 1pm - 4pm
 * remove 25ul p. 79
 + 3ul STOP soln (w EDTA)
 * run 25ul

To Page No. _____

I & Understood by me,

Polamp

Date

7/14/95

Inv nted by

[Signature]

Recorded by

[Signature]

Date

7/11/95

Project No. _____

Book No. _____

TITLE M13 PCR: 3 small products w/ Tne

84

From Page No. _____

purpose: To determine if Tne can make the 3 smallest m13 PCR products 380 bp, 768, 1356 bp, 380t using the conditions which worked for gapDH & PUC p.80 - Cheng buffer, 400 nM p, 35 cycles { 55°C annealing temp, 94°C 30" denat, 100 p9/10, 72°C 2' extension

If the products are made, we can use these conditions for m13 primer extension experiments

materials: mix for 15, 50ul rxns

[A] 150ul 5x Cheng
510ul H₂O
15ul 10mM dNTPs
15ul 50p9/10 M13mp19 RF
15ul 20uM anchor primer
705ul

[B] = 235ul A
+ 15ul 20uM 6681

[C] = 235ul A
+ 15ul 20uM 7069

[D] = 235ul A
+ 15ul 20uM 407

dilutions of Tne (5/7/95 Liz) in Tag storage buffer:

50ul stock
+ 23.3ul Tag5B
33.3ul
1:3.33x
1.5 u/ul

20ul
+ 10ul Tag5B
30ul
1:1.5
1 u/ul

20ul
+ 20ul Tag5B
40ul
1:2
0.5 u/ul

20ul
+ 20ul Tag5B
40ul
1:2
0.25 u/ul

20ul
+ 20ul Tag5B
40ul
1:2
0.125 u/ul

rxn	1	2	3	4	5	6	7	8	9	10	11	12	13
[B]	48	48	48ul										
[C]				48ul									
[D]									48ul				
1.5 u/ul Tne				2					2				
1 u/ul	2				2				2				
0.5 u/ul		2				2				2			
0.25 u/ul			2				2				2		
0.125 u/ul								2				2	

→ start rxns on ice w/ en1

Result on p.86
Only the 380t product was made.

Page N

With ss d & Underst d by me,

[Signature]

Date

7/14/95

Invent d by

[Signature]
Rec'd by
Dawlan Pomb

Dat

7/12/95

³²P 23. MIS for TFI/vent

Project No. _____
Book No. _____

Exhibit 131
Appl. No. 09/558,421

87

eN — small experiment at P. 64 except ³²P 23 instead of ³²P JAI
increasing T/O with time (P. 64-68) due to strong pause sites
visible on an agarose gel?

5 ng / λ 23mer
5x Kinase buffer
PNK
8 ATP 10 mg / μ l 7-14-95

35 μ l
10.2 μ l
1.2 μ l
5

1.75 ng
23 pmol 23mer

51.4 μ l

37°C, 30 min \rightarrow 55°C, 5 min

51.4 μ l
42.6
10.6

113 mp19 0.26 μ g / λ
0.79 nmol nt / λ
109 pmol circle / λ

11.5 pmol
circle total

200

$$\frac{23 \text{ mer}}{\text{circle}} = 2$$

70°C, 5' cool slowly

$f = 0.41$ nmol nt MIS / μ l

it will be ~237 total nmol nt / 100 μ l Rxn
same as for P. 64

To Page No. _____

Read & Understood by me,

Polansky

Date

7/14/95

Invented by

Recorded by

Date

7-15-95
RZ

From Pag No. _____

(A)

3.5 R_x100

32P 23mm - mpla (P87)

171.1

0.41 nmol it / μ l0.566 pmol acid / λ

5X Cherry (no dNTPs)

P21, 10

70

10 mM dNTP₂H₂O

1.75

100.15

C_f = 50 μ mV_p = 343(use 98 μ l / 100 μ l R.)

(1)

(2)

(3)

98 μ l

7

(A)

TFI/Vent
(specimen TFI)

5-16-95

2

(1.88 μ l total in
its only 0.94
specimen units) due
with Vent by Nir

TFI lot 31010A-502

1 μ l / λ (specimen units)1.88 μ l

(1.88 unit lot)

Vent 2 μ l / μ l

lot #17

opened 2-24-95

2

(4 units) at 1.88 is
more Vent than in TFI,
in order to get full long
products)68°C. Remove 8 μ l at 1, 2, 5, 10, 15, 20, 40, 60, 90
to 1 μ l 10X 'blue juice' in E.

Run on agarose (same as P56, 7)

To Pag N

Witnessed & Understood by me,

D. Polansky

Dat

7/14/95

Inv nted by

R c rd d by

Dat

7-14-95

Project No. _____

Exhibit 132

Book No. _____

TITLE _____

Appl. No. 09/558,421

om Page No. _____


1% agarose gel

top empty

9-16	9-16	9-16
25	30	35

 The 9-16 decreasing eny

1-8	1-8	1-8
25	30	35

 Tag 1-8 decreasing eny - run into gel 1st
The from p. 80

same photo as p.
taken from a further
distance to capture
the ~~at~~ bottom of
gel

To Page N .

Inspected & Understood by m ,
GOD OlampDate
7/14/95

Inv nted by

R c rded by

Pawlen P. Smith

Date

7/14/95

Page No. _____

heng buffer, 380bp product expected


25 cycles

30 cycles

35 cycles

units

6 5 4 3 2 1 1.5 2.5 6 5 4 3 2 1 1.5 2.5 6 5 4 3 2 1 1.5 2.5 / 100ul rxns



Tne

rTaq

- 1 unit Tne is optimal, 3.5 cycles 21.5
- 6 units rTaq gave the most product for all cycle #3

To Page No. _____

ed & Understood by me,

Polamp

Date

7/14/95

Invented by

Recorded by

[Signature]
[Signature]

Dat

7/14/95

213 PCR: A PCR buffer into Cheng

No.		purpose:
1	3g agarose. 4 mL 50X TAE. 294 mL H ₂ O 300 mL wt 214.5g w/ magnet + 20 mL EtBr	To find out what component(s) of Cheng are important for making the broad smear (vs narrow, low mult smear made w/ PCR buffer)

a/s:	A for 35 rxns = 35 μ L 200 mM anchor primer 35 μ L 200 mM 6681 primer 35 μ L 200 mM dNTP's 910 μ L H ₂ O 1015	B w/ target 507.5 μ L of A + 17.5 μ L 500g/mL m13 RF
------	--	--

*Nens 2/10/95

D	= mix of 200 mM Tricine + 10.05 mM MgOAc to add KOAc to = 10X working stock 200 μ L of 1M Tricine pH 9 (from Nens) + 10.05 μ L of 1M MgOAc + 789.95 μ L H ₂ O 1 mL	C no target DNA 507.5 μ L of A + 17.5 μ L H ₂ O
---	---	--

E	= mix of 166.6 mM Tricine pH 9, 708 mM KOAc, 8.75 mM MgOAc is 8.33X stock 166.6 μ L 1M Trine pH 9 354 μ L 2M KOAc 8.75 μ L 1M MgOAc 470.65 μ L H ₂ O
---	--

Tne (5-7-95) Liz 5 μ L diluted to 0.2 μ L w/ Taq storage buffer + 148.8 μ L Tne + 6.2 μ L Tne

xx:	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	→ 28	(12+14)					
refer	30ul B target														30ul C no target							
cl ₂	5															same as 15-28 without target						
	1.5																					
		5	5	5																		
cl		10																				
KOAc		2.35 5.88																				
			6	6	6	6	6	6	6	6	6	6	6	6								
erol			2	5	8															8	8	8
SO	8.5	0	5.75 1 .5 .75 1																			
	8.5	0	7.65	4.12	9	7	4	8.5	8.5	8.25	0.5	0.5	0.5	0.5	0.5	0.5	0					
ie	5																					
	50ul																					

To Page No.

To Page No. _____

Read & Understood by me,

2/10/95

Date

7/11/95

Invented by

Recorded by

Date

7/11/95

From Page No. _____

Buffer components

pH buffer

mM

chem

Tris pH 8.4 20

KCl 50

MgCl₂ 1.5

Tricine pH 9

MgOAc



KOAc

glycerol

DMSO

Lane x

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 No ladder



anchor primer alone
6681 primer alone

- 1) The transition from "primer dimer" (smallest band) to smear occurred when Tricine pH 9 + MgOAc 1.5 mM was substituted for Tris pH 8.4 + MgCl₂
- 2) Product was made when 85 mM KOAc was added to the Tricine MgOAc buffer - glycerol + DMSO were not required to make product
- 3) Less glycerol + DMSO resulted in more product synthesis
- 4) The smears are present in the plus and minus target lanes
- 5) The "primer dimer" is made when either primer is present alone or if both primers are present as long as PCR buffer is used

To Page 1

Witnessed & Understood by m ,

D. Olave

Date


7/14/95

Invented by

Recorded by

Date

7/14/95



Expt on p. 8

C. Comb 7/14/95

C. Comb 7/14/95

+m13
template + 2 primers
for 380bp product
*Note the enzyme
brings glycerol in
rxn. 50% glycerol
7 10, 13%

-m13
template + 2 primers
no template, only

Project No. _____

Book No. _____

TITLE _____

is 32P primer incorporated into either
the "primer dimer" or smear see
Tne on P. 76

From Page N _____

(- Target)

(+ Target)

(PSI MM RNA)

	(A)	(B)	(C)	(D) RNA
1 M Tris pH 8.4	10			20 mM
3 M KCl	8.33	8.33	-	(50 mM)
50 mM MgCl ₂	15			(1.5 mM)
1 M Tris pH 9		10	10	20 mM
MgOAC 50 mM		10.5	10.5	(1.05 mM)
Mi3 RF 50 pg/ul	-	10	10	(50 pg/50 µl PCR)
10 mM dNTPs	10	10	10	
H ₂ O	366.7	361.2	350	547.25
	410	410	-410	

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17

(A)	41															
(B)																
(C)																
6681 10 µM	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	
6681 10 µM		2		2		2		2		2		2		2		
301 10 µM	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	
16301 10 µM	2			2			2			2			2		2	

Cf = 40 on x

Tne 5-75 (L2) 5 µl

0.2 µl diluted

in Tag SB

35 cycles as per P 76

note for 3 µl
8% PAGE it's
0.08 µl 32P primer

mix 20 µl PCR with 2 µl BT, 100 mM EDTA, load 10 µl on 8% AG
as per P 76 9 150 V start ~ 2 PM (?)
moves as per AK book 7 start 4:30

To remain 30 µl PCR add 1.5 µl cycle neg stop, load 5 µl
on 8% PAGE (wells 1-17 are PCR reactions)

for sequencing Rxn use 32P 6681 - mp19 ss DNA (P 71) as per
(P 27, 4)

To Page 1

With ss d & Understood by m ,

Date

7/14/95

Inv nted by

R c rd d by

Dat

7-14-95

S Polans

^{32}P 6681:

and ^{32}P 6301

Project No. _____

Book No. _____

91

No. _____	(Can see P150, 9)				
	①	②			
100 μM 6301 (anchored)	1.57		✓	$C_f = 10$	$\mu\text{M CF}$
100 μM 6681		1.57	✓		
^{32}P γ ATP 10 $\mu\text{Ci}/\mu\text{l}$	10	10	✓		
$\gamma = 7-14.95$					
PNK 1 μl	1	1	✓		
5X Kinase buffer	3.1	3.1	✓		
	15.66 μl	15.6			
1 μM 6301 cold	78.3				
1 μM 6681 cold		78.3			
	93.98				

C_f at 1X Kinase buffer
35 mM Tris pH 7.6
50 mM KCl
5 mM MgCl₂

2 dilute hot primer with 5 parts cold primer
so contribution to PCR of P90 is:
Primer in PCR

MgCl ₂	0.83	0.033
Tris pH 7.6	5.5 mM	0.23 mM
KCl	8.3	0.33

^{32}P 6681 = mp19 for sequencing

mp19 0.26 $\mu\text{g}/\mu\text{l}$ 0.109 pmol ends/ μl 10 μl (1.09 pmol ends)

^{32}P 6681 10 μM primer 1 10 pmol primer
70°C cool slow

for 3 μl on PAGE is ~0.04 μl primer

To Page No. _____

d & Understood by me,

Date

Invented by

Date

Polaris

7/14/07

Recorded by

7-14-15

1.5 mM

MgCl₂

1.05

MgOAc


Tris

pH 8.4

Tricine

pH 9


KOAc 85 mM



PCR

Tne

³²P primers



R. J. M.

PAGE 93 OF NOTEBOOK WAS BLANK

From Page No. — 7/14/95 purpose: to determine which buffer component(s) were resp

1) for the transition seen from "primer dimer" to "smear" on p. 86 1a.

The "primer-dimer" condition was 20mm Tris 8.4

50mm KCl

1.5mm MgCl₂

The "smear" condition was

20mm Tricine pH 9

50mm KCl

1.05mm MgOAc

what comp
is causing
transition(Tris^{8.4} vs
[Mg]
MgOAc)2) Is the "smear" or "primer dimer" due to DNA contamination
materials: of the prep? runs w/o D.80.0ul 0.2^u/ul Tne (5-7-95 Lig) → 7.2ul Tne stock 5^u/ul
+ 172.8ul Tag storage buffer

! 100mm Tris 8.5 * note this is a change from the original 8.4 "dimer" condition

100ul 1M Tris 8.5 (RL)

900ul H₂O

100mm Tricine pH 9

100ul 1M Tricine pH 9 (Nen 2/10/95)

900ul H₂O5.25
10.5mm MgCl₂210ul 50mm MgCl₂ - 105ul
790ul H₂O - 895ul7.5
15mm MgCl₂300ul 50mm MgCl₂ - 150ul
700ul H₂O - 850ul5.25
10.5mm MgOAc10.5ul 1M MgOAc vs 5.25 - 10.5ul 5.25
989.5ul H₂O 1989.5ul 994.77.5
15mm MgOAc15ul 1M MgOAc - 15ul 7.5ul
985ul H₂O - 1985ul 992.5ul
2mL

e N _____

25: A = 36ul 10mM dNTPs ✓
 29.88ul 3M KCl ✓
 726.12ul H₂O (1st) ✓
 792ul

[B] template + primers

374ul [A] ✓
 17ul 20uM anchor primer ✓
 17ul 20uM 6681 primer ✓
 17ul 50 pg/ul m13RF ✓
 425ul

[C] no DNA, only enzyme

374ul [A] ✓
 51ul H₂O ✓
 425ul

use 25ul per 50ul rxn → This will give

up each rxn in duplicate

400nM primers, pH 4
 50 pg template
 50mM KCl
 200 uM dNTPs
 in 50ul rxns

components:

	small smear condition = per buffer	tricine instead of Tris	1.05 Mg ²⁺ instead of 1.5	MgOAc instead of MgCl ₂	change (Mg ²⁺) standard	broad smear condition	Tricine w/ high Mg ²⁺	Tricine w/ MgCl ₂
rxn #	1,2	3,4	5,6	7,8	9,10	11,12	13,14	15,16
Tris pH 8.4	✓		✓	✓	✓			
KCl	✓	✓	✓	✓	✓	✓	✓	✓
MgCl ₂	✓	✓						
Tricine pH 9		✓				✓	✓	✓
n MgCl ₂			✓					✓
MgOAc				✓			✓	
m MgOAc					✓	✓		

To Page No. _____

& Understood by me,

Date

Invented by

Date

Recorded by

D. Kump

7/24/95

[Signature]

7/14/95

[Signature]

Project No. _____

Book No. _____

TITLE _____

96

From Page No. _____

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
100mM Tris	8.5	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10
room M Tris	9	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10
7.5	15mM MgCl	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10
5.25	5mM MgCl	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10
7.5	15mM MgOAc	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10
5.25	5mM MgOAc	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10
mix	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul	25ul
= template + primers + dNTPs + KCl																																
= H ₂ O + dNTPs + KCl																																
"No DNA"																																
0.24ul Tris	5ul brings in ~5% glycerol																															
(5-7-95 Lig)	add enzy, vortex, quick spin (4°C) to start rxns																															
	50ul then cycle in Lab 15 9600 program 76																															
	1 min 94°C																															
	30 sec 94°C																															
	30 sec 55°C																															
	2 min 72°C																															
	4°C																															
	35 cycles																															
	STOP rxns w/ 7ul 10x STOP = 100mM EDTA p. 79																															
	8x blue juice																															
	Run 25ul on 1% agarose gel																															
T	Pag																															

With ssed & Und rstood by m ,

Dat

7/24/95

Inv nt d by

R c rd d by

Dannan Ponds

Dat

7/14/95

T Pag

N _____

Tris w/ 1.5 mM MgCl₂
 Tricine w/ 1.5 mM MgCl₂
Tris w/ 1.05 mM MgCl₂
 Tris w/ 1.05 mM MgOAc
 Tris w/ 1.05 mM MgOAc
 Tricine w/ 1.05 mM MgOAc
 Tricine w/ 1.05 mM MgOAc
 Tricine w/ 1.05 mM MgCl₂



template + primers
 no DNA
 ie no template
 no primers

ecomb

5/14/95 cc

all at 50 mM KCl

lower concentration of Mg²⁺ in Cheng vs PCR buffer is responsible for the broader smear. The Tris / tricine pH difference and MgCl₂ vs MgOAc do not affect the transition from small^{narrow} to broad smear.

is model: distance between primers
 low [Mg]
 high [K⁺OAc]
 low [Tne]
 short = small narrow smear
 long = broad smear
 ✓ by decreasing primer annealing
 ✓ by inhibiting Tne from binding primers

To Page No. _____

& Understood by me,

Date

Invented by

Date

7/21/95

Recorded by

7/14/95

Emman Eomb

claw

15 9600 - 94°C 1min
 94°C 30sec
 55°C 30sec } 35 cycles
 72°C 2min
 4°C hold

rxns by adding 7.5 ul of Stop soln = 2x Blue Juice Cf = 1x
 100mM EDTA Cf = 12.5mM EDTA

Stop soln made for future use = 50% glycerol
 100mM EDTA (10x)
 0.6x TAE
 BpB

for 5mL: 2.5mL 100% glycerol
 1mL 0.5M EDTA
 1.5mL 1x TAE = 30ul 50x TAE stock + 1470ul H₂O
 + pinch of BpB

To Page No. _____

& Understood by m ,

Polayp

Date

7/14/95

Invent d by

Recorded by

Emm

Date

7/15/95

Project No. _____

Book No. _____

TITLE

Buffers

Exhibit 134

Appl. No. 09/558,421

From Page No. _____

Buffer C, 1 L ref NB 9, pg 162

✓ 17.1 mL 1M K phos monobasic
 ✓ 7.9 mL 1M K phos dibasic
 ✓ 80 mL glycerol
 ✓ 149.12 g KCl CF = 2M
 ✓ 0.2 mL 0.5M EDTA
 ✓ 350 µl 14.3M BME - add
 ✓ qs to 1 L w/ H₂O

1M K phos monobasic
 68.045g
 qs to 500 mL
 1M K phos dibasic
 114.115g
 qs to 500 mL

1 mL of 50% Tween 20 + NP40 - post filter det 1mL + B-ME

Buffer D, 8 L ref NB 9 p 182

✓ 200 mL 1M Tris 7.5
 ✓ 1.6 mL 0.5M EDTA
 640 mL glycerol
 2.8 mL 14.3M BME - add
 ✓ 29.8 g KCl CF = 50mM
 qs to 8 L

make 0.01% NP40 + Tween 20 for dialysis
 & 0.05% for Heparin column - filter 1 L + add 1 mL 5% NP40 + Tween 20 + 350 µl.

2 L D + 700 µl
 400 µl 50% N

Buffer C, 500 mL ref NB 9 p 182

✓ 12.5 mL 1M Tris 7.5
 ✓ 0.1 mL 0.5M EDTA
 ✓ 40 mL glycerol add
 ✓ 0.175 mL β-me
 ✓ 74.5 g KCl CF = 2M
 ✓ qs to 500 mL

det 0.05% final
1mL/L

10.5 mL 50% Tween 20 + NP40 after filtering

T Pag N


Witness d & Understood by me,



Date

7/24/87

Inv nt d by



R c rd d by



Dat

7-18-87

small scale ext

and AmSO₄ (can see P163, 9)

Project No. _____

Book No. _____

101

No. 7/18/95

1ml Tag extract buffer (P.167, 3) with 0.05% each NP40 / Tween 20
50ul 10mM PMSF

7.6ul 14.3M β ME

25ul Tween 20 \rightarrow 50% mix of both \rightarrow 5ml each

25ul NP40

g cells -70°C Thorton shelf, Lab 16 - chip off 9503-15-764-D1-001K
The 5 Pol 85g

1ml Tag extract buffer mix w/ spat in 50ml
10ml 0.2g cells/ml + 10ml pipet Falcon

atc tune - on XL2020 \rightarrow (0 turn, turn small to min
to tune } then min 1-5 minimize, not over 70
stop

9x 30sec pulses in ice-H₂O bath, ~1min between
pulses - should turn browner

75°C in Falcon - H₂O bath, then cool in ice-H₂O

CF 200mM NaCl + 5% PEI to CF = 0.5% \rightarrow 16mL vol extract in grad
58.449/m cylinder

$$.2 \text{ mole} \times .016 \text{ L} = 0.0032 \text{ M}$$

30mL centrifuge tube

$$\frac{.2 \text{ g}}{58.449/\text{m}} = 0.0032 \text{ M}$$

0.187g NaCl ✓

1.45mL 1.8mL ✓
5% 16mL + 1.45mL of 5% PEI while extract is
5% stirring in c. tube, 1 drop/sec

15min, 4°C

15min, 15K 5534, 4°C (~2700xg) DNA, cell debris & heat-denatured
proteins will ppt

ant into 25mL grad cylinder:

2 = FI' fraction

13.25 mL of FI' - .2 mL = 13.05mL

d AmSO₄ to help it go into soln

\Rightarrow 2.297g ammonium

To Page No. _____

d & Understood by me,

Dorwin

Date

7/24/95

Invented by

Recorded by

Paulyn Combs

Date

7-18-95

Project No. _____

Book No. _____

TITLE _____

From Page No. _____

From Page No. _____		remove 200ul aliquots of each fraction - eppendorf on i add AmSO ₄ , spin 15 min, spin down 15K.15, remove aliquot re-measure volume	
530	AmSO ₄	vol (mL)	g AmSO ₄ slowly, stir 15'
Fr I'	0	13.05	2.297g
S 1	30%	13.8	0.414g
S 2	35	13.5	0.182g wrong 0.223g more added
S 3	40 40	13.5 13.75	0.405g 0.223g 0.426g aliquot had
S 4	45	13.5	0.419
S 5	50	13.5	0.432
S 6	55	13.4	0.442
S 7	60	13.2	
S 8	65		
S 9	70		

mix: (150 ul 0.5 M Tapes 9.3 x 2 = 300ul 0.5 M Tapes 9.3

A } 6 ul 1M MgCl₂12 ul 1M MgCl₂

{ 50 ul 3M KCl

100 ul 3M KCl

made B

3 mL

6.4 mL

20 ul mix A

412 ul mix A

6.3 ul ³²P dCTP 10mCi/mLF114/95 - 12.6 ul ³²P dCTP

60 ul 10mM dNTPs

120 ul 10mM dNTPs

405 ul 3.7 mg/mL gap activated DNA

810 ul 3.7 mg activa

2.523 mL H₂O5.046 mL H₂O

3.2 mL for 60.7 runs worth, use 48 ul/rxn

6.4 mL

1:100 in Tag dilution buffer

2ul aliquot + 198ul Tag dil bu

2ul dil + 48ul reaction mix

10', 74°C

50ul rxn

stop -1 10ul 0.5M EDTA

+ 10ul stop

spot 20ul on GEC filters

spot 20ul

T Pag N

Witnessed & Understood by m ,

Dat

Inv nted by

Dat

R c rd by

Dolan

7/24/95

Dolan

7.18.95

g N			$\left(\frac{\text{cpm}}{\text{specific activity}}\right) \left(\frac{40 \mu\text{l}}{20 \mu\text{l}}\right) (= \text{p mole})$	
			p mole	units/ μl
1	CPM1			
2	7195 c. controls			
1	no cng	71.00	7.55	
2	FI'	4679.00	498	7.47
3	SI	5411.00	576	
4	S2	5860.00	623	
5	S3	5434.00	578	
6	S4	3558.00	377	
7	S5	1394.00	148	
8	S6	299.00	31.8	→ 55% ammonium sulfate
9	mix	45046.00	} 2 μl	$\bar{x} = 45048$
10	mix	44957.00		
11	mix	45141.00		

$$\text{specific activity} = \left(45048 \text{ cpm} \right) \left(\frac{50 \mu\text{l rxn}}{2 \mu\text{l}} \right) = \frac{28.2 \text{ cpm}}{\text{p mole nt}}$$

$$\frac{28,000 \text{ } \cancel{5000} (4) \text{ p mole nt}}{10,000 \text{ p mole dCTP}}$$

$$\frac{498 \text{ p mole}}{10 \text{ n mole}} \times 3 \overset{\text{for } 10'}{=} \frac{\text{units (200)}}{2} = 7.47 \text{ units in FI'}$$

$$100,000 \text{ units} / 3.5 \text{ g cells} = 28,000 \text{ u/g}$$

Lig - 22,000 u/g

To Page N _____

Read & Understood by me,

Solomon

Date

7/24/85

Invented by

Recorded by

Carlson (comb)

Date

7.18.85

PAGE 104 OF NOTEBOOK WAS BLANK

Large scale (81.5 g cells)
The prep (can see p 176, 9)

Project No. _____
Book No. _____

105

N - except use detergent in Tag ext buffer

81.5 g

9503 -15-714 01-001R The 5 vol

326 ml

Tag ext buffer (A167, 3: with
(Room temp) 0.05% each
Tween 20 / NP40

(Cf \approx 0.2 g/ml cells
stir cells in RT buffer in beaker
strain thru cheesecloth

Minigamlin

10,000 PSI

1 pass

2 mL = FRI The, spun down cells 15 min 4°C - sup = FRI The

75°C 15 min (total time after Temp reaches 75°C)
cool fast in ice ~~stirring~~ stirring

vol = 405 ml

NaCl added to Cf 200 mM

= 4.73 g ✓

$$\frac{200 \text{ mmole}}{\text{L}} \times 0.405 \text{ L} = 81 \text{ mmole} \\ = 0.081 \text{ mole} \\ \times \frac{58.44 \text{ g}}{\text{mole}} = 4.73 \text{ g}$$

PEI for Cf = 0.5% add

45 ml 5% PEI pH 7.4 ✓ add

add dropwise, stir 15' more ✓

spin GSA 13,000 RPM 30 min ✓
ammonium = divide in 2 2 bottles

Suppl = FRI / PET

To Page No. _____

I & Understood by me,

Date

Invented by

Dat

Recorded by

7-18-95

Project No. _____

Book No. _____

TITLE _____

Exhibit 135

Appl. No. 09/558,421

From Page No. _____

vol of FRI'/PEI sup = 355 mL
 removed 1/2 mL = FRI'/PEI

$\frac{351g}{1000mL} = \frac{124.6g}{355mL} + \text{ref R. Scopes. Protein Purification p. 304}$

add 124.605 g ammonium sulfate (ground up) for 55% sat
 added slowly over ~ 15 min while stirring at 4°C

stirred 15 min, 4°C

spin down pellet 13K, 30 min GSA RC-5B, 4°C → pellet contains
 (spin in two bottles to produce 2 pellets of equal size) The DNA polymerase

saved sup at 4°C in case activity didn't come down

also 1/2 mL aliquot

respin 5K 5min to pull as much liquid off the pellets as possible
 put the 2 pellets at -70°C

METHOD 5 BANK 2

0.00 CONC % B 0
 0.00 CONC % B 0
 0.00 ALIQUOT 1
 0.00 PORT. SET 6
 0.00 PORT. SET 6
 0.00 VAL. POS 1
 0.00 VAL. POS 2
 0.00 CONC % B 0
 0.00 ALIQUOT 1

FPLC method for S200 column
 method 5 bank 2

0 conc % B 0.0
 0 conc % B 0.0
 0 mL/min 1.5
 0 port. set 6.1
 0 port. set 6.1
 0 valve pos 1.1
 0 valve pos 2.1
 400 conc % B 0.0
 400 mL/min 1.5

To Page

Witnessed & Understood by me,

Date

Invented by

Date

Record d by

E. guthrie 13th upharve
and S200

Project No. _____
Book No. _____

S200
column

107

ie No. — wash with H_2O , (HiTrap Pharaminer)
equilibrate a 10 ml blue column (used once and stored in
20% EtOH P191, 9)
at $\frac{3}{2}$ col vol/hr = 0.5 ml/min
for $2\frac{1}{2}$ hr (= 7 col vol) with buffer B

wash S200 with H_2O (180 ml vol P178, 9)
Equilibrate with buffer B
at $\frac{1}{2}$ col vol/hr
= 90 ml/hr
= 1.5 ml/min

} for S200 col will equilibrate
at 0.6 ml/min
in 500 ml in 14 hr (= 3 col vol)

S200 column

Resuspend one of the two AmSO₄ pellets of P. 106
(i.e. 0.5 of the total material from 5g cells on P105)
in buffer B (~~and~~ containing detergent) (as per P178, 9)

final vol = 0.89 ml spun out insoluble material
in microfuge 15'
remove ~20 ul to assay later AmSO₄ resuspend

Load by gravity on 180 ml S200

elute at $\frac{1}{2}$ col vol/hr (i.e. 1.5 ml/min) in buffer B

spin = 2A, 2 min/min, 1.5 ml/min, 3 ml/min 2 min/min

To Page No. _____

d & Understood by me,

Date

Invented by

Date

D. Camp

7/24/91

Recorded by

7-19-91

7-20-91

Project No. _____

Book No. _____

TITLE _____

108

From Page No. _____

unit assay on fractions 3-17 from S200 column

2ul fraction + 61.25ul Tag dil buffer \Rightarrow 1:31.62

2ul + 61.25ul Tag dil buffer \Rightarrow 1:31.62 total of 1:1000 dil

2ul of 1/1000 dil + 48ul mix A p. 102

no eng control

10min 74°C - temp varied between 72-77°C

10ul STOP 0.5M CDTA

20ul on GFC - wash 1x 10%, 3x 5% TCA, 2x GTOH

dry 30min, count

$\frac{\text{cpm}}{\text{pmol}} \left(\frac{100\text{ul}}{20\text{ul}} \right)$ specific activity = 28.2 cpm/pmol

SAM	CPM1	pmol	unit
1 no eng	41.00	4.36	546 were cloudy
Fraction 23	52.00	5.53	3
From 34	51.00	5.43	4
S200 45	45.00	4.79	5
56	153.00	16.3	4
67	721.00	76.7	7 34/5
78	634.00	67.4	8 30
89	2027.00	216	9 97 1/2
910	4597.00	489	10 220
1011	301.00	32	11 4
1112	344.00	36.6	12
1213	820.00	87.2	13
1314	208.00	22.1	14
1415	531.00	56.4	15
1516	321.00	34.1	16
1617	254.00	27.0	17

prol fr 9, 10 = 8 ml total
L = 317,250 units loaded
onto Blue Sepharose

unit def.
unit \rightarrow inc 10,000 pmole nt in 30min \Rightarrow 2027 $\frac{216 \text{ pmol}}{2 \text{ ul}} = 108 \text{ pmol} \times 1000 \times$
 $\frac{1.08 \times 10^5 \text{ pmol}}{\text{ul}} \left(\frac{30 \text{ min}}{10 \text{ min}} \right) = 3.24 \times 10^5 \text{ pmole}$
 $\frac{10,000 \text{ pmole}}{\text{unit}} = 32.4 \text{ ul} \times 300$
 $= 97,200 \text{ units}$

T Page

Witness d & Und rst d by m ,

Date

Inv nted by

Dat

[Signature]

[Signature]

R cord d by

[Signature]

7-20-97

9 'N
but separator

load fr 9, 10 off S200 (P10P) $V_f = 6 \text{ ml}$
 $1 \text{ col vol/hr} = 10 \text{ ml/hr} = 0.17 \text{ ml/min}$

Wash O/W

7 col vol buffer B total = 70 ml

(0.08 ml/min \times 70 ml in 15 hr)

collect 10 \times 7 ml fractions of wash
 7 min/fr

HOD 5 BANK 2

0.00	CONC %B	0.0
0.00	CONC %B	0.0
0.00	ML/MIN	0.08
0.00	PORT. SET	6.1
0.00	PORT. SET	6.1
0.00	VALUE POS	1.1
0.00	VALUE POS	2.1
0.00	CONC %B	0.0
0.00	ML/MIN	0.08

gradient : will scale down gradient
 for 10 ml Blue col

200 ml total gradient (20 col vol)

of 50 mM - 1 M KCl

3 col vol/hr = 0.5 ml/min

3 ml/fr

6 min/fr

1 ~~mm~~ min/min

FAIL. M=5 B=2
 11.31 RE= 11.31

C. Combs
 7/21/95

start ~ 7:45 AM

To Page No. _____

I & Understood by me,

Date

Invented by

Date

Record by

7-20-95

7-21-95

g No. _____

Blue separator

load fr 9, 10 of S200 (P10P) $V_f = 6 \text{ ml}$
 $1 \text{ col vol/hr} = 10 \text{ ml/hr} = 0.17 \text{ ml/min}$

Wash O/W

7 col vol buffer B total = 70 ml

(0.08 ml/min \times 70 ml in 15 hr)

collect 10 \times 7 ml fractions of wash
 PP min/fr

HOD 5 BANK 2

.00	CONC %B	0.0
.00	CONC %B	0.0
.00	ML/MIN	0.08
.00	PORT.SET	6.1
.00	PORT.SET	6.1
.00	VALUE.POS	1.1
.00	VALUE.POS	2.1
.00	CONC %B	0.0
.00	ML/MIN	0.08

gradient : will scale down gradient
 for TPI of P182, 9 2

for 10 ml Blue col

200 ml total gradient (20 col vol)
 of 50 mM - 1 M KCl

 $3 \text{ col vol/hr} = 0.5 \text{ ml/min}$

3 ml/fr

6 min/fr

1 mm/min

PHIL. M=5 B=2
 11.31 RE= 11.31

C Combs
 7/21/95

start ~ 7:45 AM

To Page No. _____

& Understood by me,

Date

Invented by

Date

Recorded by

Date

D. Combs

7/21/95

C. Combs

7-20-95
7-21-95

From Pag No. 7/21 Unit assay

3rd P & d CTP ref 7/14/95

dilutions in Tag dilution buffer

• Fr T' / PET
355 mL

stock

Df = 1:100

(1)

• Am SO₄ pellet
0.89 mL
resuspension

stock

1:100

Df = 1:5000

(4)

• S200 pooled peak
fractions 6 mL

stock

1:10

Df = 1:500

(7)

• Fractions from Blue
Sephacose columns

undiluted wash fractions

W5 (10)

W10 (11)

F5 (12)

peak elution fractions 10-20

stock

1:10

1:15

Df = 1:150

tube

(13) - (23)

start runs by adding 2ul of each dilution to 48ul Amix p. 102

74°C 10 min

+ 10ul 0.5 M EDTA, spot & wash & spot 20ul

Witnessed & Understood by m ,

Date

Invented by

Dat

T Pag N

Bokun

7/24/95

R c rd d by
Paulm Pomb

7/21/95

e N _____

CPMI $\frac{\text{pmole}}{\text{cpm}}$
specific activity

1100 4504.00
1200 3334.00
1400 1539.00
5000 11945.00
0.000 8088.00
10000 4497.00
500 9437.00
1000 5329.00
2000 3011.00
ash 5 108.00
ash 10 98.00
5 107.00
10 651.00
11 813.00
12 1642.00
13 3668.00
14 7866.00
15 10929.00
16 6668.00
17 6788.00
18 5668.00
19 3724.00
20 2935.00

from 13-17
2nd (of 3 ml) from each fraction $\approx V_f = 10 \text{ ml}$

analysis of blue part

Analysis against 1800 ml buffer D (P100)

continued on P. 114

To Page No. _____

I & Understood by me,

DeLano

Date

7/2/57

Invented by

[Signature]

Recorded by

Carolyn Condit


Date

7-21-57

From Page No. — Unit assay on old Tne prep (5-7-95) side-by-side the new Tne pool which was eluted from Blue agarose column on 7/21/95 p.m. NB II. Sepharos

enzyme dilutions in Tag dilution buffer:

Tne (5-7-95 Liq) - ~~later~~ previously determined to be 5 u/ml, but have lost activity over time



same dilutions also made w/ Tne from Blue Sepharose column pool

start runs w/ 2ul of each dilution into 48ul mix A ref 0.102
10' 74°C stop w/ 10ul 0.5M EDTA, spot 20ul

³²P x DCTP ref 7/14/95

SAM	CPM1	pmol	u/ml
1 1/125 2730.00	352	544	6.45
2 1/250 1638.00	202	202	7.60
3 1/500 790.00	95	95	7.50
4 1/1000 388.00			
5 1/2000 232.00			
6 1/4000 143.00			
7 1/125 9450.00			
8 1/250 5021.00			
9 1/500 2872.00			
10 1/1000 1669.00	207	31.0	
11 1/2000 899.00	107	32	
12 1/4000 503.00	56	33.8	
13 no tag 67.00			
C. Comb 7/22/95			

ave ~ 7.0 u/ml

(= 193600 units total in blue pool fr 13-17)

23.2
222 cpm/pmol
as of 7-22-95

Polansky

7/24/95

7-22-95

PAGE 113 OF NOTEBOOK WAS BLANK

from Page No. _____

check conductivity after dialysis of P111

10 μ l / ml Hbuffer D 88 μ SHep col effluent 84 μ SDialysate 94 μ S

so conductivity is good and is similar to
P183, 9 for TFI

1. Load on ~ 22 ml Hep equilibrated on
with buffer D (P100)
at 0.67 ml/min (= 40 ml/hr = 2 col vol/hr)

2. wash 1 col vol

Gradient (start ~ 9:30 AM)

50 mM - 1.05 M KCl (0.67 ml/min)
using buffer D and E (2 ml KCl)

so 0 - 50% E

600 ml total gradient vol (~ 30 col vol)

so 50 - 700 mM KCl is in 20 col vol same as
P185, 9

0.5 ml/min span = 2A

To Page No _____

Witnessed & Understood by m ,

Date

Initiated by

Date

Recorded by

2000/01/05

7/24/95

7.22.95

Page No. _____

2 L buffer F (as per 91342. PRP)

2 Tris pH 8	40 ml	✓
5 M EDTA	0.4 ml	✓
1 M DTT	2 ml	✓
10% NP40 Pierce	✓	✓
1% Tween 20 Pierce	✓	✓
lysozyme	1 L	✓
H ₂ O		

2 L

Dialyze in 1 L for 5 hr. Change to another L for 5 hr.

at assay on fractions eluted from Heparin column:

do 1 dilution (1/150) of fractions # 35-53 and a series of dilutions on the fraction with the maximum UV absorbance - #44

3 dilutions in Tag dilution buffer

stock ^{2ul} 298ul Tag dilution buffer p1000

1:150 for fractions 35-53, called H35-H53

and ~~1:15~~ 1:250, 1:500, 1:1000, 1:2000, 1:4000 of H44 - see dilutions on p 112

+ 1:8000

do make 1:125 but don't do rxn w/ it

To Page No. _____

Read & Understood by me,

Date

Invented by

Date

Polansky

7/24/91

Rec'd by

7-27-91

Project No. _____

116

Book No. _____

TITLE _____

From Page No. _____

		pmol	u/ml		
fr 44	1:250	2180.00	274	10.3	} ave 9.53 u/ml
	2:500	938.00	114	8.5	
	3:1000	560.00	55	5.8	
	4:2000	398.00	44		
	5:4000	252.00			
	6:8000	146.00			
Hep frn 1:150 dilution 35-53	no E-2	7	57.00		
	fr 35	8	1024.00		
	36	9	699.00		
	37	10	1128.00		
	38	11	1127.00		
	39	12	2104.00		
	40	13	2484.00	pmole	
	41	14	2621.00	-412	9.27
	42	15	4760.00	748	16.82
	43	16	3657.00	574	12.92
	44	17	3507.00	550	12.39
	45	18	4717.00	740	16.7
	46	19	4120.00	447	14.6
	47	20	3517.00	552	12.9
	48	21	3280.00	515	11.6
	49	22	2282.00		
	50	23	1820.00		
	51	24	1301.00		
	52	25	750.00		
	53	26	513.00		

specific activity
= 19.1 cpm/p

$$\frac{412 \text{ pmole}}{2 \text{ ul}} \times 150 \times \frac{30}{10} = 9.27 \text{ u/ml}$$

$$\bar{x} = 13.3 \text{ u/ml} \pm 2.6 \text{ u/ml}$$

pool 41-48

(Pfrn) (1.34 ml/frn)
= 10.72 ml total

pooled frn 41-48
(~10.7 ml total vol)

Dealyze into buffer E
see P 115

Recovered 2.6 ml Tne after deanalysis (in 1/2)
add 2.6 ml buffer G of 4-29-55 (sup 6, 10 ml 91342-)

T Pag No

Vitr ssed & Understo d by m ,

J. Polansky

Dat

8/24/55

Inv nt d by


R c rded by

Dat

7-2295

7-22-95 24 0.67 ml/min
 20 ml Heparin Tm
 7-22-95 0.5 mm/min
 50 mm - 1 m K20
 600 ml + Fed
 400 ml + Fed

2 ml
 1.5 ml



To Page No. _____

ed & Understood by me,

Date

Inv nted by

Date

Recorded by

7/24/95

7/24/95

7-22-95

Project No. _____

Book No. _____

TITLE _____

Exhibit 139

Appl. No. 09/558,421


From Page No. _____

This experiment is detailed on p. 119-120, VB 11

The
J-7-95 Liz

Tag

5000 units 3 2.5 2.15 1 0.5 0.25 0.125



3 2.5 2.15 1 0.5 0.25 0.125

7-22-95

Blue
Seph

Heparin

Witnessed & Und rst d by me,

[Signature]

Dat

7/24/95

Inv nted by

[Signature]

R cord d by

[Signature]

Dat

7-22-95

T Pag 1

PCR rxns with Blue Sepharose and Heparin
fractionation of Tne prep 7/22/95

g N — 7/23/95

fractionation of Blue Sepharose pool fractions and Heparin pool fraction of Tne.
Tag and Tne (5-7-95) prep will be tested in PCR alongside the new Tne.
Can we use more than 2 units of new Tne and not get a smear?

Can 5-7-95 Tne prep gives a smear w/ more than 1 unit (p. 83 NB11)

ions:

hang buffer

5 cycles w/ old program i.e. lab 15 9600 A76 94°C 1 min

Test 0.125, 0.25, 0.5, 1, 1.5, 2, 2.5, 3 units

Keep the 380bp product

50ul rxn

35x

94°C 30 sec

55°C 30 sec

72°C 2 min

4°C

cocktail w/ all components except enzyme, for 34 rxns

340ul 5x Cheng

34ul 10mM dNTPs

34ul 500g/ul M13 RF

34ul 20uM anchor primer > see p. 42 NB11


34ul 20uM 6681 primer

1156ul H₂O house distilled

1632ul → 48ul / PCR tube for 9600

enzyme dilutions in Tag SB:


Tne (5-7-95 lig) and r Tag - both 5u/ul do the following dilutions



start rxns w/ 2ul of each dilution on ice, flick, spin down

Blue sepharose pool - 23.1u/ul when normalized to 5-7-95 Tne

i.e. $\frac{5}{7} (32.3) = 23.1$ p. 112, NB11



same as Tne (5-7-95) dilutions

To Page No. _____

ed & Und rst od by me,

Date

8/11/95

Invented by

Recorded by

Carolyn Comb


Date

7/27/95

Polamp

From Page N _____

Heparin fraction #44, which had peak UV absorbance
is $6.8 \frac{\mu}{\mu l}$ when normalized to Tne (5-7-95 prep)
i.e. $\frac{5}{7} (9.53) = 6.8 \frac{\mu}{\mu l}$ p. 116 NB 11



same dilutions as for Tne (5-7-95) p.

start rxns w/ $2 \mu l$ of each dilution starting w/ the
 $1.5 \mu l$ dilutions

35 cycles, stop w/ EDTA stop soln, $7.5 \mu l$
run $25 \mu l$ on gel

Result on p. 118 NB 11 - new Tne prep (7-22-95)
is not less prone to
making a smear than the
old (5-7-95) prep. So, DNA
contamination of enz should
not be Tne's main problem.
Rather, Tne may have an
intrinsic activity that makes
it "smear" more easily than
Tag

Witnessed & Understood by me,

Date

8/1/95

Investigated by

Recorded by

Date

7/27/95

T Page 1

Fr m Page No. _____

1) Can rTaq make a smear in high Mg^{2+} ? \pm primers & temp

mix A =
 for 8 rxns
 80 μ l \rightarrow 16 μ l \rightarrow 16 μ l
 8 units
 494.4 μ l H_2O
 10x PCR buffer \checkmark
 500 μ l M13 RF target \checkmark
 20 μ M anchor \checkmark
 6681 primer, 20 μ M \checkmark
 1.31 μ l rTaq \checkmark
 1.29 μ l H_2O \checkmark
 640 μ l

mix B =
 for 8 rxns
 80 μ l \rightarrow 16 μ l \rightarrow 1.31 μ l \rightarrow 1.29 μ l
 542.4 μ l H_2O \checkmark
 10x PCR b \checkmark
 10mM dNTP \checkmark
 rTaq 5 μ l \checkmark
 640 μ l \checkmark

start rxns by adding Mg^{2+}

#	1	2	3	4	5	6	7	8	9	10	11	12	13
H_2O	17.9	17	16	12	8	4	0	17.9	17	16	12	8	4

mix A \checkmark
 = + primers
 template

80

mix B \checkmark
 no primers
 no template

80

add last just before PCR

50mM $MgCl_2$	2.1	3	4	8	12	16	20	2.1	3	4	8	12	16
100 μ l rxns													

35 cycles

stop whole rxn w/ 11 μ l stop soln w/ 10x CDTA p 79 \checkmark

Witnessed & Understood by me,

Date

8/1/95

Inv nt d by

Record d by

Paul M. Pomb

Dat

7/27/95

To Page N

N _____

res. The work better (make more product) with a hot start?

mix C for 4.5 rxns = 90ul 5x Cheng ✓

301.5 ul H₂O ✓

9ul 50pg/ul m13 RF in TC from 2/94 ✓

9ul 20uM anchor primer ✓

9ul 20uM 6681 primer ✓

9ul 10mM dNTPS ✓

427.5

7/22/95

real 19.0 uM

1:30.5x dilution

6:18.3

7/22/95

177ul Tag 5B

95ul
+ 5ul Tne 0.24ul

duplicate

hot

✓ cold start

duplicates

hot start

(15)

start

(16)

(17)

(18)

stop at 25, 30, 35

cycles

20ul

11ul stop

3ul stop

added eny
after 1' 94°C denaturation
→ H5X/2

Can human spleen genomic DNA promote smear formation - bad seed?
no template, no primers

high Mg²⁺ - short smear
2.5 rxns condition

Low Mg²⁺ - long smear condition

25ul 10x PCR buffer ✓

187.5ul H₂O ✓

5ul 10mM dNTPS ✓

7.5ul 50mM MgCl₂ ✓

225ul

25ul 10x PCR buffer ✓

189.75ul H₂O ✓

5ul 10mM dNTPS ✓

5.25ul 50mM MgCl₂ ✓

225ul

2ul D ✓

5ul H₂O ✓

5ul genomic DNA ✓

5ul Tne 0.24ul

90ul D ✓

5ul H₂O ✓

5ul Tne

90ul E ✓

2.5ul H₂O ✓

2.5ul genomic DNA ✓

5ul Tne

90ul E ✓

5ul H₂O ✓

5ul Tne

(19)

(20)

(21)

(22)

remove 10ul at 15, 20, 25, 30, 35 cycles STOP tubes + 2ul stop

Read & Understood by me,

Solay

Date

8/1/95

Invented by

Recorded by

Date

7/27/95

To Page No. _____

From Page No. _____

4) Does Tne make a smear when 1 or 2 dNTPs are missing from rxn?

4 dNTP mix, 2.5 mM each \Rightarrow 20 μ l each 10 mM dNTP stock

3 dNTP mix, A G C 2.5 mM each \Rightarrow 20 μ l 10 mM A

20 μ l 10 mM G

20 μ l 10 mM C

20 μ l H₂O

use
8 μ l of
mixes
each 100
PCR rxn
for $C_F = 2$

2 dNTP mix, G T 2.5 mM each \Rightarrow 20 μ l 10 mM G

20 μ l 10 mM T

40 μ l H₂O

mix F = 200 μ l 10x PCR buffer
for 20 rxns
1480 μ l H₂O
1680 μ l

for 9.5 rxns

G 798 μ l
28.5 μ l 50 mM MgCl₂
883.5 μ l

H 798 μ l
19.95 μ l 50 mM M
8.55 μ l H₂O
883.5 μ l

20 μ l G 20 μ l G K 20 μ l
24 μ l 4 dNTP 24 μ l 3 dNTP mix 24 μ l 2 dNTP

6.95 μ l
+ 5 μ l Tne. 24 μ l
23, 24 25, 26 27, 28

29, 30 31, 32 33, 3

100 μ l 35 cycles
11 μ l stop (p. 79)

Lab 16, 9600 method 103 1:25^{PM}

T Pag N

Witnessed & Understood by me,

Bobalung

Dat

8/1/95

Invented by

Record d by

Daniel Pomb


Dat

7/27/95

ge No. _____

2% agarose gel \Rightarrow 3.6g agarose
 300mL 1X TAE
 20uL CtBr 5989


target + primers $\xrightarrow{\text{rTaq}}$ no input DNA
 1g Cl₂ 1.05 1.5 2 4 6 8 10 1.05 1.5 2 4 6 8 10

C. Comb
7/25/95

Results


- Increasing Mg^{2+} did not cause rTaq to make a smear either in the presence or absence of input DNA

hot start $\xrightarrow{\text{Tne (7-22-95 prep)}}$ cold start \rightarrow Cheng buffer
 25 30 35 25 30 35 25 30 35 25 30 35
 3uL/boul rxn

C. Comb
7/28/95

- With 3 units Tne, the hot start rxns did not smear any less than cold start rxn. I should have used 1.36 units in order to get product instead of smear redo this expt w/ 1.36 units Tne rxn

1.5 mM MgCl₂ Tne (7-22-95) 1.05 mM MgCl₂
 + genomic - genomic + genomic - genomic
 15 20 25 30 35 15 20 25 30 35 15 20 25 30 35 15 20 25 30 35


PCR buffer
3uL/rxn \Rightarrow too muchC. Comb
7/27/95

- Addition of genomic DNA did not result in ^{more} smear over time w/ 3 units Tne

redo this expt w/ 1.36 units Tne rxn

1.5 mM MgCl₂ Tne (7-22-95) 1.05 mM MgCl₂
 + genomic - genomic + genomic - genomic
 4 4 3 3 2 2 4 4 3 3 2 2
 PCR buffer 35 cycles 3uL/rxn

3 dNTP mix = AGC present
 2 dNTP mix = GT present

C. Comb
7/27/95

- No smear made when 2 dNTPs (@ C+A) are missing, so smear is probably not made by a TdT activity.

& Understood by me,

Date

8/1/95

Invented by

Recorded by

Carolyn Comb

Date

7/27/95

To Page No. _____

Project No. _____

Book No. _____

TITLE _____

126

From Page No. — 7/20/95 Unit assay of new Tne prep (7-22-95) after dialysis into Tag storage buffer.

- old 5-7-95 Tne prep & r Tag will be done too
- 3 replicates of each dilution series, done 3x independent

cocktail: 200 μ l of Taps, $MgCl_2$, KCl mix -
 2,520 μ l H_2O w/ p1000
 60 μ l 10mM dNTPS
 405 μ l 3.7 mg/mL gapped DNA
 9.3 μ l ^{32}P dCTP ref date 7/14/95
3.2 mL

Use 48 μ l per rxn

start rxns by adding 2 μ l of enzyme dilutions on ice.

Dilutions of ^{final} Heparin pooled fraction in ~~Tag~~ ~~SB~~ dilute w/ Tag di

repeat 3x →

stock 2 μ l 248 μ l Tag dilution buffer $D_f = 1:125$

20 μ l 1:2 $D_f = 1:250$

20 μ l 1:2 $D_f = 1:500$

20 μ l 1:2 $D_f = 1:1000$

20 μ l 1:2 $D_f = 1:2000$

20 μ l 1:2 $D_f = 1:4$

Dilutions of Tne (Liz 5-7-95) and r Tag both labeled 5 μ l

repeat 3x →

stock 5 μ l 1:23 $D_f = 1:62.5$

20 μ l 1:20 $D_f = 1:125$

20 μ l 1:20 $D_f = 1:250$

20 μ l 1:20 $D_f = 1:500$

20 μ l 1:20 $D_f = 1:1000$

20 μ l 1:20 $D_f = 1:2000$

- rxns stopped w/ 10 μ l 0.5M EDTA
- 20 μ l of each rxn spotted on GFC filters →

T Page 1

With ssed & Understood by m ,

10/20/95

Dat

8/1/95

Invented by

R cord d by Paul Romb

Dat

7/27/95

SAM

CPM1

ave
pmol

ave
u/pl

C. Combs
7/27/95

Project No. _____

Book No. _____

127

1	8081.00		
2	7346.00	625	11.7
3	5871.00		
4	5139.00		
5	5341.00	441	16.5
6	4559.00		
7	3009.00		
8	2963.00	255	19.1
9	2724.00		
10	1492.00		
11	1359.00	125	18.8
12	1429.00		
13	891.00		
14	899.00	79	
15	895.00		
16	490.00		
17	402.00		
18	524.00		

The
7-22-95

19.0 u
ave

⇒ use (x ul)(19.0) = 1.36 u/100 ul
⇒ 0.0716 ul / 100 ul PCR rxn

19	6160.00		
20	6476.00	582	5.4
21	7195.00		
22	4215.00		
23	4266.00	354	6.6
24	3596.00		
7	2124.00		
8	2014.00	181	6.8
9	2055.00		
10	1160.00		
11	998.00	95	7.0
12	1024.00		
13	572.00		
14	610.00		
15	609.00		
16	361.00		
17	352.00		
18	348.00		

The
Liz

6.8 u
ave

expected
u/A

⇒ we had been using 0.2 ul / 100 ul rxn
0.2 ul is really 1.36 units, not unit
as we thought based on the
unit value 5 u/ul

19	8453.00		
20	6925.00	658	6.17
21	7075.00		
22	4769.00		
23	3803.00	387	7.25
24	4613.00		
25	2896.00		
26	2565.00	240	9.00
27	2722.00		
28	1185.00		
29	1404.00	115	8.60
30	1334.00		
31	1234.00		
32	873.00	89.7	13.40
33	953.00		
34	592.00		
35	527.00		
36	509.00		

Tag

ave
8.28

↑
expected
u/A on
real units
in rxn/50

⇒ use (x ul)(8.28 u/ul) = 1.36 u/100 ul
= 0.164 ul / 100 ul
Tag

37	56483.00		
38	57656.00		
39	56427.00		

56855 = 34.1 CPM
ave pmol

To Page No. _____

Understood by me,

Date

8/1/95

Invented by

Recorded by

Carolyn Combs

Date

7/27/95

Polamp

Project No. _____

Book No. _____

TITLE JA containing templates
(Myron Goodman assay)

from Page No. _____

"Fidel pri" (27mer)
7-11-95

Exhibit 141

Appl. No. 09/558,421

GAGACATGGCGTCCAGTCACGACCT
CTCTGTALCGCAGGGTCAGTGGACTAGTACGAGCTACT

27 bp


"Fidel Temp"
or "Fidel Temp"
(7-11-95)
(42mer)


#51351

This is not old "Fidel Temp"
of 1991, 10. This new one
less stringent primer
according to oligo people

ssDNA region is same as MB JBC (Crighthon & MG R
for + dGTP + dATP (get G-A mismatch at position 3 and)
+ dCTP for reverse

For test of dATP incorp opposite Template JA,
have all 4 dNTPs present at 200 μ M each
and look for pause one site before JA
(run on PAGE?)

G
C
A ↑
JA T A
G C (31) 
A T
T A
C G

primer (27) 

To Page No

/Witnessed & Understood by me,

Deanna Bump

Date

8/1/95

Invented by

Myron Goodman

Date

7-28-95

Recorded by

Project No. _____

Book No. _____

TITLE _____

JU vs JT in template

From Page No. _____

① ② ③ ④ ⑤ ⑥ ⑦ ⑧

2³ pri-fide Temp JT
 10 mM primer (PST)
 2³ pri-fide Temp JU
 2.5 mM primer (PST)

4 →

10 mM

~ 10 mM
total

4 →

10X Tag PCR buffer
 50 mM MgCl₂

10 10

10 10

1.5

3 3

3 3

4

10X Vent buffer

10 10

10 10

10 mM JMTPs

2 μl

→

2.5
JMS

rTag 0.0625%

2

2

Tne 0.0625%

2

2

0.125

Vent 0.0625%

2

2

0.013 pm
pol m

DeepVent 0.0625

2

2

0.013 pm
pol mH₂O

79 → 100 μl

82 →

79 →

82 →

0.013 pm
pol mso has
primer
overpo

preheat to 70°C
 Start with addition of pol.

remove 10 μl to 5 μl cycle seq stop at

0 5 10 20 40 60 90

Witnessed & Understood by me,

2001 amp

Dat

8/1/95

Inv nt d by

Rec rd d by

Dat

7-27-95

T Pag No

je N .

segment Rxn same as P 27, 4 and 90 11
using 32p pri fidel Temp with ST cond JH

To Page No. _____

ed & Understood by me,

Dat

Invented by

Date

Polay

8/1/95

Recorded by

7-28-95

Project No. _____

Book No. _____

TITLE Effect of annealing temperature on
The smear

128

From Page No. _____

- anneal at 50°C, 60°C, 70°C w/ 72°C extension - no input DNA
- take out rxn aliquots after 10, 15, 20, 25, 30, 35 cycles
- use 1.36 units Tne (5-7-95 Lig) / 100ul rxn
- test the effect of annealing temp on both the long (low Mg²⁺) and short (hi Mg²⁺) smears.

mix A = high Mg²⁺ for small smear - enough for 3.5 rxns280ul H₂O ✓

35ul 10x PCR buffer ✓

7ul 10mM dNTPs C_f = 200uM ✓10.5ul 50mM MgCl₂ C_f = 1.5mM ✓

332.5ul

mix B = low Mg²⁺ for long smear - enough for 3.5 rxn283.15ul H₂O

35ul 10x PCR buffer ✓

7ul 10mM dNTPs

7.35ul 50mM MgCl₂ C_f = 1.05mM

332.5ul

* 25 fold dilution

The stock (Lig:

50ul (real units p. 127)

3ul Tne

72ul Tag

75ul of 0.

annealing
temp

50°C

Lab 15

60°C

Lab 16

72°C

70°C

Lab 562

rxn #

1

2

3

4

5

6

596 received mix
that were made up
a different time to
the mix A+B for 1

1 ul Mix A
hi Mg

95

95

95

1 ul Mix B
low Mg

95

95

95

* ul Tne 0.24ul
5-7-95
Lig

5

100ul rxns

remove 10ul of rxn to 2ul STOP soln w/ EDTA p. 79 NB11
at cycles 10, 15, 20, 25, 30, 35

To Page 1

Witness d & Understood by me,

Dat

Inv nted by

Dat

D. Olamp

8/1/95

Record by

D. Olamp

7/28/95

g N .—

ling Lab 15 program 76

	94° 1 min
35%	94° 30 sec
	50° 30 sec
	72° 2 min
	4°C —

Lab 11 program 103

4.5 min per cycle

94° 1 min

35x { 94° 30 sec

60° 30 sec

72° 2 min

4° —

SG1 program 133 94°, min
44° 12 min 30 sec

2 temp. PCR 35x $\begin{cases} 94^{\circ} 30 \text{ scc} \\ 72^{\circ} 2 \text{ min } 30 \text{ scc} \\ 4^{\circ} - \end{cases}$

2) Agarose gel, 36 samples + primers
2) short 50° 60° 72°
3 1-10 10 1-35; 3-10 10 3-35; 5-10 10 5-35
long 50° 60° 72°
3 2-10 10 2-35; 4-10 10 4-35; 6-10 10 6-35
primers, 7/15/95 Tne, 7/22/95 Tne

Time (5-7-95 Liz) 1.36 μ /run

79	50°C					60°C					72°C						
10	15	20	25	30	35	10	15	20	25	30	35	10	15	20	25	30	35

Paterson 86-22-4
241 56-4-5

Palchun ✓ 56-28-4
TUL 56-2-5

1.5 mm
mgCl₂

1.05 mm
MgCl₂

7/23/52

Comb
28/95

To Page No._____

sed & Understood by me,

Date**Invented by****Date**

Recorded by

7/28/95

Project No. _____

Book No. _____

TITLE *unit array for 1.1x Tag*
same as P121, 9 and P52, 10

From Page No. _____

Rxn#

*ul/array**1% Tween 20/NP40*

4°C #10 (P121, 9) (no det)

1-3

2

1 *ul*

11

1.1X

4-5

3.64

Exhibit 143

Appl. No. 09/558,421

1

(2x 0.10% T/N)

7-9

2

14

(2x TFI buffer)

10-12

2

rTag (same as P121, 9) *K5 del*

13-17

2

~~1.1X~~ 1.1X *collected*

18-20

2

"new" on P34 (5.8.95)

Bontemp #11 1.1X

21-23

3.64

1-27-95

-20 5 months on P52

24-26

2

(see Rxn# 21-23 on P12)

got 54% recovery (on P53)

*% of zero time**P122, 9*

°C SAM CPM1

*ave**ul/ul*

10 1 1779.00

2029

.008

38%

10 2 1970.00

10 3 2337.00

11 4 8375.00

.033

103

11 5 8284.00

11 6 8267.00

1 7 10246.00

9774

.039

105

1 8 9851.00

14 9 9556.00

14 10 8959.00

5484

.037

107

14 11 9908.00

14 12 9584.00

14 13 10530.00

Tag 14 14 9527.00

10,119

(.04 by definition)

14 15 9706.00

14 16 9859.00

14 17 10773.00

14 18 6924.00

7063

.028

(was .025 on P53 so >100%)

14 19 7046.00

14 20 7219.00

14 21 6156.00

6257

.025

77% (P154, 9 is 0 time po (P38, 10))

14 22 6520.00

14 23 6095.00

14 24 5038.00

.019

 $\frac{.017}{.030} = 63\%$ recovered

14 25 4980.00

14 26 4755.00

14 27 82.00

14 28 119284.00

14 29 121726.00

To Page 1

Witness d & Underst d by m ,

Dat

Invented by

Date

DD Poling

8/1/95

Record d by

7-31-95

Project No. _____

Book No. _____

TITLE Hot vs Cold start PCR w/ Tne (5-

130

From Page No. — Is more specific product made and less smear in a hot vs cold start
Look at products at 10, 15, 20, 25, 30, 35 cycles

• start 100ul rxns w/ 2ul enz in order to keep [glycerol] low \rightarrow SB
 50% rxns will be 1% , do duplicate hot + cold start rxns

• materials: ~~5x stock of Tne~~

[A] mix: 100mM Tricine pH 9 \Rightarrow 100ul 1M Tricine pH 9 ✓
 5.25mM MgOAc 5.25ul 1M MgOAc ✓
 4.25mM KOAc 212.5ul 2M KOAc ✓
 682.25ul H₂O ✓
1ml

cocktail for 4.5 reactions = 90ul 5x [A] mix ✓

315ul H₂O ✓

[B]

9ul 500g/ml M13RF ✓

9ul 20uM anchor primer ✓

9ul 20uM primer 1 ✓

9ul 10mM dNTPs ✓

441ul

98ul B
 2ul Tne 0.68u/l

①

98ul B
 2ul Tne

②

98ul B
 2ul Tne

③

98ul B
 2ul Tne

④

cold start

hot start - emp a
 2:40 PM denat
 of 1st

Tne dilution: 3ul Tne (5-7-95 Li3) 0.8u/l ✓
 27ul Tne SB ✓

30ul of 0.68u/l Tne

T Pag N

Witness d & Underst d by me,

Dat

Inv nt d by

Dat

8/1/95

R cord d by

Paula Pomb

7/31/95

0001amp

g N _____

have 10ul after 10, 15, 20, 25, 30, 35 cycles + add 2ul STOP w/ 100mM EDTA
 kb 14 9600 program 10.3 = 94°C 1min

35x { 94°C 30sec
 55°C 30sec
 72°C 2min
 4°C —

% agarose gel

top

→ 23 1-10C to 1-35C thr 4-35C 2
 25 19-10C 22-35C 2

lt:

Tne (5-7-95 Liz prep) 1.30 w/rxn

Cold Start - duplicates

Hot Start - duplicates

cycles

10	15	20	25	30	35	10	15	20	25	30	35	10	15	20	25	30	35
----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

buffer = 20mM Tricine

85mM KOAc

1.05mM MgCl₂

31145

- 380bp m.3 product

Conclusion: Hot start did not result in more specific product.
 One cold start duplicate run failed, no apparent reason.

To Page No. _____

ed & Understood by me,

Date

8/1/95

Invented by

Recorded by

Paula E. Cumb

Date

7/31/95

Polansky

Project No. _____

Book No. _____

TITLE Addition of genomic DNA to a
The PCR rxn

from Page No. 123

same as p. 123 & 3 except only 1.36 μ l of Tne (5-7-95) will be used per rxn

19 90 μ l D ✓
5.5 μ l H₂O ✓
2.5 μ l genomic DNA ^{5/6/95} ASX2 ✓
2 μ l Tne 0.68 μ l p. 130

20 90 μ l D ✓
8 μ l H₂O ✓
2 μ l Tne

21 90 μ l E
5.5 μ l H₂O
2.5 μ l genomic DNA
2 μ l Tne

22 90 μ l E ✓
8 μ l H₂O ✓
2 μ l Tne

10 μ l aliquots removed after 10, 15, 20, 25, 30, 35 cycles + 2 μ l stop w/ 100 mM

result: Tne 1.36 μ l/rxn (5-7-95 Lig prep) - no pre-
1.5 mM MgCl₂ 1.05 mM MgCl₂
+ genomic - genomic + genomic - genomic
prese.

cycle # → 10 15 20 25 30 35 10 15 20 25 30 35 10 15 20 25 30 35 10 15 20 25 30 35



7/31/95 CC

Conclusion: 1) The smear is visible even after 10 cycles - we should run aliquots of earlier cycles to determine when the smear becomes visible.

2) The 1.05 mM Mg²⁺ long smears are more intense with genomic DNA present than without genomic DNA. T

may be a real result or it may be variation in rxns - repeat w/ triplicates

used & Understood by me,

Date

Invent d by

Dat

Polansky

8/1/95

Record d by

7/31/95

Paulen Pumb

Page N

Mg²⁺ titration in a short PCR rxn w/ Tne

(5/7/95 Liz)

case: To demonstrate that lowering [Mg²⁺] shifts the size of DNA products in smear from small to larger ~~in~~ in a short PCR rxn w/ Tne alone. This result has been observed w/ ~~at~~ 1.5 mM Mg²⁺ vs 1.05 mM Mg²⁺, but the range of concentrations has not been tested. Does the transition occur over a narrow or broad range of Mg²⁺? What happens w/ 85 mM KOAc + 1.5 mM Mg²⁺? We have not tested this condition before

real expt conditions:

MgOAc (mM) 0.9, 1.05, 1.1, 1.15, 1.2, 1.3, 1.4, 1.5, 1.7 (9 levels)

KOAc (mM) 50, 85 (2 levels)

template and primers for m13 380bp product present

35 cycles, 55°C annealing temp

1.36 μ /rxn rxn Tne - 5-7-95 Liz prepstart rxns w/ Mg²⁺serials: 25 mM MgOAc \rightarrow $\frac{25 \text{ mM} \cdot 100 \mu\text{L}}{1000 \text{ mM}} = 25 \mu\text{L} \text{ 1M MgOAc}$
9.75 μL H₂O

A - 50 mM KOAc for 9.5 rxns

mix B = 85 mM KOAc 9.5 rxns

19 μL 1M Tricine pH 9.0 ✓23.75 μL 2M KOAc ✓73.4.35 μL H₂O ✓19 μL 10 mM dNTPs Cf=200 μM ✓19 μL 20 μM 6681 primer ✓19 μL 20 μM anchor primer ✓19 μL 50 $\mu\text{g}/\mu\text{L}$ m13 RF ✓1.9 μL Tne (5-7-95 Liz prep) \Rightarrow 1.36 μ /rxn855 μL \leftarrow 6.8 μ / μL 19 μL 1M Tricine pH 9 ✓40.375 μL 2M KOAc ✓717.725 μL H₂O ✓19 μL 10 mM dNTP ✓19 μL 20 μM 6681 ✓19 μL 20 μM anchor ✓19 μL 50 $\mu\text{g}/\mu\text{L}$ m13 ✓1.9 μL Tne 5-7-95

855

top of gel

no tube
9 - I dropped itthere is
a tube

12 pm - 3 pm

1

2

3

4

5

6

7

8

9

10-18

18

Lab 10

9600

10.3

90 μL

1

90 μL — 1

6.4 5.8 5.6 5.4 5.2 4.8 4.4 4 3.2 same series as 1-9

MgOAc 3.4 4.2 4.4 4.6 4.8 5.2 5.6 6 6.8

add Mg, mix well & keep on ice til cycling

100 μL rxnsstop rxns w/ 100 mM EDTA, run 20 μL on 1.2% agarose gel

I & Und rstood by me,

Date

Invented by

Date

Suzanne Palacios

7/1/95

Record d by

8/1/95

Addition of genomic DNA to a Tne
short PCR rxn - 3, 6, 10, 12, 15, 20 cycle aliquot

g N 132

purpose: To determine if adding human genomic DNA to a Tne PCR rxn leads to production of a smear by an earlier cycle # than without genomic DNA. Genomic DNA ^{200ng/100ul rxn} might act like more "bad seed" material and exacerbate the smearing. This experiment was tried on p. 132. Cycle #'s 10-35 were run on the gel. There was an indication that addition of genomic DNA made the smear darker by an earlier cycle #. Now, we are repeating the p. 132 exp in triplicate and looking at even earlier cycle #'s.

summary of experimental cond: 3 rxns w/ human spleen DNA, 3 without
materials: no primers added
start rxns by adding 1.36ul/100ul rxn w/
Tne (5-7-95 Lig prep)
55°C annealing temp program 103, Lab 16
200ng of human spleen DNA/rxn

materials: mix A: for 7.5 reactions = 75ul 10x PCR buffer
569.25ul H₂O
15ul 10mM dNTPs
15.75ul 50mM MgCl₂
Tne (5-7-95 dilution
2ul Tne
15ul 5.0
20ul 0.68
Tne

top of gel			bottom of gel		
1	2	3	4	5	6
90ul					1 ✓
5.5ul		1 ✓	8ul		1 ✓
cDNA 2.5ul		1 ✓	none		1 ✓
7ul					
(0.68ul) 2ul					1
-7-95 Lig prep					

we 10ul to a tube on ice w/ 2ul STOP soln in it (100mM EDTA)
3, 6, 10, 12, 15, 20 cycles, run 10ul on 1.2% agarose gel
6-9: 3ul 6-9 went into 6-9 stop tube c = cycle #
See result on p. 134

To Page No. 134

Read & Understood by me,

Date

Inventor

Date

Erin A. Polansky

5/7/95

Recorded by

8/1/95

Cavlyn Combs

Project No. _____

Book No. _____

TITLE _____

Exhibit 146

Appl. No. 09/558,421

From Page No. _____

Tne (5-7-95 Lig prep) 1.36 μ /run, 100ul run w/ template + primers

MgOAc (mM)

0.9 | 1.0 | 1.1 | 1.2 | 1.3 | 1.4 | 1.5

50mM
KOAc85mM
KOAc

• 1.3mM MgOAc was optimal for making product - note that we have made pro w/ 1.05mM Mg^{2+} in earlier c when glycerol + ~~NaCl~~ also pro.

• 50mM KOAc is not sufficient for product formation, but 85mM KOAc is - value between 50 + 85 have not been tested

• The size of DNA products in sme does vary from small to longer as $[Mg^{2+}]$ varies from 1.5 - 0.9mM Mg

triplicate runs w/ Tne (5-7-95 Lig prep) (see P135 for reaction)
replicate 1 replicate 2 replicate 3

cycle x 3 | 6 | 10 | 12 | 15 | 20 | 3 | 6 | 10 | 12 | 15 | 20 | 3 | 6 | 10 | 12 | 15 | 20

Conclusion.

+ 200ng human
spleen genomic
DNA per 100ul rxn

There maybe a little amount of contamin DNA still (or RNA) in TNE prep - 3 cycles needed see smear.

no genomic DNA

e.emb
3/2/95

T Page 1

With ss d & Underst od by m ,

Dat

8/7/95

Inv nt d by

R c rded by

Dat

8-2-95

From Page No. _____

purpose: To determine if Mg^{2+} controls Long PCR smear size and to see if the smear is primer/template independent. Note that Long PCR rxn uses different buffering components than we've been using for si

In a single Tne enzy PCR of 330 bp product, 1.2 mM Mg^{2+} - 1.3 mM is at the center of transition from small to long smear is optimal for product formation if 85 mM KOAc present.

do: 7 levels of Mg^{2+} , \pm target and primers, 2 ratios of Tag: Tne

materials: 1 TTS Tag Long PCR system + Kalas recommendations for 7-22-95 Tne prep

dilution of Tne in Tag SB: 2 μ l Tne (7-22-95 prep, 19 μ l/ μ l)

(1) 59.64 μ l Tag SB / 1:30.82
mix

(2) 2 μ l of dilution (1)
59.64 μ l Tag SB 1:30.82

61.64 μ l of 0.024% Tne

enzyme mixes:

Final Tag (μ):	Tne (mM)	Tag (50 μ l)	Tne (0.024%)	SB
1	1	10	2.5	15.3
1	10	10	2.5	8.1
1	2	5	2.5	7.1

mix 1 with primer and template for 16 rxns: 16 μ l 10 mM dNTP

16 μ l primer mix 1

150 μ l DNA 2 mg/ml, 8 μ l H₂O, 3 μ l H₂O, 27.2 μ l H₂O

mix 1 without primer & template, 16 rxns:

16 μ l 10 mM dNTP

304 μ l H₂O

320 μ l

T Page N

With ss d & Und rst d by m,

Dancer Backup

Dat

8/7/95

Inv nt d by

R d rd d by
Dancer P m

Dat

8/12/95

Re each of these mixes again with the 1:10 TT mix
d a mix 2 (1:2 TT) 1.6 mM Mg^{2+}

To Page No. _____

ed & Understood by me,

Date

Invented by

Date

наша Родина

817195

Recorded by

8/2/95

[REDACTED]

[illegible]

With ss d & Understood by m ,

Donal a Polay

Date . . .

8/7/95

Inv nted by

R corded by

Recorded by *Paula F. Smith*

Date

8/2/95

Project No. _____

Book No. _____

TITLE

32 P end label Pri 2 --- Pri 2
for test of decreasing primer length

140

From Page No. _____

mix for 14 labeling rxns. 23 ul 32 P ATP (3.33 uM) ref. 8/4/95 (6.00 per
 23 ul 5x KINASE buffer
 4.2 ul PPK - new lot
 50.56 ul H_2O
 146 ul

rxns = 1.66 ul of each oligo - Fidel pri 2, 3, 4, 6, 8, 10, 12, 14, 16, 18
 8.34 ul mix \leftarrow 10 uM STOCKS (16.6 pmol primer total)
 12 ul in 9600 tubes
 37°C, 30' in Lab 14 9600
 55°C, 5' "
 cool to 4°C before opening tubes

add 2 ul 10 uM Fidel temp dT - 20 pmol total
 80°C, 5'
 cool to Room temp, 15 min in 9600

add 64.7 ul 10 mM Tris pH 8.0
 stored at -20°C overnight

use 12/1000 reaction for 10 uM primer

$$\frac{\text{template}}{\text{pri}} = \frac{20 \text{ pmol}}{16.6 \text{ pmol}} = 1.2$$

To Page N

With ss d & Und rsto d by m ,

Deanna Poling

Dat

8/7/95

Inv nt d by

R cord d by

Paula Lomb

Dat

8-2-95

Proj CT NO. _____

Book N. _____

TITLE Δ [tag], Δ length of primer

142

From Page No. _____

 ^{32}P fid pri. 6 - fide Temp

" 8
 " 10
 " 12
 " 14
 " 16
 " 20

(P140, 250nm primer)

Mix A

46 μ lTag μ l

cit. buffer 5B

.0016
 .007
 .04
 .2
 1
 5

50 μ l

70°C

20 min

22 min

stop with 25 μ l aq. seq stop solutionload 3 μ l on 25% PAGE

2000 V (get ~11 mA) for 3 hr (same as P155, 7)

.0032 units (tab #1) 21 pmol at in 20 min
 (based on units a max efficiency)

input primer is 1 pmol 42 mer

= ~20 pmol at incorp for 20 nt ssDNA seq

T Page 1

Witnessed & Understood by m ,

Deborah Golub

Dat

8/7/95

Inv nted by

Recorded by



Dat

8-3-95

g No. _____

23 24 25 26 27 28 29 30 31 32 33 34 35 36

10 nM primer Cf
(1 pmol primer total/50 μ l Rxn)



(for 40 Rxns)

utions:

Mix A

10 μ l

20 μ l

30 μ l

40 μ l

50 μ l

H₂O
10X PCR buffer
50 mM MgCl₂
10 mM dNTPs

1.5-40 ml	✓
200	✓
60	✓
40	✓
<u>Total 340 ml</u>	

use 40 μ l/Rxn

Cf = 1.5 mM MgCl₂

200 μ M dNTP Cf
in Rxn.

To Page No. _____

& Understood by me,

Date

8/7/95

Invent d by

Recorded by

Date

8-3-95

PAGE 144 OF NOTEBOOK WAS BLANK

Result:

1. primer length 10 and longer and extended by Tag.
2. Will test Tne vs Tag next to see if Tne does better than Tag

To Page No. _____

ed & Understood by me,

Polamp

Date

8/14/95

Invented by

Recorded by

Date

7-3-95

Project

Book No.

TITLE

Tne vs Tag for
34 pri fid. Fide 1 Template

46

From Page No.

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17
32P fid pri 6 - fidel Temp	2																
"									2								
"																	
"																	
"																	
"																	

(P140) 250 nM primer) *preheated to 70°C w/ mix A*

Tne diluted in SB 7/22/95
R.L.

0.0156 u/l		2							2								2
0.0325 u/l			2							2							
0.25 u/l				2							2						
1 u/l					2							2					

Taq 531195

0.0156 u/l					2							2					
0.0325 u/l						2							2				
0.25 u/l							2							2			
1 u/l								2								2	

Mix A (P143)

46 uL
VF = 5 uL

20 min 70°C stop with 25 uL cycle seq
stop solution

heat to 90°C, 5 min before loading

25 uL PARE run w/ P142 and 155, 7 : run 3 hr 2000 V
mix A p. 143 scaled up 1.5X = 2.31 mL H₂O

300 uL 10X PCR buffer (from Kala ✓)

90 uL 50mM MgCl₂ (made from 1 mL 21%
50% stock)

60 uL 10mM dNTPS ✓

2760 uL

50 uL 1mM
950 uL H₂O

To Page No

Witness d & Und rst od by m ,

Deena R. Poles

Dat

8/7/95


Invent d by

R corded by
Cawley Pomb


Dat

8-4-95


22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40 41 42 43 44 45 46 47 48



time dilutions: in Tag storage buffer - 12/7/94



ie 1st normalize to Tag units p. 127 19 x 5 = 11.47 ul (Tag was thought to be 5 ul)



To Pag No. _____

ed & Understood by me,

Date

Invented by

Date


Ernest Polup

8/7/95

Recorded by

Carolyn Connors

P-4-SJ



i) primer / target dependent nonspecific products, which are seen with both imunit and iomunits Tne. The products are discrete bands, mostly ≤ 1 Kb. The ~~total~~ yield of these nonspecific products, as well as the specific 7.5 Kb product increase as $[Mg^{2+}]$ increases. The highest ratio of specific to nonspecific product occurs at 1.5 mM Mg .

2) primer/target independent products. The products form an intense smear from wells down to ≤ 200 bp. As $[Mg^{2+}]$ increases, the size of products in smear decreases and the smear becomes more intense. This is the same Mg^{2+} effect that was observed in a short PCR with Tne alone.

To Page No.

Date _____


8/5/95

Results

Data on P151, 152

at the highest [pol], Tue and Toz have ~ equal
extension for 12' at primers and Toz has
greater yield for 16 nt primer at highest level

ii. In 16mer primer, T₂ again shows first appearance of ~~amplification~~ ^{sum} of but T₂ has greater yield at 12 minuts of everyone



less Tne is needed to extend a primer.
Maybe Tne is more processive
That is consistent with Tag giving go
extensive (equal or better than Tne) at high
and ~~high~~ could explain how Equil unit
of Tag, Tne for DNase I treated gaps
DNA & act differently for primer extension.

g N — sec p. 140

id-label Fid 16: 4 μ l 32 P ATP ref 8/4/95
48.08 μ l H₂O ✓4 μ l 5x Kinase buffer ✓3.32 μ l 10 μ M Fid 16 oligo ✓ (33.2 μ M, 4.46 pmol/ μ l)4 μ l 32 P ATP ref 8/4/95 (3.33 μ M, 4.46 pmol/ μ l)0.6 μ l PNK20 μ l in 9600 tube

37°C, 30 min ✓

55°C, 5 min ✓

cool to ~4°C ✓

+ 4 μ l Fidel Temp, 10 μ M 42mer (40 pmol)

50°C, 5 min

cool for 15 min to RT in PCR machine

+ 129.4 μ l 10 mM Tris pH 8.0

store at -20°C

To Page No. _____

ed & Understood by me,

Date

8/14/95

Invent d by

Recorded by

Paulyn Comb

Date

8/7/95

Polay

PAGES 150-151 OF NOTEBOOK WERE BLANK

32P 42 mer for 3' end assay
of 3' end (-) Klenow fragment

Project No. _____ Exhibit 150
Book No. _____ Appl. No. 09/558,421 61

Erdet Template 10 μ M (42 mer)	5.8 μ l	✓	5 μ mol total
3P ATP	7	✓	
5x Klenow buffer	7	✓	
PNK 1 μ l	2	✓	
H ₂ O	13.2 μ l	✓	
	35		

37°C, 30'
70°C, 5'

1 M KPO₄ pH 7.6

1 M Kmonobasic 1.3 ml
1 M Kdibasic 8.7 ml
VP = 10 ml

To Page No. _____

Read & Understood by me,

Date

8/14/95

Inventor by

Recorded by

Date

8-8-95

Project No. _____

Book No. _____

TITLE Primer extension: time course with 4 [eng], Tag vs Tne
Fid 16 primer

52

From Page No. _____

purpose: In the last primer extension experiment, ~4x less Tne than Tag was able to extend the Fid 16 primer to a full length 42mer product. p. 150. The ^{highest} level of Tag (2u), that was tested, made even more full-length product than the highest level of Tne tested (2u). Today's expt. The purpose of today's expt is to confirm the results of p. 150 and to determine whether Tne & Tag have different affinities for DNA/primer binding or if they are the same but Tne is better at extending after a stop. ^{the time course = no eng run will reveal if Tne causes a stop} sequencing rxns and the Fid 16 primer without enzyme will also be run on the 25% gel.

materials:PCR/Mg²⁺ mix to mix with enzyme prior to beginning rxns:

	conc
25 μ l 10x PCR mix ✓	1.25x
7.5 μ l 50mM MgCl ₂ ✓	1.25x
147.5 μ l H ₂ O ✓	
<u>200 μl</u>	

→ use 16 μ l ✓+ 4 μ l of eng. dilution in 5B ✓20 μ l - 1x PCR buffer = 20mM Tris 8.4, 50mM
tubes 1-8 1x MgCl₂ Cf = 1.5mM

① mix A for 10 rxns: 80 μ l 10x PCR buffer ✓
63 μ l H₂O ✓
24 μ l 50mM MgCl₂ ✓
20 μ l 10mM dNTPs ✓
40 μ l ³²P Fid 16 annealed to Fid 1 Temp p. 14
200 μ l

Keep at 70°C

enzyme dilutions in Tag storage buffer: same as on p. 147
same preps of Tag & Tne✓ 5 μ l STOP soln. from cycle sequencing kit in 9600 tubes 1-40

To Page N

Witnessed & Understood by me,

Date

Invented by

Date

D. Polansky

8/14/95

Record d by

Paula L. Paul

8/8/95

g N	1	2	3	4	5	6	7	8
ultras	1-5	6-10	11-15	16-20	21-25	26-30	31-35	36-40
2/95	1-5							
PCR	20							
20								
all								
prewarmed to								
70°C								
Tag mix	20							
1/2 Tag mix		20						
1/4 Tag mix			20					
1/8 Tag mix				20				
1/16 Tag mix					20			
1/32 Tag mix						20		
1/64 Tag mix							20	
1/128 Tag mix								20
rrns w/ 30ul of mix A prewarmed to 70°C, triterate w/ p200								
ve 10ul of rxn to a 9600 tube containing 5ul STOP at								
1min, 2min, 5min, 10min, 20min - on ice								
t to 80°C, 5min before loading 1.5ul on 25% gel								
of loading: on 25% urea gel								
ATGC noeng 31-40								
Fid16								
16ul 1.25X PCR buffer								
4ul storage buffer								
+ 80ul mix A								
removed 10ul +								
5ul stop - from cycle sequencing kit								
loaded on 2% gel in same order								
not loaded								
8ATP								
diluted								
10 ⁴ X								
to 1x10 ³ ul								
6ul								
+ 3ul								
+ 2ul H ₂ O								
+ 1ul STOP								
load 1.5ul								
To Page No.								

I & Understood by me,

Date

8/15/95

Invented by

Recorded by

Crown Comb

Date

8/15/95

10 Rxs

✓ 2.5 ml
✓ 25

Tomy PC
 combination
 B 1)
 RSA = 0.7
 5

10x PCR buff

50 mM $MgCl_2$
50 mM sterile glycerol
H₂O (sterile)

✓

✓

✓

7

50

70 mg -

350

10

32p 42mm
Cfidel Template P.61

* See p ~~44~~ 74-75
Berkton Ockerson

480

480

3² P48-mn
no

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18
---	---	---	---	---	---	---	---	---	----	----	----	----	----	----	----	----	----

A

③

45 ml —————→

distances
with T₁
dist. buff

Know exo-
lot CK041 75 u/l

at 25 u/l

5 u/l

1 u/l

Kleiner exo(t)
wt \in JPLH 130%
til 25 u/a

Klenow (+)exo
0.6 u/l
⇒ 7/2
0.006

start 1:2

Sh,

25 μ l cycle
stop sol

lost 3 pl
8% PA

To Page N

Witnessed & Understood by me,

Polamp

Dat
$$8 \overline{) 14195}$$

Invented by:

Rec rded by

Date _____

7-9-15

From:
To:
Cc:
Subject:
Date:
Priority:

Hartman, Chris
Lasken, Roger
Rashtchian, Arub
Exo minus Klenow
Monday, August 07, 1995 3:40PM
High

Project No. _____

Task N. _____

P63, 10

63

**Exo minus Klenow
Conditions**

**Task N
PERF
buffer**

oger, The unit values for the bulk exo minus klenow are as follows.

Lot No.	U/ul
CK041	75
EJP41	130

Klenow Exo(+)						Klenow Exo(+)					
CK041			EJP41			CK041			EJP41		
50	10	2	50	10	2	50	10	2	50	10	2
1.2						1.2					
0.012						0.012					

242 mer

d & Und rstood by me,

Date

Invented by

Date

To Page No. _____

Recorded by

From Page No. _____

33 correct (p138, 9)
20 μ M4.36
✓ 2.18 μ l87
~~43.6~~ p
for

M13 mp19 s DNA (+)

✓ 200 μ l21.8 pmol
total

1M Tris pH 7.5

✓ 10.6 μ l✓ 212.78 μ l 50 mM
Tris(0.00 μ mol
total)0.00 0.1025 pmol circle / μ l \Rightarrow 743 pmol nt / μ l
use 2 μ l / 50 μ l Rxn for 1.5 nmol nt / Rxn

Mix A

(33.4 Rxns)
(use 15 μ l
next time)
103 μ l ✓Taps MgCl₂ KCl
(of p' 120, 9)32P dATP 10mCi/ml 300 μ l $\frac{\text{cpm}}{\text{min}}$ 3 ✓

dATP 10mM 33.4 ✓

dCTP 10mM 33.4 ✓

dGTP 10mM 33.4 ✓

rTaq 5 μ l 66.8 μ l ✓(200 μ M)33. mp19 66.8 μ l ✓743 pmol nt / μ l 1.163 ml ✓H₂O 1.503 ml10 μ l / 50
1.5 nmol ntuse 45 μ l / 50 μ l Rxn

Witnessed & Understood by me,

S. Polansky

Date

8/14/95

Invented by

R corded by

Date

8-8-95

T Page

[illegible]


CF
 14.60
 5.00
 DNA Project No. Total dTTP
 (P. 101) Book No. dnmol TITLE Results P65
separation

0	82.60	5.4	0	18%	assuming 1/4 incorp. is 1 (is 9.1 pmol 4)
.25	137.80	9.1	12.5		
.5	172.00	11.4	25	11%	
1	303.80	20.1	50	10	
2	463.40	30.7	100	7.5	
4	725.00	48	200	6	
6	1103.80	73	300	6	1.5 pmol at m13-4 2.5% on 352 at max

200 uM dATP, dGTP, dCTP

max

0	76.00	5.0
.25	134.40	8.9
.5	201.00	13.3
1	339.20	22.4
2	482.00	31.9
4	910.20	60.3
6	1211.80	70.2
0	82.20	5.4
.25	145.80	9.1
.5	216.20	14.3
1	344.60	22.8
2	485.00	32.1
4	1062.00	70.3
6	1408.00	93.2
0	74.00	4.9
.25	158.40	10.5
.5	201.00	13.3
1	332.60	22
2	538.20	35.6
4	896.00	59.3
6	1333.20	77.3



Processivity of Taq, Tne, and Ultima

extension of 33-mer correct primer annealed to m13mp19 ssDNA
 serial enz dilutions, 2 min extension and 10 min endpt extension
 2 units - 0.0078 units in 50ul rxns, reactions started w/ 2ul enz.
 7/95

reaction cocktail for 35 rxns = 175ul 10x PCR buffer

te: 42 ³²P primer: 1 m13 circle 1347.5ul H₂O
 * mistake, see 8/11/95 52.5ul 50mM MgCl₂
 were 2ul of the labeled 35ul 10mM dNTPs
 * annealed primer was 70ul ³²P-33mer correct annealed
 diluted w/ 70ul m13 to m13mp19 - the Kinase
 rxn was done as on

" 32 P33 correct m13" of m13 ssDNA added
 0.26ug/ul m13 stock
 where more m13 added to get 1 pri/circle

zyme dilutions in Taq SB
 Taq, 5/31/95, 5ul/ul
 5ul/ul 4ul 10ul 10ul 2 fold dilutions
 16ul SB 10ul SB 10ul SB
 1:5 1:2
 1ul/ul 0.5ul/ul 0.25ul/ul 0.125 0.0625 0.0313 0.0156
 0.0078, 0.0039ul/ul

7/31/95, 11.47 ul/ul - this value is normalized to Taq P. 147
 11.47 ul/ul 5ul
 52.35ul SB then, serial dilutions made in the
 same way as for Taq
 1:11.47
 = 8.1ul/ul

Ultima 6ul/ul Lot 0643 12/31/95, Perkin Elmer
 6ul/ul 20ul SB same dilutions as for Taq + Tne, but
 Ultima units are not normalized to
 Taq units.

Read & Understood by me, Polay	Date 8/14/95	Invent d by Carolyn Conn	Date 8/9/95
-----------------------------------	-----------------	-----------------------------	----------------

Project No. _____

Book No. _____

TITLE _____

From Page No. _____

- ~~48~~ 48ul mix in a 9600 PCR tube, preheated to 70°C
- reactions were started by adding 2ul of enzy w/ P2 and triterating w/ P200
- after 2 min at 70°C, in 9600, rxns were stopped w/ 25ul Cycle sequencing Stop solution and kept at -20°C overnight prior to loading on 8% gel

Tubes	1 — 10	11 — 20	21 — 30
	Tag	Tne	T.m.a
	0.0078 → 20 units	0.0078 → 20 units	0.0078 → 20 units

The 20 units rxns were incubated for 10 min, while all the other rxns were incubated for 2 min

T Pag N.

Witness d & Underst d by m ,

Dat

8/14/95

Inv nt d by

R cord d by

Paulm Pm

Dat

8/9/95

Project No. _____

Book No. _____

TITLE Extension of 16-mer by Tag + Inc.
with Δ MgCl₂ + Δ KCl

158

From Page No. _____

2

general overview of conditions tested:

Fix MgCl ₂	Δ KCl (mM)	each condition tested w/ 0.0312u Tag + Inc 2.0u Tag + Inc in 50ul rxns at: for 20min
1.05 mM	0 25 50 85	
1.5 mM	0 25 50 85	
Fix KCl	Δ MgCl ₂ (mM)	= 44 rxns
50 mM	1, 1.2, 1.5	

for 1, 50ul rxn: 1ul 1M Tris ^{8.5} ~~8.4~~ CF=20mM * note that real PCR buffer is pH 8.4
 41.08ul H₂O
 1.42ul 3M KCl for 85mM CF
 1.5ul 50mM MgCl₂ for 1.5mM CF
 1ul 10mM dNTP CF=200uM
 * 2ul ³²P primer on Fide Temp CF=10mM
 2ul enz to start rxn
 50ul

* End-label primer as on p. 149

mix A for 50 rxns: 50ul 1M Tris, pH 8.5 ✓ 22.5 ✓
 187.5ul H₂O ✓ 843.75 ✓
 50ul 10mM dNTP ✓ 2.1015 FHC705 P63704 22.5 ✓
 100ul ³²P primer annealed to Fide Temp (✓)
 207.5ul

T Page 1

With ss d & Understood by m ,

Date

Inv nt d by

Dat

S. D. O. Lamp

8/14/95

Record d by

David M. Pombh

8/10/95

No. _____

the expt w/ varied KCl

☒ w/ 10.5 mm MgCl₂

rxns

705.5 ul A

17.85 ul 50 mm MgCl₂ ✓✓

7.65 ul H₂O ✓✓

731 ul

run + 180.6 ul B C_f = 0 mm KCl

+ 21 ul H₂O → 48 ul / rxn

4 tubes 21 ul

180.6 ul B

+ 21 ul 250.6 mm KCl C_f = 25 mm KCl / rxn

180.6 ul B

tubes

+ 21 ul 501.2 mm KCl C_f = 50 mm KCl

180.6 ul B

+ 21 ul 852 mm KCl C_f = 85 mm KCl

☒ w/ 10.5 mm MgCl₂

705.5 ul A


25.5 ul 50 mm MgCl₂ ✓

731 ul

same but use

180.6 ul ☒ for each

serial dilution of 3M KCl stock: ✓



To Page No. _____

& Understood by me,

Date

8

14/95

Invented by

Recorded by

Chun-Yu Chen

Dat

8/10/95

2 clays

From Page No. _____

for expt w/ varied $[Mg^{2+}]$:

for 14 rxns

D

= 581 μ l A11.62 μ l 3M KCl ✓9.38 μ l H_2O ✓

C = 50mM KCl/rxn

for 4.2 rxns

602 μ l180.6 μ l D+ 21 μ l 10mM Mg^{2+} /use 48 μ l/rxnCf = 1mM Mg^{2+} /rxn180.6 μ l D+ 21 μ l 12mM Mg ✓Cf = 1.2mM Mg /rxn180.6 μ l D+ 21 μ l 15mM Mg^{2+} Cf = 1.5mM Mg dil of 50mM $MgCl$ stock:50 μ l
150mM50 μ l
+ 116.5 μ l
 H_2O
15mM100 μ l
+ 25 μ l
 H_2O
12mM100 μ l
+ 20 μ l
 H_2O
10mMEnzyme dilutions in Taq storage buffer: ^(SB)Taq, 5/31/95 stock 5 μ l (not real units)5 μ l10 μ l
+ 40 μ l
SB
10 μ l10 μ l
+ 631 μ l
SB
0.0156 μ lTne, 7/22/95 stock = 11.47 μ l, normalized to Taq p. 127 & 14711.47 μ l5 μ l
+ 52.35 μ l
SB
10 μ l10 μ l
+ 631 μ l
SB
0.0156 μ l

With ssed & Und rsto d by me,

S. Polamp

Dat

8/14/95

Inv nted by

R corded by

S. Polamp

Dat

8/10/95

To Pag 1

[illegible]

Read & Understood by me,

Polarp

Date

8/14/95

Invented by

Recorded by

Date

To Page No

Project No. _____

Book No. _____

TITLE _____

162

From Pag No. _____

25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43
	✓			/			/		/		/		/		/		/	
1.5mg oxel 48				48														
25xel 48					48													
50xel 48						48	48											
55xel 48							48	48										
50xel 1mg								48			48			48			48	
1.2mg									48			48			48			48
1.5mg										48			48			48		
0.0156 Tne								2	2	2								
0.0156 Taq											2	2	2					
1 Tne	2	2	2	2										2	2	2		
1 Taq					2	2	2	2									2	2

1700V, 40mA 12 PM -
12:50 PM

T Pag N

With ss d & Underst d by m ,

Dat

Inv nt d by

Dat

D. Polansky

8/14/95

R cord d by

Paul M. Poul

8/10/95

Project No. _____

Book No. _____

TITLE _____

Redid processivity of P155
 expect expand slower [pol] trans

164

From Page No. _____

and correct error in Pri / ~~pk~~
 from 42 to (here) 20

3rd P 32 correct • m.p. 19
 P 155

2 μ l

m.p. 19 0.26 μ g / μ l

70
 70 μ l

7.6 pmol circles
 total
 now band
 Pri / circle = 1

mix A

10 x PCR buffer
 H₂O
 50 mM MgCl₂
 10 mM dNTP

70 μ l
 175 μ l
 1347.5 μ l
 52.5 μ l
 35

33. m.p. 19

0.22 pmol circles
 per 50 μ l Rxn

(preheat to 70°C
 47 μ l mix A + 2 μ l of pol
 to start \rightarrow kill with 25 μ l
 cycle seq stop etc

VP = 1370 μ l

same (eng) as p. 155 and 5 more 2 fold dilutions

Tube 15 ✓
 20 units / 50 μ l rxn

Tube 14
 2 ✓

units / rxn

Tag = 1-15

2 1 .5 .25 .125 .063 .031

.0156 .0078 .0038 0.00194 0.000

Tag = 16-29

.000484 .000242 = tube 1

To Page N

Witnessed & Understood by me,

DD Polamp

Date

8/14/95

Inv. nted by

R c rd d by

Dat

8/11/95

From: Hartman, Chris
To: Lasken, Roger
Cc: Rashtchian, Ayoub
Subject: Exo minus Klenow
Date: Monday, August 07, 1995 3:40PM
Priority: High

Project N — Exhibit 156
Book N — Appl. No. 09/558,421


P 63 10
63

**Reaction
Dickinson
Conditions**

**Tag
PCR
buffer**

nger, The unit values for the bulk exo minus klenow are as follows:

Lot No. U/ul
CK041 75
EJP41 130



& Understood by me,

Solamp

Date

8/14/95

Invented by

Recorded by

Date

8-14-95

To Page No. _____

Test supernix for [JATP],
 [JCTP]

From Page No. _____

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24

- 7 ① (-JATP)
 ② -JCTP
 ③ (-JCTP)

4 JNTPa

2.5 μ m
 5 μ m
 10 μ m
 20 μ m
 40 μ m
 60 μ m

*M_x B #11

2.5
 5
 10
 20
 40
 60
 H₂O

5
 50 μ l

- assemble on ice
- put in 5600 4°C

remove to 70°C \rightarrow 45"
 remove to 4°C

*for B
 into 220 μ m
 (at 1.1x)
 30 μ l B
 50 μ l H₂O
 4-110 μ l

With ss d & Understood by me,

JOB 21 aug

Dat

8/21/95

Invent d by

R c rded by

Dat

8-15-95

To Page 1

23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40 41 42

✓ →

43 is mix (1) ice only
ice not 70°C before EDTA

44 is mix (3) ice only

45 2# Rxn # 15 }
46 2# Rxn # 29 }
no TCA

To Page No. _____

d & Understood by me,

Date


Invented by

Date

R corded by

Polamp

8/21/95



8-11-95

Project No. _____

Book No. _____

TITLE _____

Route P 65

70

From Page No. _____

p.mol

1	103.00	3.5
2	264.00	148
3	275.00	15.5
4	638.00	41
5	1025.00	68
6	1248.00	84
7	2048.00	140
8	56.00	2.2
9	204.00	11
10	297.00	17
11	619.00	40
12	1011.00	67
13	1430.00	96
14	1879.00	128
15	94.00	2.9
16	488.00	30
17	653.00	42
18	1332.00	89
19	2200.00	150
20	3902.00	273
21	5706.00	393
22	83.00	2.1
23	457.00	28
24	643.00	41
25	1289.00	86
26	2218.00	151
27	4406.00	304
28	4082.00	282
29	53.00	
30	57231.00	412.9 cpm/pmol

412.9 cpm/pmol

1	768.00	8.1
2	854.00	32
3	926.00	53
4	928.00	53
5	979.00	68
6	1443.00	198
7	1721.00	277
8	818.00	22
9	763.00	6.8
10	971.00	65
11	931.00	54
12	1155.00	117
13	1230.00	138
14	1473.00	207
15	14.00	
16	739.00	
17	14195.00	

not very high 2H ATP BktD

$$(14195 \text{ cpm}) \quad (50 \mu\text{L Rxn} + 10 \mu\text{L EDTA}) \quad (2 \lambda \text{ spotted}) \quad (40000 \text{ pmoles Txn}) = 10.64 \text{ cH}$$

To Page

Witnessed & Und rsto d by m ,

Date

Invented by

Date

R corded by

B ok N ._____

J N .

To Page No._____

ed & Understood by me,

Date _____

Invented by

Date

Record d by

Project No. _____

Book No. _____

TITLE Determination of how to kill DNa

166

From Page No. _____

DNaase I will be used to treat the eng prep. If DNA contamin^(primer, etc) is the cause of the smear, then the DNaase I treatment may eliminate the smear. 1st we'll establish how to kill DNaase I, after it has been mixed w/ Tne, so it won't be active during a

1 DNaase I rxn will be killed w/ EDTA + heat for 5, 10 + 30 min. A second DNaase I rxn will be killed only by heat only for 5, 10, 30. After the killing treatment, the rxn will be mixed with 0x174 RF. If the DNaase I was killed, the 0x174 won't be degraded, even after a 3hr incubation.

materials: 0x174, 0.25 μ g/ μ l in 0.1mM EDTA from LTI Lot FA370

DNaase I, 1 μ l in SB = 20mM NaOAc pH 6.5
5mM CaCl₂
50% glycerol

25mM EDTA - 50ul 0.5M EDTA pH 8
950ul H₂O

200mM Tris 8.5 (note: the DNaase I buffer is 8.4)
200ul 1M Tris 8.5
800ul H₂O

20mM MgCl₂ - 200ul 50mM MgCl₂
+ 300ul H₂O

200mM KCl - 200mM 66.7ul 3M KCl
+ 933.3ul H₂O

10x DNaase I buffer = 200mM Tris-HCl pH 8.4
LTI 20mM MgCl₂
Lot EK2410 500mM KCl

- it was a bit bubbly after mixing

0.8% TAE agarose gel w/ ETBr

To Page 1

Witnessed & Underst d by m ,

[Signature]

Date

8/11/95

Invented by

Record # by

[Signature]

Date

8/15/95

Page No. _____

DNase I rxns: in 9600 PCR tubes

14 μ l H_2O
 + 2 μ l 10x DNase I buffer
 + 4 μ l DNase I; 1 μ l $C_t = 0.2$ μ l
 20 μ l
 Kill with
 heat + EDTA

(2) 14 μ l H_2O
 2 μ l 10x DNase I buffer
 4 μ l DNase I $C_t = 0.2$ μ l
 20 μ l - only heat kill this rxn

15' at RT ($= 23^\circ C$) \rightarrow during this time 1 μ g
 DNA should be digested

remove 4 μ l ($= 0.8$ units DNase I.)
 + 4 μ l 0.25 μ g/ μ l ϕ X174
 + 1.6 μ l 10x DNase I buffer
 + 10.4 μ l H_2O
 20 μ l w/ 1 μ g

series to see how long it takes active
 DNase I to degrade 1 μ g ϕ X174

- \rightarrow immediately remove 4 μ l + 1 μ l 10x loading dye w/
 100 mM EDTA = 0' Kill time
 0 time incubation w/ ϕ X174
- \rightarrow 2' later remove 4 μ l + 1 μ l LD = 0' Kill time
 2' w/ ϕ X174
- \rightarrow 15' later remove 4 μ l + 1 μ l LD = 0' Kill time
 15' w/ ϕ X174
- \rightarrow 1 hr later remove 4 μ l + 1 μ l LD = 0' Kill time
 1 hr w/ ϕ X174
- \rightarrow 3 hr remove 4 μ l + 1 μ l LD = 0' Kill time
 3 hr w/ ϕ X174.

(1) Kill the remaining 16 μ l w/ 1.36 μ l 25 mM EDTA $C_t = 2$ mM
 (1)+(2) heat to $75^\circ C$ in 9600 \leftarrow

\rightarrow 5' of heat kill, remove 4 μ l
 for rxn (1)
 + 4 μ l ϕ X174
 + 4 μ l 200 mM KCl $C_t = 50$ mM, 16 μ l
 + 1.6 μ l 200 mM Tris 8.5 $C_t = 20$ mM, 16 μ l
 + 2.02 μ l 20 mM $MgCl_2$ $C_t = 2.02$ mM
 + 4.4 μ l H_2O RT for 20 μ l
 20 μ l \rightarrow 3 hr at $25^\circ C$

To Page No. _____

sed & Understood by me,

Date

Invented by

Date

D. O. Kemp

8/21/95

Recorded by

Cawlyn Comb

8/15/95

Project No. _____

Book No. _____

TITLE _____

From Page No. _____

remove 4ul aliquots and treat in the same way after 15' and 30' killing with heat and EDTA

for rxn ②

for the DNase I rxn that was killed by only heat, take 5', 15', and 30' killing time points by removing 4ul of rxn to a tube w/

4ul ϕ X174, 0.25 μ g/ μ l
4ul 200mM XE1
12ul 200mM Tris 8.5
10.6ul 20mM MgCl₂
4.8ul H₂O
20ul

- 3hr incubation at RT
- + 3ul 10x Loading dye w/ 100 EDTA
- run 23ul on 0.8% gel

gel order

14 wells	1Kb ladder	1ug ϕ X174	5' heat + EDTA	5' heat only	15' heat + EDTA	15' heat only	30' heat + EDTA	30' heat only	OK11 ϕ X174	OK11 2ul	OK11 15' w/ ϕ X174	OK11 1hr w/ ϕ X174	OK11 3hr w/ ϕ X174	1Kb ladder
	10ug													

active DNase I
incubated w/ 1ug ϕ X174
varying amounts of time

To Page N

With ss d & Underst od by me,

Dat

8/21/95

Inv nt d by

R corded by

Robert P. P. P.

Dat

8/15/95

DPolamp

Project

Book No. _____

TITLE

init assay for -20°C sample
from J80les and also 1.1X mix sum
as p 34, 52, 80 (not among p 18)

72

From Page No. _____

note -20°C #11 arranged
on p 52

-20 #1E 27/1
-20 #4E 27/1
-20 #7E 27/1

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21

2 2 2

2 2 2

2 2 2

1.1X 5-8-95

2 2 2

1.1X Field test
(old on p 34)

2 2 2

rTag 1/25 dil
5-30-95

2 2 2 2 2

Tag Rem mix
p 120, 8

48% →

relative
to Tagp 122, 8
(time 0)

cpm

w/pl

w/λ

-20 1E { 25 5911.00
26 6883.00
27 6801.00

6531

.03

.037

81%

see p 53
where #11

-20 4E { 28 5982.00
29 5759.00
30 5205.00

5649

.026

.029

70%

of 27/1
54% and 1

-20 7E { 31 6079.00
32 5062.00
33 7422.00

6187

.028

.033

85%

renew
% #1E 4
show bottle

1.1X { 34 4974.00
5-8-95 { 35 4594.00
36 4752.00

4773

4773

.04 = .022

.022

96%

86%

.023 (p 54)

1.1X { 37 4389.00
old Field { 38 4552.00
test { 39 4971.00

4637

4637

.0213

⇒

93%

86%

rTag { 41 8601.00
42 8299.00
43 8980.00

8686 av

.04

(by

definition)

BKGD { 44 8615.00
45 78.00
46 100897.00

2X { 47 102480.00
mix { 48 102152.00

To Page 1

Witnessed & Understood by me,

D. Polansky

Date

8/21/95

Invented by

Record by

Date

8-11-95

Project No. _____

Exhibit 160

Appl. No. 09/558,421

Book No. _____ TITLE _____

From Page No. _____

specific activity of A mix = $\frac{101843 \text{ cpm}}{40,000} \left(\frac{50 \mu\text{L}}{2 \mu\text{L}} \right) = 63.7 \text{ cpm/pmol nt}$

$$\text{pmole} = \frac{\text{cpm}}{\text{S.A.}} \left(\frac{60}{20} \right) (\text{pmol})$$

① The dies at 90°C, even 5' the activity is only of the original activity. EDTA is present (ie free Mg^{2+}) Therefore, must kill DNase I at The dies a little at = w/ EDTA, maybe 10% loss of activity.

8/17/95
cc

SAM	CPM	
1	3072.00	-147
2	1533.00	-72
3	1127.00	-53
4	516.00	-24
5	198.00	-9.32
6	3423.00	-141
7	1581.00	-74.5
8	1174.00	-55
9	475.00	-22.4
10	249.00	-11.7
11	3178.00	-150
12	2007.00	-95
13	2979.00	-140
14	2332.00	-110
15	2799.00	-132
16	2601.00	-122
17	2954.00	-139
18	2798.00	-132
19	3532.00	-160
20	1251.00	-58.9
21	280.00	-13.2
22	3472.00	-144
23	2974.00	-140
24	2605.00	-123

$$140 \text{ pmole} \left(\frac{23.1}{2} \right) / 10,000 \times 3 = 0.49 \mu\text{L} \text{ (expected 0.}$$

A no killing
A 90° 5'
A 90° 10'
A 90° 30'
A 90° 1hr
B no killing
B 90° 5'
B 90° 10'
B 90° 30'
B 90° 1hr
A 75° 30'
A 75° 1hr
A 75° 2hr
A 75° 4hr
B 75° 30'
B 75° 1hr
B 75° 2hr
B 75° 4hr
Tne 90° 10'
Tne 90° 1hr
Tne 90° 2hr
Tne 75° 1hr
Tne 75° 4hr

A = 50 units Tne : 5 units DNase I
B = 50 units Tne : 1 units DNase I

Nicking assay - after treat to kill DNase I, the rxn incubated w/ 1ug ϕX174 for 3hr at 23°C. to nick any remaining DNase I. supercoiled

① note that no untreated DNA was run as a pool control for how 1ug looks. the Tne lanes serve as a control because Lig⁺ show Tne has no endonuclease.

② The 75°C 4hr treatment EDTA and the lower level DNase was best for kill the low level of DNase degrades 1ug easily

To Page 1

With ss d & Und rst od by m ,

Date

Inv nt d by

Dat

Polamp

8/21/95

Recorded by

Dunbar Pomb

8/17/95

Repeat assay of ...
with -20°C, 4°C side by side

Temp	Run#		SAM	CPM1	ave	-20/4C
1°C	1-3	(2 planned)				
20°C	4-6					
4°C	7-9					
20°C	10-12					
1°C	13-15					
20°C	16-18					
4°C	19-21					
20°C	22-24					
base	25					
base	26					
rock	27					
rock	28					
Nase	29					
Nase	30					
rock	31					
rock	32					
1 stock	33					
1 stock	34					
1/2 stock	35					
1 stock	36					
1 stock	37					
1 stock	38					
ons:						
and 4hr Tne treated with						
eI and mock rxns:						
2ul eng, 0.47u/wl P. 175						
45ul Tag dil buffer						
47ul of 0.02u/wl						
1ul Tne stock was diluted						
0.47u/wl and then to						
4ul → 2ul 3u/wl stock of Tne 5-7-95						
+ 12.4ul Tag dil buffer						
14.4ul of 5u/wl						
2ul of 5u/wl + 19.28ul Tag dil buffer / 21.28ul of .47u/wl, then diluted						

1E	1	7816.00	8440	90%
	2	9112.00		
	3	8393.00		
20	4	7582.00	7558	
	5	7182.00		
	6	7910.00		
4E	7	6265.00	6317	102%
	8	5771.00		
	9	6916.00		
7E	10	6005.00	6443	
	11	6410.00		
	12	6913.00		
11E	13	7478.00	7711	96%
	14	7917.00		
	15	7738.00		
	16	7016.00	7385	
	17	7461.00		
	18	7679.00		
	19	5526.00	5211	98%
	20	5095.00		
	21	5012.00		
	22	5396.00	5100	
	23	5050.00		
	24	4855.00		

30' 25 2887.00 } 2835 ave 70%
 26 2783.00 }
 30' 27 2672.00 } 2814
 mock 28 2956.00 }
 4u 29 3095.00
 4u 30 3280.00
 mock 31 2833.00
 36u 32 3113.00
 33 4234.00
 34 3968.00
 35 4606.00
 36 4839.00
 5u 37 3683.00
 38 4399.00
 39 472.00
 40 139229.00

2ul stock dilution
 2ul stock 5-7
 19.28ul Tag dil
 21.28ul
 2ul
 + 45ul Tag dil
 47ul of 0.02u/wl

4041x
 2835 = 70%
 4041
 after 30' 4041 treatment
 at 75°C
 2ul eng
 45ul tag

atment of Tne with DNase I: 2 rxns of 87.6 μ l H₂O
 10 μ l 10x DNase I buffer ~~60 μ l~~
 1.39 μ l Tne, 36 μ l, 5-7-95
 1 μ l DNase I, 10 μ l
 100 μ l
 (test in PCR un-killed)
 (remove 10 μ l to 400) CF = 0.5 μ l Tne
 90 μ l 1 rxn 15 min RT 0.01 μ l DNase I
 1 rxn 4 hr RT 2.5 mm MgCl₂
 15 min 1.2 % glycerol
 0.05 mm CaCl₂
 + 5.74 μ l 50 mm EDTA CF = 3 mm
 new vol = 95.74 μ l CF_{Tne} = 0.47 %
 at 75°C for 4 hr in 9600 PCR machine
 12 pm - 4 pm, put on ice

trial reactions without DNase I: 2 rxns of 88.6 μ l H₂O
 10 μ l 10x DNase I buffer
 1.39 μ l Tne, 36 μ l, 5-7-95
 100 μ l
 remove 10 μ l from each to test in PCR
 ie. no killing treatment
 90 μ l 1 rxn 15 min RT
 1 rxn 4 hr RT
 15 min
 + 5.74 μ l 50 mm EDTA CF = 0.47 μ l
 Tne
 at 75°C for 4 hr in 9600
 put on ice

mock reaction

To Page No. _____

I & Understood by me,

Date

Invent d by

Date

Polansky

8/21/95

Recorded by

Paulson Conf

8/17/95

From Page No. _____

PCR rxns: 0.5, 1, 2, 4 units Tne
 done 8/17/95
 long smear conditions and conditions to make specific
 - test DNase I treated Tne
 - test Tne that has been through a mock DNase I treat
 - test fresh Tne

mix A - to make a long smear

- for 14 rxns 1151.22ul H₂O

28ul 1M Tris 8.5 C_f = 20mM

note: in the rxn with
 4 units Tne, the
 [KCl] will be 55mM

23.38ul 3M KCl C_f = 50mM

29.4ul 50mM MgCl₂ C_f = 1.05mM

28ul 10mM dNTPs C_f = 200uM

1260ul

mix B - to make 380bp product

- for 14 rxns 1043.98ul H₂O

28ul 1M Tris 8.5

~~39.62ul~~ 39.62ul 3M KCl C_f = 85mM

dilution of fresh Tne

36.4ul 50mM MgCl₂ C_f = 1.3mM

36.4ul 5-7-95 stock

28ul 10mM dNTPs C_f = 200uM

2ul stock

151.2ul Tne 513

28ul 20uM anchor primer C_f = 400nM

153.2ul of 0.47uM

28ul 20uM 6681 primer C_f = 400nM

28ul 50pg/ul M13 RF C_f = 100pg/rxn

mix A 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24
 90 ————— 1

mix B 90 ————— 1

- H₂O 8.9 7.87 5.74 1.49 8.9 7.87 5.74 1.49 8.9 7.87 5.74 1.49 8.9 7.87 5.74 1.49 8.9 7.87 5.74 1.49

DNase treated Tne 1.1 2.13 4.26 8.51 1.1 2.13 4.26 8.51

mock treated Tne 1.1 2.13 4.26 8.51 1.1 2.13 4.26 8.51

fresh Tne 1.1 2.13 4.26 8.51 1.1 2.13 4.26 8.51

100ul rxns started on ice.

To Page N

Witness d & Understood by m ,

Date

Invented by

Dat

Polamp

8/21/95

R cord d by

Paula Smith

8/18/95

e N _____

3 in Lab 15 9600 - 1 min 94°C

35
40¹⁰ (30 sec 94°C
30 sec 55°C annealing
2 min 72°C elongation
4°C

program 76
method links 71, 75, 74

50mm KCL
1.05mm MgCl₂

85mm KCL
1.3mm MgCl₂


The treated w/ DNase I
for 30' sou the
1.0 DNase I

95 engine
units

Tne w/ DNase				mock				Fresh				Tne w/ DNase				mock				Fresh			
0.5	1	2	4	0.5	1	2	4	0.5	1	2	4	0.5	1	2	4	0.5	1	2	4	0.5	1	2	4

Mock rxn = Tne. without
DNase I taken
through all the
DNase I treatment
steps

Fresh = untreated
The used directly
from -20°C stock



8/18/95
cc

e: 85mm KCL did not prevent the small smear from forming

To Page No. _____

d & Und rstood by me,

Polamp

Date

8/21/95

Invented by

Recorded by

Emilia Combs

Date

8/18/95

From Page No. _____

Follow P. 138, 9

lot CK041

lot EJP41

no
enzyme

①

②

③

4.1

Rising (B&D system)

147.5 μ l147.5 μ l19.5 μ l

= 0.1

0.2

* see P. 75

1750

Klenow exp(-)

lot

CK041 P. 63, 75 μ l
dil in Tag SB to 4 μ l

2.5

2.0

EJP41 P. 63, 130 μ l
dil in Tag SB to 4 μ l150 μ l

2.5

150 μ l0.33 μ l Tag SB2.0 μ l2.2 μ l (BT)

41°C in 9600

remove 20 μ l to 2.2 μ l 10x BT + 100 mM EDTA
at 1, 2, 5, 30, 90, 2 hr

start 1:2

load on 2% Agarose

well #

2-7

8-13

③

①

②

To Page

Witnessed & Understood by me,

JOPolansky

Date

8/21/95

Invented by

Recorded by

Date

8-21-95

Experiments on DNase I treated Tne

Mg²⁺ titration & mixing expt.

N _____

pose: To determine if & how Tne was damaged by the DNase I treatment. p. 178 mock rxns show low polymerase activity in the PCR, so the 75°C. treatment w/ EDTA affected Tne even if DNase I was not present. A unit assay showed most of the activity was still present (~25% died). Today, we'll add mock Mg²⁺ in case there was more free EDTA than we thought. We'll also try poisoning a fresh Tne rxn w/ the treated Tne.

* titration 1.05, 1.3, 2 mM MgOAc
 smear & product conditions
 1 unit DNase treated Tne (30' treatment)
~~XXXXXXXXXX~~
 1 unit fresh Tne from 5⁴ul stock

rx mix for 12 rxn → 24ul 1M Tricine pH 9 Cf = 20mM
 30ul 2M KOAc Cf = 50mM
 24ul 10mM dNTPs Cf = 200uM
 882ul H₂O
 960ul

[B]
 80ul mix
 1.7ul 25mM MgOAc
 Cf = 1.05mM
 2.3ul H₂O
 or 90ul/rxn

[C]
 280ul mix
 18.2ul 25mM MgOAc
 16.8ul H₂O
 Cf = 1.3mM

[D]
 280ul mix
 + 28ul 25mM MgOAc
 7ul H₂O
 Cf = 2mM

mix B + 3 ul The treated w/ DNase 30' + 7ul H₂O = 100ul

* Tne at 0.33^u/ul, as determined by unit assay p. 73 NB 10 & p. 180 NB 11

mix B + 3ul Fresh Tne + 7ul H₂O

Δ diluted 5⁴ul stock to 0.33^u/ul
 5ul of 5⁴ul Tne 5-7-95
 70.8ul of Tag SB
 75.8ul of 0.33^u/ul Fresh Tne

2 mix C + 3ul treated Tne + 7ul H₂O

3 mix C + 3ul Fresh Tne + 7ul H₂O

4 mix D + 3ul treated Tne + 7ul H₂O

5 mix D + 3ul Fresh Tne + 7ul H₂O

To Page No. 180

& Und rstood by me,

olamp

Date

8/21/95

Invented by

Recorded by

Date

8/21/95

from Page No. _____

SAMP CPNI

1	2952.00
2	2739.00
3	4871.00
4	3459.00

2846

4115

68%

} The (Obase trend) 30' dil $\frac{1}{23.5}$ } Lig Tne stock 5% 5-7-15
diluted to 0.4% 1/1 (same as Obase
trend above and then diluted
 $\frac{1}{23.5}$ unit array ~~same~~ as P73, 10conclude The lost ~70% activity from
killing 4 hr killing of Obase I at 75°C (P175)product mix for Mg^{2+} titration:for 12, 100ul rxns = 24ul 1M Tricine pH 9 $C_f = 20mM$ 51ul 2M KOAc $C_f = 85mM$ 24ul 10mM dNTPs $C_f = 200uM$ 24ul 20uM 6681 primer $C_f = 400nM$ 24ul 20uM anchor primer $C_f = 400nM$ 24ul M13RF, 50pg/ul $C_f = 100pg/ul$ 789ul H_2O

960ul

Ⓔ

for 3.5 rxns

280ul mix

14.7ul 25mM $MgCl_2$
20.3ul H_2O $C_f = 1.05mM$

315ul

use 90ul/rxn

Ⓕ

280ul mix

+ 18.2ul 25mM $MgCl_2$
16.8ul H_2O

315ul

 $C_f = 1.3mM$

Ⓖ

280ul mix

+ 28ul 25mM $MgCl_2$
+ 7ul H_2O

315ul

 $C_f = 2mM$

T Pag N

itn ss d & Understo d by m ,

Date

Invent d by

Dat

J. Polak

8/28/95

R corded by

J. Polak

8/28/95

g No. _____

actions → ~~DT~~ = Tne treated w/ DNase for 30' p. 175, 1 unit, 0.33%
 FT = Fresh untreated Tne diluted to 0.33%¹, 1 unit

smear condition
 50 mm KOAc
 20 mm Tricine pH 9
 no template no primers

product condition =
 85 mm KOAc
 20 mm Tricine pH 9
 anchor primer + 6081 primer
 on mix

KOAc] mm	1.05		1.3		2 mM		1.05		1.3		2	
Time	DT	FT	DT	FT	DT	FT	DT	FT	DT	FT	DT	FT
1 unit of enzy/rxn												
1 x	1	3	4	6	7	9	10	12	13	15	16	18

20ul of each 100ul rxn was run on a 0.8% gel p. 184 →

We had also planned to do rxns with 4 units of the DNase I-treated Tne, but there was not enough of the enzy. to set up these rxns

~~Mg²⁺ was omitted from rxn by mistake~~
 using experiment to determine if the DNase I-treated Tne has
 a "poisonous" substance in it - mix untreated Tne w/ DNase-
 treated Tne smear & product conditions 1 x untreated + 0 treated

mix for smear, 6 rxns = 12ul 1M Tricine pH 9 (Cf=20mM)
 1.5ul 2M KOAc Cf=50mM
 12ul 10mM dNTPs Cf=200uM
 441ul H₂O

480ul → use 80ul/100ul rxn

80ul mix + 3ul fresh Tne (0.33%¹) + 17ul H₂O
 80ul mix + 3ul " + 1.52ul treated Tne + 15.48ul H₂O
 80ul mix + 3ul " + 3ul " + 14ul H₂O
 80ul mix + 3ul " + 6ul " + 11ul H₂O
 80ul mix + 3ul " + 12ul " + 5ul H₂O

T Page No. _____

ed & Understood by me,

Polansky

Date

8/28/95

Inv nt d by

Recorded by

Carolyn Lewis

Date

8/31/95

From Page No. 181

2nd try at the mixing expt. - I forgot to add MgOAc to the rxns on

smear buffer: Cf = 20mM Tricine pH 9.0
50mM KOAc
1.05mM MgOAc
200μM dNTPs

no primers or template
mix 1 unit of fresh Tne w/
0, 0.5, 1, 2, 4 units of
mock-treated Tne. p. 17.

for 6 rxns = 12ul 1M Tricine pH 9 ✓
15ul 2M KOAc ✓
25.2ul 25mM MgOAc ✓
12ul 10mM dNTPs ✓
415.8ul H₂O ✓
480ul → use 80ul mix per 1, 100ul rxn

product buffer: Cf = 20mM Tricine pH 9
85mM KOAc

1.3mM MgOAc
200μM dNTPs
400nM 6681 primer
400nM anchor primer
100pg/rxn M13mp19 RF template

mix 1 unit fresh Tne w/
0, 0.5, 1, 2, 4 unit
of mock treated Tne
p. 1


A unit assay was done on the ^{mock} treated Tne on Fri 8/18/95 ^{p. 75 N}
We'll assume a concentration of 0.33^u/ul, which was the ^{an}
concentration on 8/18. The base treated Tne had not lost any more acti
at 4°C over 3 days, as shown by unit assay on 8/21/95 p. 180. It's likely

for 6 rxns = 12ul 1M Tricine pH 9 ✓ the mock didn't loose an

25.5ul 2M KOAc ✓
31.2ul 25mM MgOAc ✓
12ul 10mM dNTPs ✓
12ul 20mM 6681 ✓
12ul 20mM anchor ✓
12ul 50pg/w M13 RF in TE ✓
363.3ul H₂O ✓
480ul, use 80ul mix / 100ul rxn

To Page N

Witnessed & Understood by m ,




Dat

8/28/95

Inv nted by

R cord gby



Dat

8/22/95

g No. _____	19	20	21	22	23	24	25	26	27	28	29
mix /	80ul				1						
ct mix						80ul				1	
Tne / 4ul	3ul (unit)									1	
K-treated 0.33%ul	0	1.52	3	4	12	0	1.52	3	4	12	
✓	17	15.48	14	11	5	17	15.48	14	11	5	
	100ul										

94°C 1 min

94°C 30 sec

55°C 30 sec

72°C 2 min

4°C - hold

Method 76, Lab 15 9600 Method 103 Lab 16 9600

35
cycles

Preparation of fresh Tne: 5ul of 5%ul Tne 5-7-95 Lim stock

70.8ul Tne SIB

75.8ul of 0.33%ul fresh Tne

3% TAE gel w/ ctBr, 20ul of each rxn was run on gel

To Page No. _____

Read & Understood by me,

J. Lang

Date

8/28/95

Invent d by

Recorded by

C. Lang

Date

8/22/95

Project No. _____

B4


Book No. _____

TITLE _____

From Page No. _____

Results:

mg ²⁺ titration		Mixing expt.	RNase treatment of Tne
p. 179		p. 182	p. 185
DT = Tne treated 30' with DNase I		1 unit Fresh Tne	+ = Tne treated w/ RNase
p. 175		+ increasing units of mock treated Tne (no DNase)	- = Tne not treated w/ RNase
FT = Fresh, untreated Tne		smear cond. product cond.	smear product
50 mM KOAc		units of mock treated	
85 mM KOAc			
no template/primer			
1.05 1.3 2 1.05 1.3 2			
1 unit eng →		0.5 1 2 4 0.5 1 2 4	+ - + -

CS2
8/23/95Conclusions from mg²⁺ titration expt:

A unit assay was done on the DNase I treated Tne (see p. 180) and fresh Tne. Equal units of the DNase I treated and fresh Tne were used in the mg²⁺ titration PCR reactions. The treated Tne did not make any specific product, while the Fresh Tne did. We conclude that the 75°C incubation with EDTA (to kill the DNase I) damage Tne. Therefore we do not yet know if treating Tne w/ DNase can eliminate the "bad seed" DNA and prevent a smear.

If there was some residual DNase I activity, the product & smear could be degraded. 2 pieces of data argue against the active DNase explanation. 1) our nicking assay p. 174 shows that 1 µg of Ø174 was not degraded in 3 hr and was only nicked a little (~10%) at the DNase I killing treatment.

2) in a previous expt (p. 177) the mock-treated Tne showed the same low activity in a PCR as the DNase I treated Tne. So, the DNase was not responsible for the low amount of smear products made.

Tne + EDTA did not at 94°C. maybe adding mg²⁺ back can overcome the damage to Tne. p. 174

Further expts to try to measure thermostability of the 75°C w/ EDTA treated Tne, by cycling the Tne to 94°C before doing the unit assay → purify the DNase I-treated Tne away from DNase or "poison"-gel

Witnessed & Understood by me,

Date

8/28/95

Invented by

Recorded by

David P. Pambly

Date

8/23/95

To Page No. _____

Expt.: RNase A + 1 + treatment of Tne. Is RNA the bad seed that primes the smear.

Project No. _____

185

PCR 1 is 1+2

PCR 3 is 5+6

See result on p. 184

PCR 2 is 3+4

PCR 4 is 7+7

at 3 → RNase

#1, ~~2~~

oul worth of buffer w/ 100ul worth of dNTP, target, primers

near V V 1ul Tricine

V V 1.25ul 2m KOAc

V V 2.1 ul 25mM MgOAc

V V 2 ul 10mM dNTP

V V 2ul 6681, 20mM

V V 2ul anchor

V V 2ul m13, 50pgul

1 RNase TI } both diluted 10 fold
1 RNase A } in 10mM Tricine

35.65 ul H₂O

50ul

TI 1460u/ul

A 10mg/mL = 10ug/ul

1 ul → 40ug RNA/mL

dil 10x
ul
ul

near V V 1ul Tricine #2

V V 1.25ul 2m KOAc

V V 2.1ul 25mM MgOAc

tail → 2ul Tne Liz 5u/A

1 ul RNase I

1 ul RNase A } dil 10x in 10mM Tricine

ul H₂O

46.65ul

50

want 1 unit
→ dil to 0.5u/ul in SB⁺

15' 37°C → mix → PCR

same w/o RNase

#3, 4

ESC
8/22/95

d & Und rstood by me,

Polansky

Date

8/28/95

Invented by

Recorded by

Paula Combs

Date

8/22/95

Project No. _____

Book No. _____

TITLE _____

From Page No. _____

expt 3 RNase:

for product

1 ul Tricine 1M

✓ ✓ 2.13 ul 2M KOAc

✓ ✓ 2.6 ul 25mM MgOAc

✓ ✓ 2 ul dNTP 10mM

✓ ✓ 2 ul GGS 20mM

✓ ✓ 2 ul anchor 20mM

✓ ✓ 2 ul m13 50pg/ul

1 RNase

1 RNase

✓ ✓ 34.27 H₂O

50 ul

✓ ✓ 1 ul Tricine

✓ ✓ 2.13 2M KOAc

✓ ✓ 2.6 25mM MgOAc

2 ul Tric 0.5M

✓ ✓ mix 40.27 H₂O

1 RNase

~~it~~

✓ of same w/o RNase

#5

8/22/95
CS

#6

#7

#8

To Page 1

Witnessed & Understood by me,

S. Polak

Dat

8/28/95

Invented by

Record d by

Paula Pauli

Dat

8/22/95

eN 184

conclusions from the mixing expt. on p. 184:

The purpose of the mixing expt was to see if 0.5, 1, 2, 4 units of mock treated Tne could poison a PCR with 1 unit of fresh Tne. The mock treated Tne received the 4hr 75°C EDTA treatment but did not contain any DNase I.

The mock-treated Tne did not poison the ability of fresh Tne to make a ^{10¹⁰A} smear under the standard smear buffer conditions of 50mM KOAc, 20mM Tricine, 4.05mM MgOAc. The differences in smear intensity are probably just representative of ^{typical} variation in smear intensity. Therefore the mock-treated Tne does not contain a poison that is effective at the levels tested.

The mock-treated Tne also did not poison fresh Tne's ability to make a specific product (0.5u mock w/ 1u fresh still made product). The smears seen with 1, 2, 4u mock probably are the result from having too many total units. Using more than 1 unit Tne/100ul run always results in a smear. The unit of mock was not enough to make a smear. 1 unit mock may not be exactly the same as 1 unit of treated (which doesn't show any PCR activity in the Mg²⁺ titration expt) and the window for ~~to~~ activity is probably very narrow. - More controls & expt should be done to confirm the absence of a poison.

conclusions from RNase expt - Tne treated with RNase A + RNase T₁ and then used directly in a PCR. If RNA is the "bad seed" RNase might cure the formation of a smear and increase product yield. The RNase treatment had no effect on either smear formation or product yield. We conclude that RNA is not priming the smear reaction.

To Page No. _____

I & Understood by me,

Polansky

Date

8/28/95

Invented by

Recorded by

C. E. E. E.

Date

8/22/95

Project No. _____

184

Book No. _____

TITLE _____

From Page No. _____


Results:

Mg^{2+} titration
p. 174
DT = Tne treated 30' with DNase I p. 175
FT = fresh, untreated Tne
smear cond. 50 mM KOAc
no template/primer
product cond. 85 mM KOAc

Mixing expt.
p. 182
1 unit Fresh Tne
+ increasing units of mock treated Tne (no DNase)
smear cond. units of mock treated

RNase treatment of Tne
p. 185
+ = Tne treated w/ RNase
- = Tne not treated w/ RNase

1 unit eng →
1.05 1.3 2 1.05 1.3 2 0.5 1 2 4 0.5 1 2 4 + - + -
smear product



CS
8/23/95

Conclusions from Mg^{2+} titration expt:

A unit assay was done on the DNase I treated Tne (see p. 180) and Fresh Tne. Equal units of the DNase I treated and fresh Tne were used in the Mg^{2+} titration PCR reactions. The treated Tne did not make any specific product, while the Fresh Tne did. We conclude that the 75°C incubation with EDTA (to kill the DNase I) damage Tne. Therefore we do not yet know if treating Tne w/ DNase I can eliminate the "bad seed" DNA and prevent a smear.

If there was some residual DNase I activity, the product & smear could be degraded. 2 pieces of data argue against the active DNase I explanation. 1) our nicking assay p. 174 shows that 1 µg of p. 174 was not degraded in 3 hr and was only nicked a little (~10%) after the DNase I killing treatment. 2) in a previous expt (p. 177) the mock-treated Tne showed the same low activity in a PCR as the DNase I treated Tne. So, the DNase was not responsible for the low amount of smear products made.

Tne + EDTA died at 94°C. maybe adding Mg^{2+} back can overcome the damage to Tne p. 174

Further expts to try: measure thermostability of the 75°C w/ EDTA treated Tne, by cycling the Tne to 94°C before doing the unit assay → purify the DNase I-treated Tne away from DNase or "poison"-gel

Witnessed & Understood by me,

Date

Inventor by

Date

[Signature]

8/28/95

Recorded by
[Signature]

8/23/95

To Page No. _____

ig N _____ 2.2 X reactions, $V_p = 0.1$ each

Rxn mix (B8D)

Cf at 1X

2 HPO ₄ pH 7.6 P.61	33 μ l	16.5 μ l ✓	50 mM
2 NTP ₂ 10 mM	130.2 μ l	6.6 ✓	200 μ M each
2 SA nuclear free depleted/stratagene cat 30004157	4.4 μ l	2.2 ✓	0.1 mg/ml
50% glycerol	180.8 μ l	90.4 86.9 ✓	13.17% (includes contributed by Klenow ex.)
1 M Mg Cl ₂	4.60	2.31 ✓	7 mM
2 must go in after 4 is diluted & ppt of Mg ₂ PO ₄)			
3 correct. mp 19	50.8 μ l	26.4 ✓	note 5 μ l Kb 29.5 μ l rxn contributes
39.9: 4 λ / 50 μ l rxn vol			($\frac{5}{300}$) 50% glycerol
1.12 pmol primer			= 0.833%
1.24 pmol circle			to Cf
1.24 pmol primer = 2			∴ overall Cf
H ₂ O	360.14	183.6 ✓	is 140% glycerol in rxn
$V_p =$	649	324.5	

To Page No. _____

sed & Understood by me,

[Signature]

Date

8/28/95

Invented by

Record d by

Date

7-28-95

From Page No. _____

10X PCR buffer
10 mM MgCl₂✓ 40
✓ 56

for 4 Rxns (Cf 100)

(Cf = 1100 CPM/pmol &
200 μ M each)H₂O✓ 260
356use 55 μ l / 100 μ l Rxn

(1) (2) (3)

A

89 μ l

→

✓

(53.6 u)

The. Ciz 36 μ l

2.6

✓

The PL7-225J
19 μ l (see P 27, 11)

5

✓

95 u

rTag EKBTI
del to 36 μ l

2.6

✓

95.6

Tag 5B

2.4

2.4

✓

heat reaction 1' 94
then lower heat to 7
then start reaction

dCTP, dGMP, dTTP*

6

✓ 100 μ lby addition
of dCTP, dGMP
mix

X mix

3rd PCTP 16
100 μ lremove 10 μ l to 5 μ l cycle seq stop solution
at 1 2 5 15 30 60 90 min10 mM dNTPs 8
✓ 64

run of 70 PAGE

with dDA, dDT, and PFT (Reactions #22, 23, 24)
at top of P7P) reloaded here as
number 22, 23, 24 also. 22 is no Eng

To Page 1

With ss d & Understood by me,

S. Polansky

Date

8/28/95

Invented by

Recorded by

Date

8-28-95

86

see P155-157 in Althman

Project No. _____

Book No. _____

TITLE

Procedures

Incl FY

vs

vs AlthmanFrom Page No. 600x11Mix A procedures(for 37 rxns, 40 μ l / 5 μ l rxn)

from P.75

32 P 33 run correct. mp19

✓ ✓

130 μ l

7.8 pmol circles

circle / primer = 2

3.9 pmol primer

= 0.6 pmol circle /

H₂O

✓ ✓

1109.5 μ l

for 37 rxns

10 x PCR buffer

✓ ✓

147 (1x at 40 μ l)use 40 μ l / rxn50 mM MgCl₂

✓ ✓

55.5 μ l (1x at 50 μ l)

Ct = 1.5 mM Mg

10 mM dNTPs

✓ ✓

37

(1x at 50 μ l)1x at 50 μ l rxn = 240 μ l = 1.48 μ l

3.8% Tris-Liz

Tris-FY

(Althman)

well H

A

B

sol dil

.00005 μ l

#1

12

23

for 40 μ l rxn

.0001

2

13

24

.0002

3

14

25

.0004

4

15

21

.0008

5

16

27

.0016

6

17

28

.0032

7

18

29

.0064

8

19

30

.0128

9

20

31

.0256

10

21

32

.0512

11

22

33

50 μ l

n same as

put 2 μ l pol into 8 μ l of 1.25 x PCR buffer
 preheat to 70°C 1 min in 9600

start with 40 μ l of Mix A procedures (also preprogrammed to 70
 stop at 2 min with 25 μ l cycle seq stop sol

number "0" is Mix A procedures 40 μ l
 Tag SB 2
 H₂O 8

cycle seq stop 50 μ l
 25 μ l

To Page 1

Witnessed & Understood by me,

Date

Invented by

Dat

Record d by

Fidel Pri (P54)
100 μ M
(its a 27mer)

5 μ l ✓

(500 pmol total)

γ -³²P ATP 10 mCi/ml
(3.33 μ M)

10 ✓

(33 pmol total)

5X Knaal buffer
PNK 10⁴ U
H₂O

10 ✓

1 ✓

24 ✓

50 μ l

37°C, 30' → 70°C, 5'

[JT]

25 ✓

5

[JU]

25 ✓

5

Fidel Temp (JT)

Fidel Temp dU

P54 100 μ M

10 mM Tris pH 8

247.78

50 μ l ✓

277.8

50 μ l

500

(1000 pmol
Temp pri = 2)

Cf = 900 nM primer

↓

90°C 2 min

↓

cool slow

use 5 μ l / 50 μ l extension reaction
for Cf = 90 nM primer

Project No. _____

Book No. _____

TITLE _____

2

rom Page No. _____

(14 Rxns)

Cheng

Top PCR buffer

14 Rxns

(2 ml, 100 p. 20 except
 it has JNTP & at 1 ml each at 5x)

5 x Cheng (+JNTP)

70 μ l

✓

its 0.05 mM Mg

10 x PCR buffer
 50 mM MgCl₂
 10 mM JNTPs

70

✓

14.7

✓

1.05 mM

14

✓

Human DNA spleen
 10 μ M 1/2 267 bp 1/2

7 14

7 14

✓

14 28

14 28

✓

28

552.3 28
 545.3 503.3

✓

400 nt
 run

H₂O

534 532

672

inc 47 μ l / Rxn

672

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24 25

I

II

no wt (Lig 1-7.5)

0.25 μ l 2

6

1

2

0.4
 2.5 μ l of
 5' end of
 3' end of

0.5 2

6

1

2

1 2

2.5

10

2

1.5 2

2

4.67

2

2 2

2

2

3

✓

The 25 μ l 5.4 μ l SB

0.25 5.5 5.5 2

2

0.5 5.5 5.5 2

2

1 1.2 10.0 2

2

1.5 1 5.27 2

2

2

2 1.5 5.5 2

2

2

3 2 4.27 2

2

2

4 2.5 8.37 2

2

2

50 μ l

To Page No

Witnessed & Understood by me,

Date

Invented by

Date

Rec rd d by

506 01 amp

9/7/95

[Signature]

7-1-95

PCR with OAPDH, gloom,

Cheng vs PCR temp

for the WT and Δ SY
see P 78, 11 for wt conditions

Proj # N _____
B ok No. _____

ag No. _____



94°C 1 min

94°C 30 s
55°C 30 s
72°C 2 min } 35 cycles

4°C hold

(get BS + EDTA added $\leq 30'$ after finish
Cheng buffer Tag PCR buffer
2 units
in Tag PCR buffer

The
units



← 267

Result: only get expected 267 bp product for the Δ SY if
PCR buffer used consistent with lower ionic strength
of PCR buffer helping the less successive deleted pol
(Cheng is 55 mM KOAC compared to 50 mM)

To Page No. _____

sed & Understood by m , S. Polak	Date 9/7/95	Invented by 	Date 9-1-95
		Recorded by 	

Project

Book No.

TITLE

PCR with Tne Δ5 FY
 267 bp - 7.5 bp products

14

from Page No.

Mix A

10x PCR buffer

100

✓ for 20 PCR

5x dNTP

200 μl

✓

50 mM NaCl

21 μl

✓

Human DNA 20 ng/μl

20 μl

✓

H₂O

810.5

25.5 μl

✓

7.4%

Tne Δ5 FY

8.5 μl

✓

980 μl

mix A

1 2 3 4 5 17 P 9

48 μl

267 bp 16/12 h/g

BDNF 10 μM primers

1.366 kb primers

2.0

2.82

4.1

5.5

6.166

7.5

94°C, 1'

94°C 30S

55°C 30S

72°C

40°C

35 cycles

2 min (#1-9) or 7 min (#10-18)

start 11:35

(get BJ + EDTA in ≤ 30' after finish in case 3' ext is a problem)

start 7:40
 for 7 min elong need ~ 5 hrs
 done ~ 2:30
 need ~ 2 1/2 hrs for 2 min e
 (done ~ 11:10)

Witnessed & Understood by me,

Bob O'Camp

Date

9/7/95

Invented by

Record by


Date

9-6-95

T Page No

No. _____

2' 7'
 elongation elongation
 267 200 1.36 2.8 2.8 5.5 6.1
 292



267	+	
200	+	(for 7 min along)
1.36	+	
2	+	(for 7 min)
2.8	-	
4.1	-	
5.5	-	
6.1	-	
7.5	-	

To Page No. _____

ed & Understood by me,

Dat

9/7/95

Invented by

Recorded by

Date

9-8-95

D. Polans

From Page No. _____

10x Tag PCR buffer = 200 mM Tris 8.4, 500 mM KCl

10x ultimate buffer is 100 mM Tris pH 8.8, 10 mM KCl

I will try 0 10 25 50 mM KCl CF
(for now I'll stay with 20 mM CF Tris 8.8

(A)

1 M Tris pH 8.8	20 μ l	✓	(cf=20mM)	(for 10x 10)
50 mM MgCl ₂	21 μ l	✓	(cf=1.05mM)	
4 dNTPs	20	✓	(200 μ M ea)	
Human DNA 80 ng/l	5 μ l	✓		
9.4 μ l Tris 0.5FY	8.5 μ l	✓	(4 μ l/50 μ l as per P)	
H ₂ O	86.5 μ l	✓	(24 μ l/100 μ l PCR)	
	940			use 94 μ l/100 μ l P

① ② ③ ④ ⑤ ⑥ ⑦ ⑧
94 μ l

(A)

KCl

200 mM	5	5	5	5	5	5	5	5	cf= 1
500 mM		5			5		5		2
1 M			5			5		5	5

primers

267 bp 16/12 h/g 1

1.366 kb

H₂O

5 1
100 μ l

cycle as per P 94 with 2 min elongation

start 7:40 AM done

To Page

Witnessed & Understood by me,

0801amp

Date

9/7/95

Invented by

Recorded by

Date

9-7-95


No. _____

Result

note 267 bp product does
just begin to appear
at 35 cycles and
highest K_m Cl

level of 10 mM
concluded:

- 1) higher K_m Cl₂
helps in contrast
to expectation that
it would inhibit
this distributive
form of Tne.
- 2) target DNA
is limiting here
since high yields
were obtained for
4x more DNA
p 93 and 94



apparently need $\geq 1 \mu\text{l}$ of H/S ONA (pang)/ 50 μl Rxn
need only 0.25 μl Tne

To Page No. _____

d & Underst od by me,

Date

Invented by

Date

Recorded by

Wanda Boleyn

9/7/95

9-7-95

Prepare Tth sol. for shipping
to Roach

From Page No. _____

They want 10,000 units and some SB as LTI expect
include 100 mM KCl

Therefore

Tth (formerly thought to be TFI)
4-30-95 (see P8 for units)
4.33 u/ μ l

2.5 ml

(10825 u
total)* 2 M KCl in 50% glycerol
and 20 mM Tris pH 8

0.132 ml

 $V_f = 2.632 \text{ ml}$ final units = $\boxed{4.1 \text{ units}/\mu\text{l}}$ * KCl
glycerol
H₂O
in Tris pH 8

0.298 g ✓ (MW = 74.55)

1 ml ✓

to 2 ml ✓

20 μ l ✓ $V_f = 2 \text{ ml}$

To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Dorcas A. Polansky

9/9/95

Recorded by

9-9-95