

Project No. _____
B k N . _____

Clone The mutants in pTTQ19

39

Tag No. _____

Plan: Delete 5'→3' exo activity of Tne pol mutants (F→Y)
and 3'→5'/F-Y ex.

Digest pTTQ19 with Sph + SmaI. Clone Tne mutants
at sphI and SmaI site. The digest with sphI / SmaI (HindIII).
will put Tne pol. in frame and thus, there was need for
my other manipulation.

pTTQ19 : 5μl (0.5μg)

TE 20μl

React 3μl

SphI/SmaI 1μl/1μl

30' / 37°C → freeze until other
two are ready

pUC19 FY : 60μl DNA
6μl React 2
2μl HindIII

pUC35FY (H2) : 60μl
6μl React 2
2μl HindIII

After 30 at 37°C add 5μl 1mM dNTP mix + 1μl (5U)
Klenow — Incubate for 5' on ice. Add EDTA 2μl
H2O → S10H not → see next page.

-10	RBS	met asn ser arg gly ser val asp leu gln pro ser leu ala leu ala	ATG AAT TGC CCG GGA TCC GTC GAC CTG CAG CCA AGC TTG GCA CTG GCC	EcoRI	SmaI	BamHI	PstI	HindIII	pTTQ8					
-10	RBS	met ser leu ala ala gly arg arg ile pro gly asn ser leu ala	ATG AAT TGC CCG GGA TCC GTC GAC CTG CAG CCA AGC TTG GCA CTG GCC	EcoRI	SmaI	BamHI	PstI	HindIII	pTTQ9					
-10	RBS	met asn ser ser ser val pro gly asp pro leu glu ser thr cys arg his ala ser leu ala leu ala	ATG AAT TGC CCG GGA TCC GTC GAC CTG CAG CCA AGC TTG GCA CTG GCC	EcoRI	SstI	KpnI	SmaI	BamHI	XbaI	Sall	PstI	SphI	HindIII	pTTQ18 ✓
-10	RBS	met ser leu his ala cys arg ser thr leu glu asp pro arg val pro ser ser asn ser leu ala	ATG AAT TGC CCG GGA TCC GTC GAC CTG CAG CCA AGC TTG GCA CTG GCC	EcoRI	SstI	KpnI	SmaI	BamHI	XbaI	Sall	PstI	SphI	HindIII	pTTQ19 ←
-10	RBS	met asn leu ile thr asn ser ser ser val pro gly asp pro leu glu ser thr cys arg his ala ser leu	ATG AAT TGC CCG GGA TCC GTC GAC CTG CAG CCA AGC TTG GCA CTG GCC	EcoRI	SstI	KpnI	SmaI	BamHI	XbaI	Sall	PstI	SphI	HindIII	pTTQ181
-10	RBS	met thr met ile thr asn ser ser ser val pro gly asp pro leu glu ser thr cys arg his ala ser leu	ATG AAT TGC CCG GGA TCC GTC GAC CTG CAG CCA AGC TTG GCA CTG GCC	EcoRI	SstI	KpnI	SmaI	BamHI	XbaI	Sall/AccI	PstI	SphI	HindIII	pUC18

otide sequences of the promoter and polylinker regions of the pTTQ vectors and pUC18. Sequence extending from the -35 region of the lac or tac promoter to the dist.
polylinker is given for pTTQ8, 9, 18, 19 and 181. The comparable region of pUC18 is also shown. Unique cloning sites in the polylinker, the -35 and -10 regions of th
d the RBS are shown.

Handy

8/1/95

Recorded by Debra Watson

11/1/95

40

Project No. _____

Book No. _____

TITLE

The clones in pTTQ19 (Δ' -5'-exo)

From Page No. _____

pUC Tne FY HindIII \rightarrow blunt ended fragment \rightarrow dissolve in 17 μ l TpUC Tne 35FY HindIII \rightarrow blunt ended fragment \rightarrow dissolve in 17 μ l TSph I digestionFY35FY

DNA

17 μ l17 μ l

React 6

2 μ l2 μ l


Sph I

1 μ l1 μ l

30 min / 37°C

 \rightarrow

Run gel.

#1 & #2 \rightarrow pTTQ19 sph / sma#3 \rightarrow pUC Tne FY#4 \rightarrow pUC Tne 35FY

purify Vector & insert as a mixture by gene clean in 1

(a) pTTQ19 + 2.0 Kb (Tne FY)

H₂O.

(b) pTTQ19 + 2.0 Kb (Tne 35FY)

ligation15 μ l DNA mix4 μ l 5X buffer1 μ l 5U ligase

Ligate for 15 min at room temp.

Transform DH10B. Plate 10% & 90% culture - 30°C / ON.

To Page No.

Witnessed & Understood by me,

Date

8/2/95

Invented by

Date

7/20/95

The clones in pTTQ19. (Δ 5'-exo)

Project No. _____
 Book N. _____

Page No. _____

Result of transformation:


PTTQ Tne 35/FY	10%	105
	90%	TNTC
PTTQ Tne FY	10%	115
	90%	TNTC.

Left the plates in the bench over the weekend.

Inoculate 6 clones from Overnight growth at 30°C. for mini prep (2 mL EG/Amp).

Standard mini prep!
 Digest clones with PTTQ19.

Listed as PTTQ Tne 35FY in App. 1 was dissolve in 150 μl TE. The Sst I will be coming from



Please see Mary Long's note book # 57 (LTI Book # 3959, p183)

Signed & Understood by m , 	Date 8/1/95	Invented by Deb A. Baltage	Date 7/21/95
	To Page No. _____		