

Project No. \_\_\_\_\_

B k N . \_\_\_\_\_

Clone The mutants in pTTQ19

39

Tag No. \_\_\_\_\_

Plan: Delete 5'→3' exo activity of Tne pol mutants (F→Y)  
and 3'→5'/F-Y ex.

Digest pTTQ19 with Sph + SmaI. Clone Tne mutants  
at sphI and SmaI site. The digest with sphI / SmaI (HindIII).  
will put Tne pol. in frame and thus, there was need for  
my other manipulation.

pTTQ19 : 5μl (0.5μg)

TE 20μl

React 3μl

SphI/Sma 1μl/1μl

30' / 37°C → freeze until other  
two are ready

pUC19 FY : 60μl DNA  
6μl React 2  
2μl HindIII

pUC35FY (H2) : 60μl  
6μl React 2  
2μl HindIII

After 30 at 37°C add 5μl 1mM dNTP mix + 1μl (5U)  
Klenow — Incubate for 5' on ice. Add EDTA 2μl  
mmM → EtOH ppt → see next page.

-10	RBS	met asn ser arg gly ser val asp leu gln pro ser leu ala leu ala	ATG AAT TCC CCG GGA TCC GTC GAC CTG CAG CCA AGC TTG GCA CTG GCC	EcoRI	SmaI	BamHI	PstI	HindIII	pTTQ8					
-10	RBS	met ser leu ala ala gly arg arg ile pro gly asn ser leu ala	ATG AGC TTG GCT GCA GGT CGA CCG ATC CQC GGG AAT TCA CTG GCC	EcoRI	PstI	BamHI	SmaI		pTTQ9					
-10	RBS	met asn ser ser ser val pro gly asp pro leu glu ser thr cys arg his ala ser leu ala leu ala	ATG AAT TCG ACC TCG GTA CQC CCG GAT CCT CTA GAG TCG ACC TCC AGG CAT GCA AGC TTG GCA CTG GCC	EcoRI	SstI	KpnI	SmaI	BamHI	XbaI	Sall	PstI	SphI	HindIII	pTTQ18
-10	RBS	met ser leu his ala cys arg ser thr leu glu asp pro arg val pro ser ser asn ser leu ala	ATG ACC TTG CAT GGC TCC AGC TCG ACT CTA GAG GAT CQC CCG GAT CCG AGC TCC AAT TCA CTG GCC	EcoRI	SphI	PstI	Sall	BamHI	XbaI	SmaI	SstI			pTTQ19
-10	RBS	met asn leu ile thr asn ser ser ser val pro gly asp pro leu glu ser thr cys arg his ala ser leu	ATG AAT TTG ATT ACG AAT TCC AGC TCG GTA CQC GGG GAT CCT CTA GAG TCG ACC TCC AGG CAT GCA AGC TTG	EcoRI	SstI	KpnI	SmaI	BamHI	XbaI	Sall	PstI	SphI	HindIII	pTTQ181
-10	RBS	met thr met ile thr asn ser ser ser val pro gly asp pro leu glu ser thr cys arg his ala ser leu	ATG ACC ATG ATT ACG AAT TCC AGC TCG GTA CQC GGG GAT CCT CTA GAG TCG ACC TCC AGG CAT GCA AGC TTG	EcoRI	SstI	KpnI	SmaI	BamHI	XbaI	Sall/AccI	PstI	SphI	HindIII	pUC18

otide sequences of the promoter and polylinker regions of the pTTQ vectors and pUC18. Sequence extending from the -35 region of the lac or tac promoter to the dist.  
polylinker is given for pTTQ8, 9, 18, 19 and 181. The comparable region of pUC18 is also shown. Unique cloning sites in the polylinker, the -35 and -10 regions of th  
d the RBS are shown.

Handwritten signature

8/1/95

Recorded by Debba Bhatnagar

11/1/95

40

Project No. \_\_\_\_\_

Book No. \_\_\_\_\_

TITLE The clones in pTTQ19 ( $\Delta'$ -5'-exo)

From Page No. \_\_\_\_\_

pUC Tne FY HindIII  $\rightarrow$  blunt ended fragment  $\rightarrow$  dissolve in 17  $\mu$ l T

pUC Tne 35 FY HindIII  $\rightarrow$  blunt ended fragment  $\rightarrow$  dissolve in 17  $\mu$ l T

Sph I digestion


FY

35 FY

DNA 17  $\mu$ l  
 React 6 2  $\mu$ l  
 Sph I 1  $\mu$ l

17  $\mu$ l  
 2  $\mu$ l  
 1  $\mu$ l

30 min / 37°C  $\rightarrow$  Run gel.



#1 & #2  $\rightarrow$  pTTQ19 sph / sma  
 #3  $\rightarrow$  pUC Tne FY  
 #4  $\rightarrow$  pUC Tne 35 FY

purify Vector & insert as a mixture by gene clean in 1 H<sub>2</sub>O.  
 (a) pTTQ19 + 2.0 Kb (Tne FY)  
 (b) pTTQ19 + 2.0 Kb (Tne 35 FY)

ligation

15  $\mu$ l DNA mix  
 4  $\mu$ l 5X buffer  
 1  $\mu$ l 5U ligase

Ligate for 15 min at room temp.

Transform DH10B. Plate 10% & 90% culture - 30°C / ON.

To Page No. \_\_\_\_\_

Witnessed & Understood by me,

*[Signature]*

Date

8/2/95

Invented by

*[Signature]*

Date

7/20/95

The clones in pTTQ19. (Δ 5'-exo)

Project No. \_\_\_\_\_  
 Book N. \_\_\_\_\_

Page No. \_\_\_\_\_

Result of transformation:


PTTQ Tne 35/FY	10%	105
	90%	TNTC
PTTQ Tne FY	10%	115
	90%	TNTC.

Left the plates in the bench over the weekend.

Inoculate 6 clones from Overnight growth at 30°C. for mini prep (2 mL EG/Amp).

Standard mini prep!  
 Digest clones with PTTQ19.

Listed as PTTQ Tne 35FY in App. 1 was dissolve in 150 μl TE. The Sst I will be coming from



Please see Mary Long's note book # 57 (LTI Book # 3959, p183)

Used & Understood by m , H. J.	Date 8/1/95	Invented by Deb A. Baltage	Date 7/21/95
-----------------------------------	----------------	-------------------------------	-----------------

To Page No. \_\_\_\_\_