

DEI + $(NH_4)_2SO_4$ ppt.

Page N _____

an to continue with purification - following the ~~old~~ ^{new} protocol as for wild type Tne - p.108.

$$11 - 6.8 \text{ mL} \quad (.05)(6.8) = 2M \times \quad x = 174 \mu\text{L of } 2M \text{ KCl}$$

$$5' \text{exo} - 4.8 \text{ mL } 3.8 \quad (.05)(4.8 + x) = 2M \times \quad x = 97.4 \mu\text{L of } 2M \text{ KCl}$$

$$(.40)(6.8 + x) = 10\% \times$$

$$291 \mu\text{L} = x \quad x = 291 \mu\text{L } 10\% \text{ DEI}$$

$$(.40)(3.9 + x) = 10\% \times$$

$$143 \mu\text{L} = x \quad x = 143 \mu\text{L of } 10\% \text{ DEI}$$

Make each a Anal 50 mM KCl slowly add ~~DEI~~ a 10% DEI sol'n to a Anal [3] of .4%. vortex - let shake 30 minutes @ 4°C. spin in 2 mL eppendorf in micro-centrifuge 20 minutes @ 4°C - Save Supernatant.

60% $(NH_4)_2SO_4$ fractionation

$$\text{TY1} \quad \frac{36 \text{ g solid}}{100 \text{ mL}} = \frac{x}{4.8 \text{ mL}} \quad 2.45 \text{ g}$$

$$3'5' \text{exo} - \frac{36 \text{ g}}{100 \text{ mL}} = \frac{x}{3.5 \text{ mL}} \quad 1.26 \text{ g}$$

vortex - let shake 30 min @ 4°C
spin in 55-34 - 20,000 x g -
Decant + Save Supernatant - pellet's

To Page No. _____

Used & Understood by m ,

Man Jorgo

Date

6/20/55

Invented by

E. Jorgo

Date

06/16/55

Recorded by

From Page No. _____

Bump Heparin with 5 M NaOH - wash w/ H₂O extensive
 Equilibrate w/ Buffer A

Buffer A - Heparin -

Buffer B - Heparin

25mM Tris pH 7.4

10% glycerol

5mM Bme

.1mM PMSF

.1mM EDTA

10mM KCl

conductivity - 1.2mS

A.S.

25mM Tris pH 7.4

10% glycerol

5mM Bme

.1mM PMSF

.1mM EDTA

1.5M KCl

TY-1 - Dissolve Pellet in 10mL of Buffer A

4.5mS - conduct

Add 30 mL additional of Buffer A

2.1mS - conduct

Load 9 35mL on 2mL TBSO Heparin @ .75mL/min
 collect flow through material - wash to base line -

Gradient Program - D - 100% B @ .5mL/min - 20mL linear g
 wash 100% B - 10mL - @ .5mL/min
 collect 500 μ L fractions -

To Page 1

Witnessed & Understood by me,

Date

Initiated by

Date

Mam Longo

6/20/25

Recorded by

6/15/95

Page N _____

6/15/95

MIX Rxn

Stock

For 20 mL

SMTAPS

1 mL

50 mM $MgCl_2$ 800 μ L

2M KCl

500 μ L

1M DTT

200 μ L

10 mM dNTPs

400 μ L

act. Salmon testes

5 mL

(2.1)

- 1.1 mL dCTP
vial

20 mLs

Aliquot 500 μ L / tube store in -20°C freezer - yellow tubes -

To Page No. _____

Read & Understood by me,

Date


Invent d by

Date

May Long

6/20/95

Recorded by



6/16/95

Project No. _____

Book No. _____

TITLE

Heparin - FY-1

From Page No. _____

06/15

SAM

CPM1

EX

FY-1

1	115552.0052
2	53328.0054
3	9146.0056
4	4556.0058
5	1260.0059
6	3744.0060
7	1028.0061
8	574.0062
9	536.0063
10	346.0064
11	730.0065
12	438.0066
13	348.0067
14	21268.0068
15	668.0069
16	372.0070
17	866.0071
18	74836.0072
19	146.0073

Pool 49-SS dialyze O/N in Queso Buffer A

my 6/20/95

60

70

80

90

100

Pharmacia LKB Biotechnology

24 µl Rxn
 1 µl pack
 Sample -
 incubate @
 in 8' - qu
 w/ 10 µl of S
 EDTA - SP
 20 µl on 6
 wash
 5' 1x 10' TCI
 3' 3x 5' T.
 2x 5' to
 dry + cou
 econoflow

Pool - 49
 dialyze O/
 in again
 Queso Buf
 See p. 144

11/11/95
 11/11/95

With ssed & Understood by m ,

Date

Inv nted by

Date

T Page No.

Mar Longo

4/20/95

R c rded by

6/20/95

S. Kym

06/16/95

Hepain 3-5 cko mutant

Project No. _____

Book N. _____

143

age No. _____

06/15

0806/N50/98109/8.93

The 3-5 cko mutant

SAM

CPM1

1	266.00	20
2	324.00	23
3	1126.00	24
4	24684.00	26
5	33768.00	28
6	111394.00	30
7	78652.00	32
8	29724.00	34
9	8666.00	36
10	54.00	38
11	2912.00	40
12	1402.00	42
13	13900.00	44
14	212.00	46

Pool-
25-35

24 μ l mix
1 μ l fraction
Sample -
incubate @
74°C 8 min -
quench w/
10 μ l g-SM
EDTA -
Spot 20 μ l
on 6F1C
uban -
1x 10% TCA
17-Pi

3x 5% TCA
2x EtOH
dry +
count -

Pool-25-35
dialyze 4hrs
in QLSO
Duffer A-7
see p. 143

my 6/20/95 27

56/15/95

technology

Code No. 18-1001-44

To Page No. _____

sed & Understood by me,

Date

Inv nt d by

Date

May 16/95

6/20/95

Recorded by

6/16/95