

Page N _____

11-8-94 received plates from Kala S. for pUC19 PCR
with Top + Deep Vent and no Mn^{++} or dNTP bias.
picked 20 white colonies and 2 Blue for 2 ml LB+100 μ g
overnights at 30°C.

12-9-94

miniprep as per (p 41, 4)
using 0.5 ml cells.


resuspended DNA in 50 μ l TE + RNase

digest each (# 21, 22 ^{are from} the 2 blue colonies)
as follows (tube 23 is 0.2 ml of 0.8 μ g/ml pUC19)

1) buffer 4 (10x)	11.6	✓	✓	✓	569
2) H ₂ O	2	✓	✓	✓	94
3) R1 10 μ l	5 μ l				
4) H ₂ O	0.5	✓	✓	✓	23.5
5) H ₂ O	0.5	✓	✓	✓	
6) H ₂ O	0.1	✓	✓	✓	4.7
7) H ₂ O	0.5	✓	✓	✓	14.1
8) H ₂ O	20 μ l				

37°C, 2 hours

Result:
all are
full length
loc 2



loc 2 gives 2 bands

Issued & Understood by me,

Deena B. B. B.

Date

1/6/95

Invented by

Recorded by

To Page No. _____

Date

7-94
12-9-94