

80

Project No. _____

Book No. _____

TITLE 23 mer degradation: V, OV, Tr.c
buffers: Cheng vs. Vent vs. KlenTag

From Page No. _____

① ② ③ ④ ⑤ ⑥ ⑦ ⑧ ⑨ ⑩ ⑪

Cheng buffer 5X	20	→									
10X KlenTag buffer *			10	→							
Vent buffer						10	→				
Tog storage buffer			2	2	-						
Mg OAc 12 mM	9.5	→									
Mg SO ₄ 100 mM			1.2	μl	→						
glycerol 50%										16	→
DMSO 100%											
32P 23mer **	3	μl	→								
Vent pol 0.05 μl	2			2		2			2		
Deep Vent 0.05 μl		2			2			2			2
Tne 0.5 μl			2			2			2		
H ₂ O	65.5	→	71.8	71.8	73.8	75	→			69	→

Preheat to 70°C, start by addition of DNA pol
remove 10 μl to 5 μl cycle reg stop mix at 10, 20, 30 min
well #1 is 23mer uncut
Vp = 100 μl

T Page

Witness d & Underst od by me,

Date

11/29/94

Invented by

Record d by

Dat

11-4-94

ag N . _____
 (14) (15)

- ✓
- ✓
- ✓
- ✓ ← (note KlenTag system relies on Tag storage buffer for glycerol and Tween-20/NP40 - for Tne, it is diluted in Tag storage buffer so no supplement is needed for vent and Deep vent dilution is in storage buffer (with Triton and 50% glycerol)
- ✓ (1.2 mM Mg OAc Cf)
- ✓ (1.2 mM Mg SO₄ Cf)
- ✓ Cf = 8% glycerol
- ✓ Cf = 2% DMSO
- ✓

2 } dilute in vent/Deep vent storage/dilution buffer (its 2.1% Triton)
 2 (dilute in Tag storage buffer)
 → ✓ Cf = .002% Triton will include 2 µl Tag storage buffer next time (similar to TFL storage buffer with 0.5% Tween/NP40)

3 µl, 0.66 ^{pmol} 13.5 µl (8.91 pmol)
~~23 µl~~ 16.8 µl (25.1 pmol)
 6 pmol H₂O 24.7 µl
 5.5 µl
 0.36 pmol primer

** for 72p 23 mer mix 0.66 pmol 13.5 µl
 plus 16.8 µl cold 5' 3' 23 mer plus
 24.7 µl H₂O. Cf = 5.5 µl and specific activity is reduced 4x 2x
 * 10x KlenTag is 500 mM Tris HCl pH 7.0 160 mM (NH₄)₂SO₄ and no Mg SO₄

Cf = 360 mM primer

To Page No. _____

sed & Understood by m , Benedicta Polansky	Date 11/29/94	Invented by 	Dat 11-4-94
Recorded by			