

<sup>32</sup>P primer for 14/ Vent  
Human spleen DNA

ag N — <sup>32</sup>P 2633 (into the anchor primer)  
follow P. 53 except use more <sup>32</sup>P ATP

~26<sup>th</sup> primer would ATP is 100% effective in labeling

reagent	2633	159 μM	1	25	675	0.25 μl	33.75
<sup>32</sup> P γ ATP	6000 Ci/mmol						
10 mCi/μl	10-21-94						
(1.67 μM ATP)							
5X Kinase buffer							
PNK 50 <sup>u</sup> /μl							

(159 μM primer) (41.8 μM ATP)

Any down 11C6 ladder 10 μl H<sub>2</sub>O 1 μl <sup>32</sup>P dGTP 15' 37°C 1 μl EDTA

37°C 30 min → 5' 55°C → add

spin col same as P154, 7, and 145, 3


dilute <sup>32</sup>P 2633 with 100 μl H<sub>2</sub>O (V<sub>p</sub> = 133 now)

spin in microfuge in "micron 3"


(amicon # 4240?) - after all went in, put

add 200 μl more H<sub>2</sub>O and spin again

remove volume that did not enter filter



invert filter



10-24-94

Had a problem: filter kept peeling back on micron 3. Maybe g force was too high on Beckman microfuge "E" model will skip separation of free ATP.

<sup>32</sup>P 2633 is diluted only 33.75 fold for C<sub>f</sub> = 4.71 μM

To Page No. \_\_\_\_\_

Reviewed & Understood by me, Deanna Pokany	Date 10/24/94	Invented by [Signature]	Date 10-19-94 10/24/94
		Recorded by	