

10.5 kb

age No. _____

10/95: To try another miniprep - 10.5 kb fragment contained in pDELTA1 - and amplify

used overnight cultures grown in the presence of Tet^r + Kan^r saved at 4°, over 60% and.

Alkaline lysis protocol. all resuspended in 25 µl of TE extracted from 3 x 1.5 ml culture

Did an enzyme titration } - amplification done only with miniprep DNA. no plasmids were kept yet.

included was just 1 rx with Tag & 20

Vol 50 µl. 200 µM dNTP 4, 3, 2 mM Mg 0.4 µM primers { 1, 2, 5 U enzyme } 0.1 0.1 0.5

cycling 2 step 95°, 3' (95°, 45" / 55°, 30" / 72°, 5') 25 cycles

separated primers with buffer B containing 2 mM Mg added supplement Mg accordingly.

	2 mM	3 mM	4 mM
sed. 2 µl of miniprep template	0	5	10
con. unknown.	50	45	40
separated 20x premix	50	50	50

	100	20	4	4	40	5	5	50
x buffer								
dNTP								
P1								
P2								
temp								
enzyme								
Mg								

5 µl rx as needed. enzyme 1 µl added 5 + 2 = 7 later det. con.

To Page No. _____

Read & Understood by me,

Date

Invented by

Date

Recorded by

1/10/95


A. S. S. S. S.

From Page No. _____

1	2		2
3	4	10	3
5	6		4
7	8		2
9	10	20	3
11	12		4
13	14		2
15	16	50	3
17	18		4

19 Tag 20 2 ml

1		2		5				
2	3	4	2	3	4	2	3	4

Result:

cycling has to be optimized. - Lab of sequencing
 2 mM Mg didn't work in any of the sets. - It
 worked earlier in 25 µl 20 µl Tag.

- amount of template?

- Increasing the enzyme didn't seem to work.
 so does Mg

- Tag alone at 20 / 50 µl didn't work

- get fresh enzyme.

To Page

Witnessed & Understood by me,

Date

1/25/95

Invented by

Recorded by

Dr. Subramaniam

Date

1/10/94