

Project No. _____

Book No. _____

TITLE _____

118

12/1/94


From Page No. _____

Purpose: Since 1U of Tag - titration with Mg didn't work, increase the amount of enzyme to 2U.

- Both new + old dV primers were tried.

- Expt was done under same conditions as I.

old dV 0 1.5 2.5 2.5 2.5 3.0 mM Mg



new dV

Result:

- mispriming persists

- 2U better than 1U of

- 2U / 1.5 mM Mg ✓

- again new dV for cloning work - can give up.

* Try again with new dV primers + increasing amount Tag + Mg.

To Page 1

Witnessed & Understood by me,

Date

Invented by

Date

[Signature]

12/1/94

R c rded by

[Signature]

12/1/94