


Project No. \_\_\_\_\_  
Book No. \_\_\_\_\_

TITLE Completion of RNase Assay -

122


From Page No. \_\_\_\_\_

Take samples from -20°C freezer - spin in micro centrifuge  
15 minutes -  
decant ethanol - air dry pellets -  
Add - ml of RNA blue juice - heat 30 sec at 90°C  
Run out on 16%  
sequencing gel -  
400 volts -



Conclusion -  
Appears to be  
RNase free! Next  
time use more RNase -  
Only used half of  
recommended amount  
used 1ug v- recommended  
2ug.

Bradford on Pools



Witnessed & Understood by me,

May Longo

Date

4/13/95

Invented by

E. Flynn

Recorded by

Date

04/04/95

To Page N

# Exonuclease Assay - The Pool

Project N \_\_\_\_\_

Book No. \_\_\_\_\_

123

Page N \_\_\_\_\_

C. NO. 30042 SOP.

Tube	Rxn Mix 4.5	Enzyme Units	$\mu$	H <sub>2</sub> O
1	↓	0.0	-	5 $\mu$
2		10.0	←	4 $\mu$ 50 $\mu$ /ul
3		20.0	1	
4		10	2	
5		15	3	
6		20	4	
7		0		5 $\mu$ dil'n buffer

Rxn Mix	16 rxns -
10x PCR	80
50mM MgCl <sub>2</sub>	80
S' ds sub	16 pmol 32 $\mu$ 5 pmol/ $\mu$
B' ds sub	16 pmol 32 $\mu$ 5 pmol/ $\mu$
H <sub>2</sub> O	494
	720

Heat @ 37°C for 1 hour - 1-7  
Heat @ 72°C for 1 hour - 8-14

see page - 124 for data

To Page No. 124

Used & Understood by me,

- Mary Longo

Date

4/15/95

Investigated by

E. Flynn

Recorded by

Date

04/06/95

From Page No. \_\_\_\_\_

Rxn Mixture - in 8 rxns. -

(all tubes on ice before use -)

15264-03  
.34  $\mu$ g/ $\mu$ l10x PCR buffer -  
50mM MgCl<sub>2</sub> -  
→  $\phi$ X174  $\Phi$  DNA -  
Autoclaved H<sub>2</sub>O40  $\mu$ l ✓  
40  $\mu$ l ✓  
8  $\mu$ g (23.5  $\mu$ l) ✓  
256.5360  $\mu$ l

Endo mix

H<sub>2</sub>ODiluted enzyme 50/ $\mu$ 

1

45

5

2

45

1

2 units - 2  $\mu$ l

3

45

4

5 units - 1  $\mu$ l

4

45

3

10 units 2  $\mu$ l

5

45

2

15 units 3  $\mu$ l

6

45

1

20 units 4  $\mu$ l

7

45

5 Dil Buffer ✓

Incubate @

72°C

in 3 hours -

or

37°C

5.5 hours

Tag

Double Strand. Assay -

25260-027

EF 1702

.33  $\mu$ g/ $\mu$ l

10x PCR buffer

40 ✓

50mM MgCl<sub>2</sub>

40 ✓

-  $\phi$ X174  $\Phi$  RF8  $\mu$ g (24.2  $\mu$ l) ✓Autoclaved H<sub>2</sub>O

255.8

360 -

Endo

H<sub>2</sub>ODil. Enzyme 50/ $\mu$ 

1

45

5

2

4  $\mu$ l of 50/ $\mu$ l

2

45

1

5

1  $\mu$ l

3

45

4

10

2  $\mu$ l

4

45

3

15

3  $\mu$ l

5

45

2

20

4  $\mu$ l

6

45

1

-

T Pag 1

7

45

5 Dil Buffer ✓

Witnessed &amp; Understood by me,

Date

Invented by

Date

Manny Longo

4/13/95

Rec rded by

04-01-95

Pag No. \_\_\_\_\_

Spin samples down add 5  $\mu$ l of Blue Juice -  
Run out on 1.2% Agarose gel -

1 2 3 4 5 6 7 8 9 10 11 12 13 14



SS Endo

DS Endo

1 2 3 4 5 6 7 8 9 10

H<sub>2</sub>O 2 5 10 15 20 B

8 9 10 11 12 13 14

H<sub>2</sub>O 2 5 10 15 20 B

C = 100  
100

at 100 - 45  
10  
45

conv  
3.4  
10%

Endo looks good - however DS Endo - shows conversion to linear but this is also present in the buffer only lane - could just be a contaminant in the Dil'n buffer -

Dil'n Buffer used - from A.G. flasks from the 4°C Deep cooler - orange tip -

Conclusion: - free of SS Endo nuclease - possible <sup>some</sup> DS endo nuclease but control w/ buffer only shows significant conversion to linear so believe that the dil'n buffer in or has <sup>DS endo</sup> activity not the enzyme prep.

To Pag No. \_\_\_\_\_

Assessed & Understood by m ,

Date

Invented by

Date

Recorded by

Mary Tonger

4/13/95

S. Figure

5/10/95