

Project No. \_\_\_\_\_


Book No. \_\_\_\_\_

TITLE Completion of RNase Assay -

122


From Page No. \_\_\_\_\_

Take samples from -20°C freezer - spin in micro centrifuge  
15 minutes -  
decant ethanol - air dry pellets -  
Add -  $\mu$ l of RNA blue juice - heat 30 sec at 90°C  
Run out on 16%  
sequencing gel -  
400 volts -



Conclusion -  
Appears to be  
RNase free! Next  
time use more RNA -  
Only used half of  
recommended amount  
used 1  $\mu$ g v- recommended  
2  $\mu$ g.

Bradford on Pools



Witnessed & Understood by me,

May Longo

Date

4/13/95

Invented by

E. Flynn

Recorded by

Date

04/06/95

To Page N



# Exonuclease Assay - The Pool

Project N \_\_\_\_\_

Book No. \_\_\_\_\_

123

Page N \_\_\_\_\_

C. NO. 30042 SOP.

Tube	Rxn Mix 4.5	Enzyme Units	$\mu$	H <sub>2</sub> O
1	↓	0.0	-	5 $\mu$
2		10.0	←	4 $\mu$ 50 $\mu$ /ul
3		20.0	1	
4		10	2	
5		15	3	
6		20	4	
7		0		5 $\mu$ dil'n buffer

Rxn Mix	16 rxns -	
10x PCR	80	
50mM MgCl <sub>2</sub>	80	
S' ds sub	16 pmol	32 $\mu$ 5 pmol/ $\mu$
B' ds sub	16 pmol	32 $\mu$ 5 pmol/ $\mu$
H <sub>2</sub> O	494	
	720	

Heat @ 37°C for 1 hour - 1-7  
Heat @ 72°C for 1 hour - 8-14

see page - 124 for data

To Page No. 124

Used & Understood by me,

May Longo

Date

4/15/95

Investigated by

E. Hyman

Recorded by

Date

04/06/95



Project No. \_\_\_\_\_

Book No. \_\_\_\_\_

TITLE Endo Assay - 18038 QCP-T

From Page No. \_\_\_\_\_

Rxn Mixture - fn 8 rxns -

(all tubes on ice before use -)

15264-03  
.34 µg/µl10x PCR buffer -  
50mM MgCl<sub>2</sub> -  
→ ΦX174 Ⓢ DNA -  
Autoclaved H<sub>2</sub>O40 µl ✓  
40 µl ✓  
8 µg (23.5 µl) ✓  
256.5

360 µl

Endo mix

H<sub>2</sub>O

Diluted enzyme 50/µl

1	45	5	-
2	45	1	2 units - 2 µl
3	45	4	5 units - 1 µl
4	45	3	10 units 2 µl
5	45	2	15 units 3 µl
6	45	1	20 units 4 µl
7	45	5	Dil Buffer ✓

Incubate @ ~~72°C~~ 37°C fn 3 hours -  
5.5 hours

Tag

Double Stranded Assay -

25260-027

EF 1702

.33 µg/µl

10x PCR buffer  
50mM MgCl<sub>2</sub>  
- ΦX174 Ⓢ RF  
Autoclaved H<sub>2</sub>O40 ✓  
40 ✓  
8 ~~24.2~~ 24.2 ✓  
255.8

360 -

Endo

H<sub>2</sub>O

Dil. Enzyme 50/µl

1	45	5	-
2	45	1	2 4 µl of 50/µl
3	45	4	5 1 µl
4	45	3	10 2 µl
5	45	2	15 3 µl
6	45	1	20 4 µl
7	45	5	Dil Buffer ✓

Witnessed &amp; Understood by me,

Date

Invented by

Date

Manny Longo

4/13/95

Rec rded by

04-04-95

T Page 1



Pag No. \_\_\_\_\_

Spin samples down add 5ul of Blue Juice -  
Run out on 1.2% Agarose gel -

1 2 3 4 5 6 7 8 9 10 11 12 13 14



SS Endo

DS Endo

1 2 3 4 5 6 7 8 9 10

H<sub>2</sub>O 2 5 10 15 20 B

8 9 10 11 12 13 14

H<sub>2</sub>O 2 5 10 15 20 B

C = 100  
100

at 100 - 45  
10  
45

conv  
3.4%  
10.1%

Endo looks good - however DS Endo - shows conversion to linear but this is also present in the buffer only lane - could just be a contaminant in the Dil'n buffer -

Dil'n Buffer used - from A.G. flasks from the 4°C Dec. cooler - orange bp -

Conclusion: - free of SS Endo nuclease - possible <sup>some</sup> DS endo nuclease but control w/ buffer only shows significant conversion to linear so believe that the dil'n buffer in or has <sup>DS endo</sup> activity not the enzyme prep.

To Pag No. \_\_\_\_\_

Assessed & Understood by m ,

Date

Invented by

Date

Recorded by

Mary Tonger

4/13/95

S. Figure

5/10/95